


Sep 18, 2025

🌐 Genotyping Mice by PCR and Copy Number Standards

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5krynv1b/v1>



Thom Saunders¹, Honglai Zhang¹, Jennifer Leo¹, Sivakumar Jeyaraja¹, Eden A. Dulka¹, Zachary T. Freeman¹

¹University of Michigan

Thom Saunders: Orcid ID: 0000-0003-2015-101X

Sivakumar Jeyaraja: Orcid ID: 0000-0001-6614-5489

Eden A. Dulka: Orcid ID: 0000-0002-7724-8923

Zachary T. Freeman: Orcid ID: 0000-0003-1291-382X

Transgenic Animal Mod...



Thom Saunders

University of Michigan

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5krynv1b/v1>



External link: <https://michmed.org/z45g9>

Protocol Citation: Thom Saunders, Honglai Zhang, Jennifer Leo, Sivakumar Jeyaraja, Eden A. Dulka, Zachary T. Freeman 2025. Genotyping Mice by PCR and Copy Number Standards. [protocols.io https://dx.doi.org/10.17504/protocols.io.8epv5krynv1b/v1](https://dx.doi.org/10.17504/protocols.io.8epv5krynv1b/v1)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 28, 2025

Last Modified: September 18, 2025


Protocol Integer ID: 225753

Keywords: genomic mouse dna, sensitivity of the genotyping assay, experimental artifact pcr assay, genomic copy, genotyping mice, genotyping assay, mice by pcr, transgene dna, crispr founder, amount of transgene dna, transgenic founder, other genomic dna, inserted dna, copy number standards the purpose, pcr, important because transgenic founder, pcr screen, use of copy number standard, copy number standard, identification of transgene, insensitive assay, new primer pair, standardized copy, specific primer, transgene, copy standard, primer

Abstract

The purpose of this protocol is to describe the identification of transgenes in genomic mouse DNA. The use of copy number standards allows one to determine the sensitivity of the genotyping assay. This is important because transgenic founders or CRISPR founders are often mosaic and may carry less than one genomic copy of the inserted DNA when tail tip or other genomic DNA is assayed. PCR screens must be designed to detect transgene DNA at 0.1 copy level or lower. Insensitive assays are known to be the reason that transgenic founders were not detected on numerous occasions. To avoid this experimental artifact PCR assays need to be standardized copy standard template prepared by mixing non-transgenic tail DNA with a known amount of transgene DNA. When sensitive and specific primers are used then the occurrence of false negative results during genotyping can be eliminated. If primers do not meet sensitivity requirements or if they produce more than one band on a gel with copy standards then new primer pairs are ordered and tested until a suitable pair is identified.

Image Attribution



Copy sensitivity test with C468 KO spanning PCR primers on synthetic template spiked into wild type genomic DNA. Note that the C468 KO primers can detect as little as 0.01 copies of the transgene template when it is present in genomicDNA.

- Lane 1: Invitrogen 1 Kb Plus Ladder.
- Lane 2: 100 transgene copy standard.
- Lane 3: 50 transgene copy standard.
- Lane 4: 10 transgene copy standard.
- Lane 5: 1.0 transgene copy standard.
- Lane 6: 0,1 transgene copy standard.
- Lane 7: 0.01 transgene copy standard.
- Lane 8: 0.001 transgene copy standard.
- Lane 9: 0.0001 transgene copy standard.
- Lane 10: Wild type DNA (zero copy transgene).
- Lane 11: Negative control. No template.
- Lane 12: Negative control. No template.
- Lane 13: Invitrogen 1 Kb Plus Ladder.

Image provided by Jennifer Leo.

Safety warnings

 All activities involving animals require approval by an IACUC or other relevant regulatory body prior to executing this protocol.

Ethics statement

All animal activities described were reviewed and approved by the University of Michigan IACUC.

Procedure: Prepare Copy Number Standards

- 1 Assumption: the length of the inserted transgene DNA is 5480 bp.
- 2 Assumption: the Haploid content of a mammalian genome is 3×10^9 bp.
- 3 Assumption: 2 μ g of mouse tail DNA are available.
- 4 Since the transgenic founder mice are hemizygous or mosaic calculate the amount of transgene DNA that would occur in 1 microgram of genomic DNA in a one copy mouse:
- 5 Divide the number of base pairs in the transgene by 3×10^9 bp or
 $5,480 \text{ bp} \div 3 \times 10^9 \text{ bp} = 1.83 \times 10^{-7}$
- 6 For a mass of 1 microgram of DNA you need to add
 1.83×10^{-7} micrograms of transgene DNA to 2 micrograms of DNA to make a copy standard of one transgene copy per hemizygous genome.
- 7 Thus, for a 10 copy standard 18.3 pg transgene DNA is mixed with 2 μ g tail DNA
Thus, for a 1 copy standard 1.83 pg transgene DNA is mixed with 2 μ g tail DNA
Thus, for a 0.1 copy standard 0.183 pg transgene DNA is mixed with 2 μ g tail DNA
Thus, for a 0.01 copy standard 0.0183 pg transgene DNA is mixed with 2 μ g tail DNA

Suggestions for Primer Design

- 8 Primer design suggestions are taken from the Mouse Blastocyst DNA Extraction and PCR Amplification protocol.
DOI [dx.doi.org/10.17504/protocols.io.5qpvodr5bg4o/v1](https://doi.org/10.17504/protocols.io.5qpvodr5bg4o/v1)
- 9 Procedure: Primer Design for Specific and Sensitive PCR Assays.
- 10 Use Primer-Blast to pick primers <http://www.ncbi.nlm.nih.gov/tools/primer-blast>
- 11 Adjust Primer Parameter default settings.
Minimum primer melting temperature: change to 60°C.
Optimal primer melting temperature: change to 63°C.



Maximum primer melting temperature: change to 66°C.

Minimum primer melting temperature difference: change to 1°C.

Adjust Primer Pair Specificity Checking Parameters

Click box to turn on "Enable search for primer pairs specific to the intended PCR template"

Set Search Mode to "Automatic"

Set Database to "Genomes for selected eukaryotic organisms (primary assembly only)"

Set Organism to "Mus musculus (taxid:10090)" for mouse.

Open Advanced Parameters

Change Primer Minimum size to 14

Change Primer Optimum size to 29

Change Primer Maximum size to 31.

- 12 nota bene. Stratman et al. reported primers of 27-30 nucleotides made up of 50-60% GC content will produce a 100-500 bp PCR product that uniformly detects genomic DNA with single copy sensitivity.

Stratman JL, Barnes WM, Simon TC. 2003. Universal PCR genotyping assay that achieves single copy sensitivity with any primer pair. Transgenic Res. 12:521-522.

Isolate Mouse Tail Tip DNA (or other DNA source)Untitled section

- 13 Prepare Reagent: TNES Buffer

A	B	C
Constituent	Quantity	Final Concentration
1 M Tris-HCl, pH 7.5	1.0 ml	10 mM Tris, pH 7.5
5 M NaCl	12.5 ml	400 mM NaCl
500 mM EDTA	20.0 ml	100 mM EDTA
10% SDS	6.0 ml	0.6% SDS
Mol. Bio. Water	60.5 ml	
Total Volume	100 ml	

Produces 100 ml of Buffer. Store at Room Temperature.

14 Prepare Reagent: Proteinase K Stock Solution.

A	B	C
Constituent	Quantity	Final Concentration
Proteinase K	10 mg	
Mol. Bio. Water	1 ml	10 mg/ml

Produces 1 ml of 10 mg/ml Proteinase K. Store at -20°C.

15 Prepare Reagent: Saturated NaCl.

A	B	C
Constituent	Quantity	Final Concentration
NaCl	350.6 g	
Mol. Biol. Water	1000 ml	6 molar

Store at 37°C to aid NaCl solubility, salt crystals will be apparent.

16 Obtain tail biopsies from 2- to 3-week-old mice in 1.5 ml microtubes per your approved animal use protocol. Hold mouse firmly at base of tail with one hand, with the other cut off 0.5 of the tail tip with a scalpel or single edge razor blade. Ear punches from older mice can be used in place of tail tips. Keep tissues frozen at -80°C if they will not be immediately processed.

17 Add 600 microliters of TNES and 35 microliters Proteinase K.

18 Incubate overnight (8-24 hr.) at 55°C.



- 19 Add 166.7 microliters 6M NaCl. Shake vigorously for 15 sec.
- 20 Remove supernatant to new 1.5 ml microtube and add 1 – 2 volumes cold 95% ethanol.
- 21 Insert closed end glass capillary into tube and twirl to spool DNA onto the glass capillary. If DNA does not spool then spin down in a microcentrifuge at top speed for 5 minutes and work with DNA pellet.
- 22 Rinse spooled DNA by dipping glass capillary in a tube filled with 70% EtOH and allow to air dry 5-10 min.
- 23 Place capillary with DNA into a microtube containing 100-500 µl of TE buffer (10mM Tris, pH 8.0, 1mM EDTA). TE volume depends empirically on pellet size.
- 24 Heat at 65°C for 10 minutes to aid dissolution of DNA.
- 25 Quantitate DNA and store at 4°C until needed.
- 26 Method adapted from: Miller SA. Dykes DD. Polesky HF. 1988. A simple salting out procedure for extracting DNA from human nucleated cells. Nucleic Acids Research. 16:1215.

PCR Tests of Primers for Sensitivity and Specificity.Untitled section

- 27 Use 200 ng of the transgene DNA copy standards DNA described above as a template in 25µl PCR reactions as described below.
- 28 Prepare Reagent: Prepare PCR Cocktail

	A	B	C
	Constituent	Quantity For 1 Reaction	Quantity For 20 Reactions
	Mol. Biol. Water	14.0 µl	280.0 µl
	10X PCR Buffer	2.5 µl	50.0 µl



	A	B	C
	10 millimolar Nucleotides	4.0 µl	80.0 µl
	Primer A (10pmol/microliter)	1.25 µl	25.0 µl
	Primer B (10pmol/microliter)	1.25 µl	25.0 µl
	Taq Polymerase (0.3125 U/µl)	1.0 µl	20.0 µl
	Total Volume	24.0 µl	480.0 µl

29 Aliquot 24.0 µl of cocktail into PCR tubes.

30 Add 1.0 microliter of DNA to each PCR tube.

nota bene: if tail DNA concentration is less than 200 ng/µl then add 200 ng of DNA in a volume up to 3 µl.

31 Final volume of PCR reaction is 25.0 µl.
Overlay with mineral oil and place into thermal cycler.

32 Amplify 30 to 35 cycles.
The first step in the cycle is denaturation for 30 seconds at 98°C.

33 Analyze on a 1 to 1.5% agarose gel in TBE, add 0.5 µg/ml of ethidium bromide to the gel prior to casting the gel. Visualize DNA on a UV transilluminator.

34 Remember to include a negative control of water only with PCR reagents.

35 Remember to include a DNA 0.1 copy standard as a positive control.
The first time primers are used they should be tested on a range of copy standards to demonstrate their sensitivity. Sensitive primers will detect as little as 0.001 transgene copies per genome.
If more than one band is produced discard the primers and order new primers.



If the primers can not detect 0.1 copies of the transgene in the copy standard discard the primers and order new.

36 For every sample run a parallel PCR with internal positive control primers that amplify genomic DNA from any mouse on the face of the earth. This will ensure that the DNA sample is "amplifiable" and free of PCR inhibiting contaminants. For example, amplification of beta globin indicates that the tail DNA is free of PCR inhibitors. If it is positive for beta globin and negative for the expected PCR product expected from sensitive and specific primers this indicates the animal is not transgenic. Parallel PCR reactions or an internal positive control prevents false negative results, and positive animals will not be accidentally discarded as a consequence of PCR inhibitors in the DNA preparation.

37 Beta globin primers are used to amplify the endogenous mouse beta globin gene. Amplification will thus demonstrate that the DNA sample is an "amplifiable" substrate for PCR.

5' CCA ATC TGC TCA CAC AGG ATA GAG AGG GCA GG 3' 32 mer

5' CCT TGA GGC TGT CCA AGT GAT TCA GGC CAT CG 3' 32 mer

Size of expected amplification product is 494 base pairs.

38 Thermal cycler profile:

Heat at 98°C for 60 seconds.

Cycles 1 through 35:

94°C 30 seconds

60 °C 90 seconds

72° C 120 seconds


Final treatment:

72°C 10 min

4° C Hold

Results of Copy Standard PCR

39 The PCR products are analyzed by agarose gel electrophoresis.



Copy sensitivity test with C467 5' transgene specific primers on synthetic template spiked into wild type (wt) genomic DNA. Note that the C467 5' transgene primers can detect as little as 0.001 copies of the transgene template when it is present in genomic DNA.

- Lane 1: Invitrogen 1 Kb Plus Ladder.
- Lane 2: 100 transgene copy standard.
- Lane 3: 50 transgene copy standard.
- Lane 4: 10 transgene copy standard.
- Lane 5: 1.0 transgene copy standard.
- Lane 6: 0.1 transgene copy standard.
- Lane 7: 0.01 transgene copy standard.
- Lane 8: 0.001 transgene copy standard.
- Lane 9: 0.0001 transgene copy standard.
- Lane 10: Wild type DNA (zero copy transgene).
- Lane 11: Negative control. No template.
- Lane 12: Negative control. No template.



Protocol references

Konkel DA, Tilghman SM, Leder P. 1978. The Sequence of the Chromosomal Mouse beta-globin Major Gene: Homologies in Capping, Splicing and Poly(A) Sites. *Cell* 15:1125-1132.

Miller SA, Dykes DD, Polesky HF. 1988. A simple salting out procedure for extracting DNA from human nucleated cells. *Nucleic Acids Research*. 16(3):1215.

Ralser M, Querfurth R, Warnatz HJ, Lehrach H, Yaspo ML, Krobisch S. 2006. An efficient and economic enhancer mix for PCR. *Biochem Biophys Res Commun*. 347:747-751.

Stratman JL, Barnes WM, Simon TC. 2003. Universal PCR genotyping assay that achieves single copy sensitivity with any primer pair. *Transgenic Res*. 12(4):521-2.

Ye J, Coulouris G, Zaretskaya I, Cutcutache I, Rozen S, Madden T. 2012. Primer-BLAST: A tool to design target-specific primers for polymerase chain reaction. *BMC Bioinformatics*. 13:134.