



Mar 21, 2025

🌐 Genotyping the SNP rs1143634 of the gene Interleukin-1 β with Real-Time Polymerase Chain Reaction (PCR) in dengue patients

DOI

<https://dx.doi.org/10.17504/protocols.io.6qpvr9ow3vmk/v1>



Raquel da Silva Carvalho¹, Vanessa Rafaela Milhomen Cruz Leite², Jéssica arletto de Sousa Barros¹, Irmtraut Araci Hoffmann Pfrimer¹

¹Postgraduate in Environmental and Health Sciences, Pontifical Catholic University of Goiás (PUC-GO), Center for Immunological Studies and Research (NEPY), Goiânia, Goiás;

²Laboratory of Genetics and Molecular Biology - Institute of Biological Sciences II - Federal University of Goiás, Goiânia, Goiás

Jéssica arletto de Sousa Barros: These authors contributed equally to this work;

Irmtraut Araci Hoffmann Pfrimer: These authors contributed equally to this work



Raquel da Silva Carvalho

Pontifícia Universidade Católica de Goiás

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr9ow3vmk/v1>



Protocol Citation: Raquel da Silva Carvalho, Vanessa Rafaela Milhomen Cruz Leite, Jéssica arletto de Sousa Barros, Irmtraut Araci Hoffmann Pfrimer 2025. Genotyping the SNP rs1143634 of the gene Interleukin-1 β with Real-Time Polymerase Chain Reaction (PCR) in dengue patients . **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6gpvr9ow3vmk/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: November 14, 2024

Last Modified: March 21, 2025

Protocol Integer ID: 112134

Keywords: Dengue; Genotyping; IL-1 β , single nucleotide polymorphism, dengue patient, dengue case, dengue, interleukin, inflammatory response, regulating inflammatory response, immune response, chronic inflammation, adaptive immunity, immune, stimulation of chronic inflammation, snp rs1143634, such as the snp rs1143634

Abstract

Interleukin 1 (IL-1) plays a role in regulating inflammatory responses and is crucial for both innate and adaptive immunity. The secretion of IL-1 β can influence various biochemical, physiological and pathological functions in the body, ranging from polarization towards a Th1 response and intensification of pro-inflammatory mediators in the acute phase, to stimulation of chronic inflammation and carcinogenic potential in some cases. In addition, dengue cases are alarmingly high every year, reaching thousands of cases in Brazil and Latin America. The increase in severe cases of dengue may be related to genetic factors, such as the SNP rs1143634 (G>A) in the interleukin 1- Beta (IL-1 β) gene. Although synonymous, this SNP can influence the immune response, potentially aggravating the infection.

Materials

1. PCR 96- Well Plate (**DNASE FREE**)
2. SEALANT
3. PIPETTE (V°: 2/ 10/ 20 MICROLITERS)
4. TIP (V°: 2/ 10/ 20 Microliters **DNASE FREE** with Filter)
5. DNA SAMPLE - Extracted and Quantified
6. MASTER MIX (THERMO FISHER)
7. TAQMAN GENOTYPING ASSAY (THERMO FISHER)
8. NUCLEASE-FREE WATER

INITIAL PREPARATION

1 Before starting

During the extraction of human DNA, attention must be paid to possible sources of contamination, as well as to the manufacturer's recommendations.




Figure 1: Extraction of human DNA by PureLink™ Genomic DNA Mini Kit (Invitrogen), for future genotype detection.

2 GENERAL ASPECTS

1. Definition:

The polymerase chain reaction (PCR) is a technique used to detect RNA and DNA, with the test focusing on DNA. This technique, also known as quantitative PCR or quantitative real-time PCR, is a method of simultaneously quantifying and amplifying DNA. Real-time PCR is based on the detection and quantification of a fluorescent compound (during the early stages of PCR), based on the basic principle that the significant increase in the amount of PCR product correlates with the initial amount of the control. This technique requires only 3 picograms of material, which is around 1000 times less genetic material than conventional assays.

2. **Objective:** Quantify and amplify a specific DNA sequence (SNP IL-1 β rs1143634), identifying the presence of the SNP.

3  20240603_162117 (1).mp4 9.3MB

PERSONAL PROTECTIVE EQUIPMENT

- Nitrile gloves
- Disposable cloak/sleeved jacket
- Goggles
- Surgical mask
- Gown

4 REAGENTS NEEDED

PCR 96- Well Plate (**DNA'SE FREE**)

1. SEALANT



2. PIPETTE (V°: 2/ 10/ 20 MICROLITERS)
3. TIP (V°: 2/ 10/ 20 Microliters **DNA'SE FREE** with Filter)
4. DNA SAMPLE - Extracted and Quantified
5. MASTER MIX (THERMO FISHER)
6. TAQMAN GENOTYPING ASSAY **C_9546517.10** (THERMO FISHER Scientific)
7. NUCLEASE-FREE WATER

4.1



A	B
REAGENTS	VOLUME PER REACTION
MASTERMIX	10ul
TAQMAN PROBE	1ul
GENOMIC DNA	4ul
NUCLEASE- FREE WATER	5ul
TOTAL VOLUME	20 ul

Table 1: Volume of reagents per reaction for Real-Time Polymerase Chain Reaction to detect genotypes.

5 CYCLING PANEL

PCR was carried out using Roche's LightCycler 480 II thermal cycler.





Roche's LightCycler 480 II thermal cycler

	A	B	C	D
	STAGE	TEMPERATURE	TIME	CYCLE
	PRE-READING	60 C°	30s	1
	ENZYME ACTIVATION	90 C°	5m	40
	DENATURATION	60 C°	30s	1
	RINGING/EXTENSION	60 C°	30s	1
	POST-READING	60 C°	30s	1

DESCRIPTION OF THE STEPS

- 6 **Equipment:** 96-well plate
Functions: To carry out the reaction.
Note: The plate has a specific barcode, so it cannot be reused.
- 7 **Process:** Mix each solution well and centrifuge or Multi Mix (Vortex) for up to 1 minute to remove air bubbles.





- 8 **Equipment:** 96-well optical plate
Functions: Mix the volumes of each reagent in each well of the plate up to the pre-set volume.
- 9 **Equipment:** Close the plate with the sealant, then centrifuge for another 5 minutes to remove the air bubbles until well homogenized. 
- 10 **Equipment:** Store the plate with the samples for up to 72 hours. 

ACTIONS IN CASES OF NON-COMPLIANCE

- 11
 - After identifying the non-conformity, describe it in detail in the laboratory log. Describe the errors and the actions taken to correct them. 
- 12
 - Call technical assistance if the equipment is faulty or out of calibration.

INTERPRETING GENOTYPING RESULTS

- 13 After real-time PCR, when the curves of each sample are analyzed, the genotypes are recognized as Wild (GG), being those that have not mutated, Heterozygous (GA), with an exchange of nitrogenous bases in one chromosome and finally Mutant (AA) with total exchange. In the Roche thermal cycler they are identified by color, Wild (GG) by blue, Heterozygous (GA) by red and Mutant (AA) by green.
This data will be analyzed by statistics using the R studio software.






Figure 3: Result after real-time PCR.



Protocol references

1- OLIVEIRA, R: PCR em tempo real: métodos e aplicações. Universidade de Aveiro, 2010.

<http://hdl.handle.net/10773/7230>

2- MENEZES, M;LIMA, L; MARTINELLO, F: Laboratory diagnosis of SARS-CoV-2 by real time reverse transcriptase-polymerase chain reaction (RT-PCR). Universidade Federal de Santa Catarina (UFSC). 2020

DOI: 10.21877/2448-3877.20200006

3-OLIVEIRA M;WATANABE, A; CESAR, D; CANDIDO, J; LIMA, N; MOREIRA, O:TESTES DIAGNÓSTICOS PARA O SARS-COV-2: UMA REFLEXÃO CRÍTICA. Quím Nova [Internet]. 2022Jun;45(6):760–6. Available from:

<https://doi.org/10.21577/0100-4042.20170895>