


Aug 15, 2025

🌐 GeoMx Digital Spatial Profiler Immuno-oncology (IO) Proteome Atlas (IPA) Assay for Human Pancreas

DOI

<https://dx.doi.org/10.17504/protocols.io.5jyl8qyy8l2w/v1>



Alexandra Cuaycal¹, Dongtao Ann Fu¹, Heather Kates¹, Samuel Ewing¹, Jing Chen¹, Maigan Brusko¹, Martha Campbell-Thompson¹, Clayton E Mathews¹

¹University of Florida

Alexandra Cuaycal: Postdoctoral Associate

Dongtao Ann Fu: Molecular Pathology Core, RRID:SCR_016601

Heather Kates: Data Management Analyst III

Samuel Ewing: Biological Scientist

Jing Chen: Associate Research Scientist

Maigan Brusko: Assistant Professor

Martha Campbell-Thompson: Professor

Clayton E Mathews: Professor

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Alexandra E. Cuaycal

University of Florida

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8qyy8l2w/v1>

Protocol Citation: Alexandra Cuaycal, Dongtao Ann Fu, Heather Kates, Samuel Ewing, Jing Chen, Maigan Brusko, Martha Campbell-Thompson, Clayton E Mathews 2025. GeoMx Digital Spatial Profiler Immuno-oncology (IO) Proteome Atlas (IPA) Assay for Human Pancreas. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.5jyl8qyy8l2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

This protocol was used to spatially profile targeted proteins in human pancreas through the HuBMAP program.

Created: June 12, 2025

Last Modified: August 15, 2025

Protocol Integer ID: 219976

Keywords: pancreas, FFPE, GeoMx, Proteome atlas, IPA, NGS, Nanostring, morphology, insulin, cytokeratin, CD31, spatial proteomics, fluorescence, digital spatial profiler, HuBMAP, RNA sequencing, UF Diabetes Institute, cell morphology marker, tissue spatial information, geomx digital spatial profiler immuno, oncology proteome atlas, nanostring geomx digital spatial profiler, proteome atlas, labeled antibody, rna probe, florida molecular pathology core, defined antibody panel, morphology marker, hybridized rna probe, nanostring geomx, protein quantification via bulk rna, antibody, protein quantification, bioinformatics, antibody panel, antigen, geomx dsp instrument use, oncology, bulk rna, antibodies against antigen, following geomx ngs guideline, marker, paraffin sections from human pancrea, geomx ngs guideline, human pancreas sections within region

Funders Acknowledgements:

NIH

Grant ID: 1R01DK122160

NIH

Grant ID: P01AI042288

NIH

Grant ID: U54DK127823

NIH

Grant ID: 5R01DK123329

Disclaimer

For detailed information about this protocol, please visit Nanostring University:

<https://university.nanostring.com/geomx-dsp-spatial-proteogenomic-protocol>

Abstract

Purpose: This protocol details processing fixed paraffin sections from human pancreas for the Nanostring GeoMX digital spatial profiler (DSP) immuno-oncology proteome atlas (IPA). This technology is based upon a defined antibody panel of 575 markers. Each antibody is conjugated to an RNA probe. The RNA probes are designed with UV cleavable linkers and DSP barcodes. Fluorescent-labeled antibodies against antigens that define specific cell types were used as morphology markers. Morphology markers are used to identify cell types of interest with nuclear counterstaining permitting tissue spatial information. Sample regions of interest (ROI) are defined based on the three morphology markers and each ROI is illuminated with UV light thereby cleaving hybridized RNA probes. Cleaved probes are aspirated via microcapillary and deposited into a unique well of a 96-well collection plate. The collection plate is used for next generation sequencing (NGS). **Scope:** This document was written following GeoMX NGS guidelines (university.nanostring.com) with minor modifications for the University of Florida Molecular Pathology Core. The entire workflow is described for manual slide staining and GeoMX DSP instrument use. Steps for library construction, sequencing, bioinformatics, and data analysis are in outline form only. Library preparation and sequencing are outsourced to the University of Florida ICBR facility. **Expected outcome:** Targeted-protein quantification via bulk RNA-sequencing will be obtained from fixed human pancreas sections within regions of interest and subareas defined by cell morphology markers.

Image Attribution

Human pancreas stained with panCytokeratin, Insulin, and CD31 with Syto83 nuclear marker.

Guidelines

Maintain the DSP instrument according to vendor recommendations including weekly wash buffer changes. Maintain an RNase-free workspace. Validate morphology marker antibodies using assay conditions with antigen retrieval and assay incubation steps but use mock IPA antibody probe solution. Prepare all solutions using DEPC-water and RNase-free supplies. Clean equipment with RNase AWAY, rinse with DEPC-treated water, and air-dry.

Materials

Chemicals (store at  Room temperature except where noted):

Deparaffinization reagents (xylene or substitute, ethanols 100%, 95%)

10xPBS pH 7.4 (PBS, FisherScientific BP399-1 or comparable)

20x SSC (DNase RNase free, Sigma, S6639)

10x Citrate Buffer, pH 6.0 Antigen Retrieval Solution (Sigma, C9999-100ML)

10% neutral buffered formalin (NBF, FisherScientific, 50-980-502, EMS 1574004)

Formamide (deionized or molecular grade, FisherScientific, AM9342) store at  4 °C

DEPC-treated water (DEPC water, FisherScientific, BP561-1)

RNase AWAY (FisherScientific, 7003PK)

Tween 20 (FisherScientific, BP337-100)

TRIS-base (Sigma, 10708976001 or equivalent)

Glycine (Sigma G7126 or equivalent)

Paraformaldehyde, 4% in PBS (ThermoFisher, J61899.AK)

Solution preparations:

1xPBS (PBS, 1 liter) - add 100ml 10xPBS to 900ml DEPC water, store at  4 °C .

1xAntigen Retrieval Solution (100ml) - add 10ml 10xCitrate Buffer, pH 6.0 to 90ml DEPC water for final 1x Citrate Buffer, pH 6.0.

2xSSC (1 liter) - add 100ml 20x SSC and 900ml DEPC water.

2xSSC/T (100ml) - add 10ml 2x SSC, 1ml 10% Tween-20, and 89ml DEPC water.

4xSSC (1 liter) - add 200ml 20x SSC and 800 ml DEPC water.

2x SSC/50% formamide - add 100ml each 4x SSC and 100% formamide.

Morphology markers:

Primary antibodies are conjugated to fluorophores for cell types of interest based on fluorophores compatible with the DSP instrument (AF488, AF532, AF594, AF647 or similar). Use up to 3 conjugated primaries with suitable nuclear marker. Nuclear markers are SYTO 83 (AF532, S11364, FisherScientific, 5mM in DMSO). Antibodies used: Insulin, Mouse anti-all, AF488 (NICBTACLS) FisherScientific 53-9769-82. Cytokeratin (PanCK, Mouse anti-all, AF594 NAE/1/AE/3) Novus Biologicals NBP2-33200AF594. CD31, Mouse anti-hu, AF647 (NJC/70A) Abcam ab215912.

	A	B	C	D
Antibody Name:	Insulin	PanCK	CD31	
Fluorophore:	AF488	AF594	AF594	
Clone:	(ICBTACLS)	(AE-1/AE3)	(JC/70A)	
Source:	mouse	mouse	mouse	



	A	B	C	D
	Dilutions:	1:200	1:400	1:40
	Company:	Invitrogen	Novus	Abcam
	Cat #:	50-112-4641	NBP2-33200AF594	Ab215912

Morphological markers used for human pancreas.

Lab Supplies:

Superfrost Plus slides (FisherScientific, 22-035813)

Glass desiccator for slide drying and storage

Kimwipes

RNase-free beakers and cylinders to prepare buffers

Coplin jars or slide trays

Pipettes (5-1000, 12-channel P20) and filter tips (DNase/RNase free)

Reagent reservoir, 20ml, sterile (FisherScientific, 8094)

Microcentrifuge tubes (DNase/RNase free)

Microcentrifuge rack

Adhesive PCR plate foils (ThermoFisher, AB0626).

Nanostring GeoMX DSP IPA reagents and kits:

GeoMX collection plate (4pk, 100473)

GeoMX instrument buffer kit (100474)

GeoMx Protein Slide Prep Kit (121300312)  4 °C

GeoMX Pro Code Pack Z and Y  -20 °C (121400207)

GeoMx Human IO Proteome Atlas Kit (121300160)  -80 °C .

Equipment:

Oven for drying slides

Leica XL slide autostainer or alternative

TintoRetriever pressure cooker (BioSB, MM 0023 - 090121 or alternative)

Waterbath (FisherScientific, FSGPD02)

Slide incubation chamber (black or opaque)

Nanostring GeoMX DSP instrument

PCR thermocycler

Centrifuge with rotor for 96-well PCR plates (at least 2000 x g)

Minicentrifuge for quick spins.



Safety warnings

- ⚠ Follow Environmental Health and Safety protocols for all hazardous chemicals (xylene, formalin, formamide).

Slide preparation

3h 50m

- 1 Cut paraffin block to fresh sample. Section at $\sim 5\ \mu\text{m}$ and place sections centrally on Superfrost Plus slides as shown in Fig. 1. Tissue samples should have maximum width and length of ~ 14 and $36\ \text{mm}$, respectively, for inclusion within the GeoMX instrument slide holder gasket. Tissue outside the gasket region should be carefully removed without creating folds or leaks to gasket seal. Air dry mounted slides overnight before use or storage. Recommended storage is for < 2 weeks in a desiccator at $4\ ^\circ\text{C}$.




Figure 1. Slide dimensions for GeoMx DSP. Tissue is placed in green scan area.

- 2 Bake slides at $60\ ^\circ\text{C}$ for 00:30:00 to 03:00:00 immediately before use.
- 3 Dewax and rehydrate sections using slide tray stations or a Leica XL slide autostainer. Do not allow tissue sections to dry once rehydrated and avoid touching tissue when drying around edges.
 - 3.1 Xylene 3 x 00:05:00. CitriSolv or another xylene substitute can be used.
 - 3.2 100% Ethanol 2 x 00:05:00.
 - 3.3 95% Ethanol 2 x 00:05:00.

3h 30m

5m

5m

5m



3.4 ddH₂O 2 x  00:05:00 .

5m

Antigen retrieval

1h 45m

4 Prepare Antigen Retrieval Solution (1x Citrate Buffer, pH 6) in Coplin jar.

5 Transfer slides to Antigen Retrieval Solution jar and put into BioSB pressure cooker.

6 Secure the pressure cooker lid and run on high pressure and high temperature for  00:15:00 .

15m

7 Allow pressure cooker to release pressure before opening the BioSB pressure cooker, and cool down for  00:25:00 without the lid.

25m

8 Immediately transfer slides to 1x TBS-T.

9 Wash slides in TBS-T 1x  00:05:00 .

5m

10 Slides can be stored in TBS-T up to  01:00:00 .

1h

Blocking and antibody incubations

17h

11 Draw a closed hydrophobic barrier around the tissue with hydrophobic pen (Vector).

12 Block with  200 μ L of Buffer W per slide for  01:00:00 at  Room temperature in a covered humidity chamber.

1h

13 During blocking, thaw the Protein Probe (antibody) mixes, mix by flicking and spin down. Do not vortex.

14 Prepare a working antibody solution by diluting the Protein Probe (antibody) mixes and Morphology Markers (MM, antibodies for protein markers) into Buffer W according to Table 1.

15 Keep working solution on ice until ready to use:

	A	B	C	D	E	F	G	H
		Pro Core Ab mix (IPA tube 1 of 2)	IPA Module (IPA tube 2 of 2)	MM 1: Ins	MM 2: PanCK	MM 3 CD31	Buffer W	Total volum e
	Channel	-	-	AF488	AF594	AF647	-	
	Dilution (v/v)	-	-	1:200	1:400	1:40	-	
	Volume per sample	8.4 uL x n	52.5 uL x n	1.05 uL x n	0.53 uL x n	5.25 uL x n	142.3 uL x n	210 uL x n

Table 1. Working antibody mix solutions for IPA slide prep (n = number of slides).

16 Remove Buffer W from slides (one at a time) by gently tapping each slide on a clean paper towel, then using a kimwipe to remove excess without touching the tissue.

17 Add  200 µL of working antibody solution to each slide. Cover the humidity chamber and incubate overnight at  4 °C .

16h

Post-fixation for pancreas

40m

18 Wash slides in 1x TBS-T 3 x  00:05:00 .

5m

19 Incubate in 4% PFA for  00:30:00 at room temperature.

30m

20 Wash slides in 1x TBS-T 2 x  00:05:00 .

5m

Nuclear staining

20m



21 Prepare enough of a 1:6000 solution of SYTO83 in 1x TBS-T.

22 Add  200 μL of the nuclear stain solution to each slide and incubate for  00:15:00 at  Room temperature .

15m

23 Wash slides in 1x TBS-T 2 x  00:05:00 .

5m

24 Slides can be stored in 1x TBS-T wash solution up to 1 hour at  Room temperature or up to 6 hrs at  4 °C

25 Scrape off hydrophobic pen barrier on a slide with a razor prior to GeoMs DSP run.

GeoMX DSP run

2h

26 Load slides in DSP within 6 hours of last wash for optimal results.

27 Slides must be stored in dark to avoid cleavage of DSP tags by UV light.

28 Start run by selecting New/Continue Run under Data Collection in the GeoMx DSP Control Center. Include a non-template control (NTC).

29 Clean the bottom of slide with 70% ethanol and place slide in slide holder with slide label towards front of instrument.

30 Complete loading slides then lower slide tray clamp. Cover sections with 6 ml Buffer S. Record slots for each slide.

31 Load a 96-well aspiration plate in DSP. Check solution containers are filled and waste container is empty.

32 Follow prompts in the GeoMx DSP Control Center to identify slides.



- 33 Name slides according to protocol
<https://dx.doi.org/10.17504/protocols.io.n2bvj6dnblk5/v1>. Include UF HuBMAP sampleID and morphology marker abbreviations in order of channels followed by date (MMDDYY). Example: P1-21_DNA_PanCK_INS_CD31_121024.
- 34 Upload unzipped configuration files to DSP kits used for Morphology, Core, and Module kits used in slide preparations (nanosting/dsconfigfiles). Config files must match names of assays used on slides.
- 35 Define scan area for each sample using x- and y-sliders and scan. The minimum tissue slide filters particles and small tissue areas from scan. The sensitivity slider adjusts intensity of tissue identification.
- 36 Draw regions of interest (ROI) using geometric, segmented or other tools. Maximum ROI size is ~650 μm diameter circle.
- 37 Segment ROIs for AOI (area of illumination) for cell types of interest based on positive or negative selection by morphology markers. Minimum of ~ 50 nuclei/AOI is recommended for sufficient probe reads.
- 38 Fluorescently labeled antibodies were used as morphology markers to detect cell types of interest. Anti-Insulin is used to designate beta-cells in the endocrine compartment, anti-CD31 is used to define endothelium, and anti-pan cytokeratin (PanCK) is used to define ductal epithelium in the exocrine compartment. Acinar cells are defined by the lack of staining for CD31 and PanCK in the exocrine ROI.
- 39 Exit scan workspace and approve ROI to start ROI collection.  02:00:00 . 2h
- 40 After ROI/AOI collections are completed, follow prompts to remove slide holder and collection plate.
- 41 Remove buffer from each slide with pipette and dispose. Open clamps and unload slides. Store slides in 2x SSC at  4 °C up to 7 days.
- 42 Clean slide holder as directed and replace in instrument. Follow prompts to close instrument door.
- 43 Fill in the GeoMX Seq Code plate to be used with the plate barcode in the Readout Group Information grid. Click update button and finalize and download readout package buttons to USB drive.

Slide reuse



- 44 Protein assay slides can be restained following UV exposure for 3 minutes to strip tags from bound probes followed by washes using with 2xSSC/T. Proceed with additional staining.

NGS library and sequencing

- 45 Transfer collection plate to the University of Florida ICBR Gene Expression and Genotyping Core or other NGS core facility for library preparation and Illumina sequencing. Refer to Nanostring user manual for additional information.
- 46 Library will be prepared after indexing PCR with pooling and purification.
- 47 Library will be evaluated for quality and fragment size distribution using a BioAnalyzer 2100 (Agilent).
- 48 Core facility will perform Illumina sequencing.

Data analysis

- 49 Core facility will process Illumina FASTQ sequencing files and transfer back to laboratory.
- 50 Share FASTQ files and annotation file with data analyst who will analyze RNAseq data using R scripts for GeoMx DSP workflow (Bioconductor - GeoMxWorkflows). Data analysis may also be performed with GeoMX DSP software on auxiliary server.
- 51 Manage annotations (tags). Add columns to further define regions and cell types from morphology markers.
- 52 Run various testing including QC, data scaling, normalization, background, statistical tests, and ratios.
- 53 Analyze in spatial context using spatial deconvolution software and scRNAseq data.
- 54 Complete analysis and upload data to HuBMAP HIVE.