



May 13, 2025

GFP-Clu-tails purification from Escherichia coli cells

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1>

Andreas Bracher¹, F Ulrich Hartl¹

¹Department of Cellular Biochemistry, Max Planck Institute of Biochemistry



Patricia Yuste-Checa

MPI Biochemistry

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1>

Protocol Citation: Andreas Bracher, F Ulrich Hartl 2025. GFP-Clu-tails purification from Escherichia coli cells. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.n2bvj3w55lk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 26, 2024



Last Modified: May 13, 2025

Protocol Integer ID: 94463

Keywords: ASAPCRN, tails purification from escherichia coli cell, fusion protein gfp, fusion protein, tails purification, escherichia coli cell, purification, gfp, protein, tails from escherichia, clu

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000282

Abstract

This protocol details how to purify the fusion protein GFP-Clu-tails from Escherichia coli.

Attachments



GFP-Clu-tails purifi...

1.7MB

Materials

Buffers

Binding buffer: PBS + [M] 20 millimolar (mM) imidazole

Elution buffer: PBS + [M] 250 millimolar (mM) imidazole




His₆-Ubiquitin-GFP-Clu-tails expression and cell lysis

- 1 Express His₆-Ubiquitin-GFP-Clu-tails in E. coli BL21 (DE3) codon+RIL cells cultured in  4 L LB Medium  Overnight at  18 °C with  0.25 millimolar (mM) IPTG. 
- 2 Centrifuge culture and keep pellet. 
- 3 Re-suspend pellet in  70 mL volume of ice-chilled Binding buffer, add Complete protease inhibitor cocktail (Roche) and  1 millimolar (mM) phenylmethylsulfonyl fluoride (PMSF).
- 4 Lyse cells by ultrasonication in ice bath (15 cycles of  00:00:25 ultrasonication with  00:01:35 intermittent cooling). 
- 5 Clear lysate by centrifugation at  22000 rpm in a JA25.50 rotor at  4 °C .

Ni²⁺-chelating Sepharose metal affinity chromatography

- 6 Load supernatant onto a  8 mL Ni-chelating Sepharose (Cytiva 17-0575-01) column previously equilibrated with binding buffer by gravity flow at  4 °C .
- 7 Wash the column with 5 CV of ice-chilled Binding buffer. 
- 8 Wash the column with 2 CV of ice-chilled PBS +  50 Molarity (M) imidazole. 
- 9 Elute His₆-Ubiquitin-GFP-Clu-tails protein with 6x  5 mL of ice-chilled Elution buffer. Collect fractions of 5 mL volume. Bright green color indicates presence of His₆-Ubiquitin-GFP-Clu-tails. Store fractions  On ice .
- 10 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.




Protease cleavage and removal of protease

- 11 Pool fractions containing His6-Ubiquitin-GFP-Clu-tails peak (here fractions 4 and 5).
Cleave off His6-Ubiquitin with His-tagged Usp2 during ⌚ 04:00:00 🌡️ On ice . 4h
- 12 Exchange protein buffer to Binding buffer with a HiPrep™ 26/10 Desalting column (Cytiva 17-5087-01) equilibrated with Binding buffer at 🌡️ 4 °C using bright green fluorescence as a marker for protein containing fractions.
- 13 Remove uncleaved material and His-tagged Usp2 by metal affinity chromatography as above. Collect flow-through and wash fractions with bright green fluorescence.
- 14 Concentrate to 🧴 1.5 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at 🌡️ 4 °C .

Size exclusion chromatography on HiLoad 16/600 Superdex-200 (Cytiva 28-9893-35)

- 15 Apply concentrate on a HiLoad 16/600 Superdex-200 (Cytiva 28-9893-35) column equilibrated with PBS. Develop the column at $4\text{ }^{\circ}\text{C}$ and collect 3 mL fractions.
- 16 Analyze eluted fractions by SDS-PAGE and Coomassie blue staining.



- 17 Merge fractions 16 and 17. Concentrate to 0.7 mL volume by ultrafiltration using 10 kDa cut-off spin concentrator at $4\text{ }^{\circ}\text{C}$, aliquot and flash-freeze purified GFP-Clu-tails in liquid nitrogen for storage at $-70\text{ }^{\circ}\text{C}$.

Note

Concentrations were determined by absorbance at 488 nm (eGFP!) using absorbance coefficients of $61,000\text{ M}^{-1}\text{ cm}^{-1}$ or $1.45\text{ L g}^{-1}\text{ cm}^{-1}$.

Approximate yield: From 4 L of culture around 2.5 mg of GFP-Clu-tails were obtained.