

Sep 23, 2023

## GFP pull down assay

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwd6x2lmk/v1>

Elias Adriaenssens<sup>1</sup>

<sup>1</sup>Sascha Martens lab, University of Vienna, Max Perutz Labs - Vienna



Elias Adriaenssens

Sascha Martens lab, University of Vienna, Max Perutz Labs - ...

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.e6nvwd6x2lmk/v1>

**Protocol Citation:** Elias Adriaenssens 2023. GFP pull down assay. protocols.io  
<https://dx.doi.org/10.17504/protocols.io.e6nvwd6x2lmk/v1>

**License:** This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** June 27, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 84086

**Keywords:** GFP pull down assay, ASAPCRN, assay this protocol, gfp, assay, protocol

**Funders Acknowledgements:**

Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-000350

Marie Skłodowska-Curie MSCA Postdoctoral fellowship

Grant ID: 101062916

## Abstract

This protocol describes GFP pull down assay.

## Attachments



[754-1922.pdf](#)

43KB

## Materials

### Materials

- GFP-Trap agarose beads (Chromotek)
- dH<sub>2</sub>O
- Protein Loading dye

### Bead assay buffer

|  | A               | B      |
|--|-----------------|--------|
|  | Tris-HCl pH 7.4 | 25 mM  |
|  | NaCl            | 150 mM |
|  | DTT             | 1 mM   |



## GFP pull down assay

2h 5m

- 1 Mix GFP-tagged TBK1 with  20  $\mu\text{L}$  of equilibrated GFP-Trap agarose beads (Chromotek) at a final concentration of  1 micromolar ( $\mu\text{M}$ ). Make sure to wash beads 2x in dH<sub>2</sub>O before washing with bead assay buffer to equilibrate the beads adding the protein to the beads.
- 2 To this end, wash  20  $\mu\text{L}$  of beads twice with dH<sub>2</sub>O and equilibrate with bead assay buffer. 
- 3 Resuspend the beads in  40  $\mu\text{L}$  bead assay buffer, to this add GFP-TBK1 at a final concentration of  5 micromolar ( $\mu\text{M}$ ). 
- 4 Incubate the beads with GFP-TBK1 for  01:00:00 at  4 °C at a horizontal tube roller.   

- 5 Wash the beads three times to remove unbound GFP-tagged bait protein. 
- 6 Prepare protein master mixes with prey protein in bead assay buffer at the following concentrations:
  - mCherry-OPTN (1  $\mu\text{M}$ ),
  - mCherry-NDP52 (1  $\mu\text{M}$ ),
  - GST-NAP1 (1-10  $\mu\text{M}$ ).
- 7 Add the protein master mixes to the beads and incubate for  01:00:00 at  4 °C at a horizontal tube roller.   
 
- 8 Wash the beads three times to remove unbound proteins, remove any supernatant from the beads and resuspend the beads in  60  $\mu\text{L}$  of 1x Protein Loading dye, and heat-inactivate at  95 °C for  00:05:00.   
 
- 9 Analyze the samples by SDS-PAGE and Coomassie staining as described above. 