



Sep 07, 2023

Gibson Assembly Cloning

 [The Journal of Biological Chemistry](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2lyjwyqlx9/v1>

Claire Chiang^{1,2}, Suzanne Pfeffer^{1,2}

¹Stanford University School of Medicine;

²Aligning Science Across Parkinson's (ASAP) Collaborative Research Network

Public-lab



Suzanne R Pfeffer

Stanford University School of Medicine

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2lyjwyqlx9/v1>

External link: <https://doi.org/10.1016/j.jbc.2025.110679>

Protocol Citation: Claire Chiang, Suzanne Pfeffer 2023. Gibson Assembly Cloning . **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2lyjwyqlx9/v1>

**Manuscript citation:**

Adhikari A, Tripathi A, Chiang CY, Sherpa P, Pfeffer SR (2025) Allosteric regulation of the Golgi-localized PPM1H phosphatase by Rab GTPases modulates LRRK2 substrate dephosphorylation in Parkinson's disease. The Journal of Biological Chemistry 301(10). doi: [10.1016/j.jbc.2025.110679](https://doi.org/10.1016/j.jbc.2025.110679)

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 06, 2023

Last Modified: May 31, 2024

Protocol Integer ID: 87471

Keywords: ASAPCRN, gibson assembly cloning, gibson assembly, gibson reaction, appropriate antibiotic, reaction into dh5 alpha, µl reaction

Funders Acknowledgements:

Aligning Science Across Parkinson's Disease

Grant ID: ASAP-000463

Abstract

Gibson assembly requires a vector backbone and one or more inserts that have been PCR amplified. The inserts should have 15-20 base pairs of overlap at ligation sites, so primers used to amplify the inserts should contain a tail overhang with this homology. The vector backbone should be linear.

A Gibson reaction uses a master mix that can be homemade or commercially purchased. Add the appropriate amount of master mix, then your vector and insert at 2-3 fold molar excess for a 20 µl reaction. Exact ratios will need to be optimized based on the size of your vector and inserts and how many inserts you have. Incubate this reaction at 50°C for 60 minutes. Then transform this reaction into DH5 alpha (or competent cell of choice) and plate on LB agar with appropriate antibiotic.



Materials

PEG-8000

Tris-HCl pH 7.5

MgCl₂

DTT

dNTPs [NEB, N0447S]

NAD

10U/μl T5 exonuclease [NEB, M0663S]

2U/μl Phusion polymerase [Thermo, F530S]

40U/μl Taq Ligase [NEB, M0208S]

Sterile water

Competent cells (DH5α) [Thermo, 18258012]

PCR machine

LB, LB agar

antibiotic of choice



Prepare Gibson master mix

- 1 Make 5 ml of 5x reaction buffer
 - 1.1
 - 25% PEG-8000 1.25 g
 - 500 mM Tris-HCl pH 7.5 2.5 ml
 - 50 mM MgCl₂ 0.25 ml of 1M
 - 50 mM DTT 0.25 ml of 1M
 - 1 mM each dNTPs 0.05 ml each of 100mM
 - 5 mM NAD 0.5 ml of 50 mM
 - 2 Prepare master mix
 - 2.1 Master Mix
 - 320 µl 5X Reaction Buffer (above)
 - 0.64 µl 10U/µl T5 exonuclease
 - 20 µl 2U/µl Phusion
 - 160 µl 40U/µl Taq Ligase
 - 700 µl Water
- 3 Use master mix as 1.5x. Snapfreeze in  15 µL aliquots and store at  -80 °C .

Gibson assembly

- 4 Prepare linearized vector backbone by restriction enzyme digestion or PCR amplification.
- 5 Prepare inserts by PCR amplification. Ensure that each insert has 15-20 base pairs of overlap with other inserts/vector at ligation sites, added by primer design.
- 6 Take a 15 ul aliquot of 1.5x Gibson master mix on ice. Add 50-100 ng of vector with 2-3 fold molar excess of PCR inserts, and sterile water if necessary, to make total reaction volume up to  20 µL .
- 7 Mix well and let incubate in PCR machine for  01:00:00 at  50 °C

1h



Plasmid preparation

1h

- 8 After incubation is complete, transform into competent cells (DH5α) and plate onto LB agar with appropriate antibiotic. Let this grow at  37 °C  Overnight .
- 9 Pick a single colony to grow in 5 ml LB culture at  37 °C  Overnight and purify plasmid to check sequencing.