

Aug 27, 2019

## Gibson Protocol

 In 1 collection

DOI

[dx.doi.org/10.17504/protocols.io.6t5heq6](https://dx.doi.org/10.17504/protocols.io.6t5heq6)

N.J. Hillson<sup>1,2</sup>

<sup>1</sup>JBEI; <sup>2</sup>LBNL



Michael J Fero

---

OPEN  ACCESS



**DOI:** [dx.doi.org/10.17504/protocols.io.6t5heq6](https://dx.doi.org/10.17504/protocols.io.6t5heq6)

**Protocol Citation:** N.J. Hillson 2019. Gibson Protocol. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6t5heq6>

**Manuscript citation:**

D. G. Gibson, "Enzymatic assembly of DNA molecules up to several hundred kilobases," Nature Methods, vol. 6, no. 5, pp. 343-345, 2009.

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use these protocols in our group and they work.**

**Created:** August 26, 2019

**Last Modified:** August 27, 2019

**Protocol Integer ID:** 27229

1 Prepare 6 ml of 5X ISO Buffer in a 15 ml falcon tube as follows:

3 ml 1 M Tris-HCl pH 7.5  
+ 150 ul 2 M MgCl<sub>2</sub>  
+ 240 ul 100 mM dNTP mix (25 mM each: dGTP, dCTP, dATP, dTTP)  
+ 300 ul 1 M DTT  
+ 1.5 g PEG-8000  
+ 300 ul 100 mM NAD  
+ \_\_\_\_\_dH<sub>2</sub>O to  
6 ml  
Store at -20 C in 320 ml aliquots.

2 Prepare 1.2 ml of Gibson assembly master mix as follows:

320 ul 5X ISO Buffer  
+ 0.64 ul 10 u/ml T5 exonuclease\*  
+ 20 ul 2 u/ml Phusion polymerase  
+ 160 ul 40 u/ml Taq ligase  
+ \_\_\_\_\_dH<sub>2</sub>O to  
1.2 ml  
Store at -20 C in 15 ml aliquots.

\*This is optimized for 20-150 bp sequence homology overlaps

3 Thaw a 15 ml aliquot of the Gibson assembly master mix and keep on ice until use.

4 Measure the DNA concentration (ng/ml) of each assembly piece.

5 Add 100 ng of the linearized vector backbone and equimolar amounts of the other assembly pieces to the thawed 15 ml master mix in a 20 ml total volume assembly reaction mixture as follows:

linearized vector backbone (100 ng)  
+ each additional assembly piece (to equimolar with backbone)  
+ 15 ml Gibson assembly master mix  
+ \_\_\_\_\_dH<sub>2</sub>O to  
20 ml

6 Incubate the assembly reaction at 50 C for 60 minutes, and then place on ice.



- 7 Transform 5 ml of the assembly reaction into 100 ml of competent *E. coli* and/or run a diagnostic agarose gel to check for successful assembly.