


May 15, 2019 Version 2

gmcSCRB-seq protocol V.2

DOI

dx.doi.org/10.17504/protocols.io.2aegabe



Johannes JWB Bagnoli¹, Christoph Ziegenhain², Aleksandar Janjic¹, Lucas Esteban Wange¹, Beate Vieth¹, Swati Parekh³, Johanna Geuder¹, Ines Hellmann¹, Wolfgang Enard¹

¹Ludwig-Maximilians-Universität München; ²Karolinska Institute Stockholm;

³Max Planck Institute for Biology of Ageing



Aleksandar Janjic

Ludwig-Maximilians-Universität München

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.2aegabe

Protocol Citation: Johannes JWB Bagnoli, Christoph Ziegenhain, Aleksandar Janjic, Lucas Esteban Wange, Beate Vieth, Swati Parekh, Johanna Geuder, Ines Hellmann, Wolfgang Enard 2019. gmcSCRB-seq protocol. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.2aegabe>

Manuscript citation:

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 24, 2019

Last Modified: May 15, 2019

Protocol Integer ID: 22566

Keywords: molecular crowding, scRNA-seq, SCRB-seq

Abstract

gmcSCRB-seq is an alternative lysis protocol to mcSCRB-seq (**Publication** and **Protocol**). gmcSCRB-seq uses Guanidine Hydrochloride in the lysis buffer and requires an additional clean up step. We recommend using this alternative lysis protocol in cases where cells are difficult to lyse or where the DNA appears degraded following the pre-amplification step. However, we have observed reduced sensitivity with gmcSCRB-seq compared to mcSCRB-seq and would therefore only recommend using it in cases where the latter does not produce high quality libraries.

Attachments



[mcSCRBseq_oligodT.tx..](#)

⌵

44KB

Guidelines

- For troubleshooting help, feel free to join our [mcSCRB-seq Slack channel](#), leave your question in the comments section, or message us directly.
- The complete list of reagents and plastic ware with order numbers can be found in the 'Materials' section.
- Make sure all steps involving single-cell lysate and RNA before reverse transcription are carried out swiftly.
- Size selection of libraries is optional, but has in our experience improved reliability in cluster densities when sequencing.
- All primer sequences are listed below:

Oligo	Vendor	Purification	Concentration	Sequence
barcoded oligo-dT (E3V6NEXT)	IDT	TruGrade	2 µM	Biotin- ACACTCTTTCCCTACACGA CGCTCTTCCGATCT[BC6] [UMI10][T30]VN
TSO unblocked (E5V6NEXT)	IDT	HPLC	100 µM	ACACTCTTTCCCTACACGA CGCrGrGrG
PreAmp (SINGV6)	IDT	Desalted	10 µM	Biotin- ACACTCTTTCCCTACACGA CGC
3' enrichment primer (P5NEXTPT5)	IDT	HPLC	5 µM	AATGATACGGCGACCACCG AGATCTACACTCTTTCCCTA CACGACGCTCTTCCG*A*T *C*T
i7 Index Primer (N7XX)	IDT	TruGrade	5 µM	CAAGCAGAAGACGGCATAC GAGAT[i7]GTCTCGTGGGCT CGG

Find the cell barcode sequences in the attached text file (Abstract tab).

Materials

	NAM E	CATA LOG #	VEN DOR
	EDTA 0.5M	E788 9	Sigm a Aldric h
	Ethan ol, absol ute	9065. 4	Carl Roth
	Exon uclea se I (20 U/μl)	EN05 82	Ther mo Fishe r Scien tific
	Exon uclea se I React ion Buffe r (10x)	EN05 82	Ther mo Fishe r Scien tific
	IGEP AL CA- 630	I8896	Sigm a Aldric h
	Maxi ma H- Rever se Trans cripta se (200 U/μl)	EP07 53	Ther mo Fishe r Scien tific
	Maxi ma RT Buffe r (5x)	EP07 53	Ther mo Fishe r Scien tific
	Polye thyle ne glycol 8000	8951 0	Sigm a Aldric h
	PBS 7.4	10010 -23	Gibco - Ther

			mo Fisch er
Phusi on HF Buffe r	B051 8		New Engla nd Biola bs
Protei nase K	9034		Takar a
Sera- Mag Spee d Bead s	65152 10505 0250		Ther mo Fishe r Scien tific
Terra PCR Direct Poly mera se Mix	6392 71		Takar a
Ultra Pure DNase/ RNase- Free Distill ed Water	10977 -049		Invitr ogen - Ther mo Fishe r
Sodiu m Azide 99.5 %	S200 2- 100G		Sigm a Aldric h
Sodiu m Chlori de 5M	S515 0-1L		Sigm a Aldric h
Trizm a hydro chlori de soluti on 1M pH 8.0	T269 4		Sigm a Aldric h
Bioan alyze r High Sensi tivity DNA	5067 -462 6		Agile nt Tech nolog ies



	Analy sis Kits		
	MinEl ute Gel Extra ction Kit	2860 6	Qiage n
	Nexte ra XT DNA Librar y Prepa ration Kit	FC- 131- 1096	illumi na
	Quant -iT PicoG reen dsDN A Assay Kit	P758 9	Invitr ogen - Ther mo Fishe r
	dNTP s (25 mM each)	R018 2	Ther mo Fishe r Scien tific
	Alumi nium seals for cold stora ge	391- 1275	
	Adhe sive PCR film seals	AB05 58	Ther mo Fishe r Scien tific
	twin.t ec 96- well DNA LoBin d Plate s	0030 1295 04	Eppe ndorf
	twin.t ec 384- well DNA	0030 1295 47	Eppe ndorf



LoBin d Plate s		
0.5 ml PCRcl ean tube DNA LoBin d	0030 1080 35	Eppe ndorf
1.5 ml PCRcl ean tube DNA LoBin d	0030 1080 51	Eppe ndorf
5.0 ml PCRcl ean tube DNA LoBin d	0030 10831 0	Eppe ndorf
15 ml PCRcl ean tube DNA LoBin d	0030 1222 08	Eppe ndorf
E-Gel EX Agaro se Gels, 2%	G402 002	Invitr ogen - Ther mo Fishe r

Before start

Wipe bench surfaces with RNase Away and keep working environment clean.



Preparation of lysis plates

- 1 Prepare **Lysis Buffer** for the number of plates needed.

Note

Note that 96-well plates are the preferred setup for gmcSCRB-seq, as the well cleanup (Step 9-15) would be more difficult with 384-well plates.

Reagent	96-well plate
8M Guanadine Hydrochloride	343.75
2-mercaptoethanol	5.5
NEB HF Phusion buffer (5x)	1.1
H2O	199.65
Total	550

- 2 Prepare 96 well plate(s) containing 5 μ L **Lysis Buffer** per well.

Note

Lysis plates should be prepared shortly before use, but can be stored for up to 1 week at room temperature. Prior to use, double check to make sure salts have not fallen out of solution or that the lysis buffer has not evaporated.

 5 μ L Lysis Buffer

Sample Collection

- 3 Sort 1 cell to each well of a 96 well plate containing 5 μ L **Lysis Buffer**.
- 4 Immediately seal the plate with an aluminium cold storage seal.



- 5 In a cooled centrifuge, spin down the plate for 30 sec @ max. speed and place immediately on dry ice.

🧊 -80 °C Store plates containing single-cell lysates in a -80 °C freezer for up to 6 months.

Initial Clean-up

- 6 Thaw plates briefly (up to 1 min) at room temperature
- 7 Spin down (30 sec @ 1000 rcf) in a centrifuge pre-cooled to 4 °C.

🧊 4 °C

- 8 Prepare **Clean-up Beads:**

Reagent	Amount
PEG 8000	11 g
NaCl, 5M	10 mL
Tris-HCL, 1M, pH 8.0	500 µL
EDTA, 0.5M	100 µL
IGEPAL, 10% solution	50 µL
Sodium Azide, 10% solution	250 µL
UltraPure Water	up to 49 mL
Total	49 mL

-Add all ingredients into a 50 mL falcon tube, but do not add the total amount of water until after PEG is completely solubilized

-Incubate at 40°C and vortex regularly until PEG is completely dissolved

- Resuspend bead stock carefully (Sera-Mag Speed Beads)

- Pipette 1000 µL of bead suspension into a 1.5 mL tube

- Place on magnet stand

-Remove supernatant

-Add 1000 µL 10 mM Tris-HCl, pH 8.0, 1 mM EDTA (TE), and resuspend beads

-Place on magnet stand

-Remove supernatant

-Repeat wash one more time

-Add 900 µL 10 mM Tris-HCl, pH 8.0, 1 mM EDTA (TE), and resuspend beads




-Add to PEG solution above and mix well.

Note


Beads should be prepared ahead of time and can be stored at 4 °C or room temperature.


- 9 Add 10 uL of **Clean-up Beads** to each well (1:2 ratio).

 10 µL Clean-up Beads

- 10 Seal the plate and vortex to mix the bead-lysate mixture. Briefly spin down the plate (<200 rcf, <10 sec) so all of the mixture is at the bottom of the well but the beads are still in solution.

- 11 Incubate for 5 minutes at RT.


 00:05:00 binding of the cDNA onto the beads


 20 °C Room temperature

- 12 Place on magnet stand until clear

- 13 Discard supernatant

- 14 Wash twice with 50 uL 80% ethanol (while on magnet) and discard supernatant

 50 µL 80% ethanol (freshly prepared)

- 15  00:05:00 air dry beads

Note

Depending on lab temperature and humidity, drying times can vary.

Reverse Transcription

- 16 Prepare **Reverse Transcription Mix** as follows:

Reagent	96-well plate
UltraPure Water	88 μ L
PEG 8000 (50 % solution)	165 μ L
Maxima RT Buffer (5x)	220 μ L
dNTPs (25 mM each)	44 μ L
TSO E5V6NEXT unblocked (100 μ M)	22 μ L
Maxima H Minus RT (200 U/ μ l)	11 μ L
Total	550 μL

Note


- RT MM can be prepared while plate is incubating or while drying on magnet.
- If ERCCs will be used, decrease the amount of H₂O and add appropriate amount of ERCCs.
- Caution: Reverse Transcription Mix with PEG needs to be mixed carefully!

 4 °C Keep Reverse Transcription Mix on ice

- 17 Once drying is complete, add 4 μ L H₂O to each well.

 4 μ L H₂O


- 18 Add 5 μ L Reverse Transcription Mix to each well.

 5 μ L Reverse Transcription Mix

Note

If a robot (eg. Formulatrix Mantis) is used, make sure to calibrate correctly to the viscous solution.

- 19 Add 1 μ L of barcoded oligo-dT primer [2 μ M] (E3V6NEXT adapter) to each well.

 1 μ L barcoded oligo-dT primer [2 μ M] (E3V6NEXT adapter)



- 20 Seal plate with a PCR seal, vortex briefly and spin down (30 sec @ 1000 rcf) in a centrifuge pre-cooled to 4 °C.
- 21 In a thermocycler with heated lid, incubate:
- 42 °C 90 min
- 8 °C ∞

cDNA Pooling & Purification


- 22 Vortex to mix the bead-cDNA mixture. Briefly spin down the plate (<200 rcf, <10 sec) so all of the mixture is at the bottom of the well but the beads are still in solution.
- 23 Prepare **Bead Binding Buffer**

Reagent	Amount
PEG 8000	11 g
NaCl, 5M	10 mL
Tris-HCL, 1M, pH 8.0	500 µL
EDTA, 0.5M	100 µL
IGEP AL, 10% solution	50 µL
Sodium Azide, 10% solution	250 µL
Ultra Pure Water	up to 50 mL
Total	50 mL



- Add all ingredients into a 50 mL falcon tube, but do not add the total amount of water until after PEG is completely solubilized
- Incubate at 40°C and vortex regularly until PEG is completely dissolved


- 24 Pool all wells (including beads) of one plate into a 2 mL tube and add 960 µL (ratio 1/1)
Bead Binding Buffer


 960 µL Bead Binding Buffer

Note

Vortex to mix the sample and bead binding buffer.

25

 00:05:00 binding of the cDNA onto the beads


 20 °C Room temperature

- 26 Place on magnet stand until clear

- 27 Discard supernatant

- 28 Wash twice with 2 mL 80% ethanol (while on magnet) and discard supernatant

 2 mL 80% ethanol (freshly prepared)

- 29  00:05:00 air dry beads

Exonuclease I Treatment

- 30 Elute cDNA in 17 µL UltraPure Water & transfer to new tube

 17 µL UltraPure Water

Note

Avoid transferring beads as these can inhibit the downstream PCR.



31 To the 17 μ L cDNA, add:

🧪 2 μ L Exonuclease I Buffer (10x)

🧪 1 μ L Exonuclease I (20 U/ μ L)

Full length cDNA amplification

32 In a thermocycler with heated lid, incubate:

🔥 37 °C 20 min (Exol digest)

🔥 80 °C 10 min (Heat inactivation)

🔥 8 °C ∞ (Store)

33 Prepare **PreAmplification Mix** as follows:

Reagent	1x
Terra direct Buffer (2x)	25 μ L
SINGV6 Primer (10 μ M)	1 μ L
Terra polymerase (1.25 U/ μ L)	1 μ L
UltraPure Water	3 μ L
Total	30 μL

34 Add 30 μ L of **PreAmplification Mix** directly to the Exonuclease I digested sample.

🧪 30 μ L PreAmplification Mix

cDNA purification & quantification


35 In a thermocycler with heated lid, incubate as follows:


Step	Temperature	Time	Cycles
Initial Denaturation	98 °C	3 min	1x
Denaturation	98 °C	15 sec	13-21x
Annealing	65 °C	30 sec	
Elongation	68 °C	4 min	
Final Elongation	72 °C	10 min	1x
Store	8 °C	∞	

Note

Cycle number highly depends on the input amount and should be optimized depending on the specific cell type used in the experiment. For ES cells, 13-15 cycles are sufficient.


- 36 Mix PreAmplification PCR with 40 μ L **Clean-up Beads** (1/0.8 ratio)


 40 μ L Clean-up Beads

- 37  00:05:00 binding of the cDNA onto the beads

- 38 Place on magnet until clear and discard supernatant

- 39 Wash twice with 150 μ L 80% ethanol (while on magnet) and discard supernatant

 150 μ L 80% ethanol (freshly prepared)

- 40  00:05:00 air dry beads

- 41 Elute cDNA in 15 μ L UltraPure Water & transfer to new tube

 15 μ L UltraPure Water


- 42 Quantify the cDNA using the Quant-iT PicoGreen dsDNA assay kit or equivalent Qubit following the manufacturer's protocol. Use 1 μ L of clean cDNA for quantification.

Expected result

cDNA concentration should be > 1 ng/ μ L, depending on cell type and cycle number

- 43 Optional: Quality check the cDNA using the Agilent 2100 Bioanalyzer with High Sensitivity DNA Analysis Kits.

Expected result




Tagmentation, Library PCR & Indexing

- 44 Prepare **Tagmentation Mix** and dispense 19 μL to a new 96-well plate.

Reagent	1x
Tagment DNA Buffer (2x)	10 μL
Amplicon Tagment Mix (Tn5)	5 μL
UltraPure Water	4 μL
Total	19 μL


🧴 19 μL Tagmentation Mix

- 45 Dilute cDNA to 0.8 ng/ μL and add 1 μL to each reaction.

 1 μL cDNA (0.8 ng/ μL)

46 In a thermocycler with heated lid, incubate as follows:

 55 °C Tagmentation

 00:10:00 Tagmentation

47 To stop the reaction, add 5 μL NT buffer to each reaction and mix by pipetting up and down.

 5 μL NT Buffer


 00:05:00 Incubation at room temperature

48 Prepare **3' Enrichment PCR Mix** as follows and add 24.5 μL to each tagmentation reaction.

Reagent	1x
NPM PCR Mix	15 μL
P5NEXTPT5 (5 μM)	0.5 μL
UltraPure Water	9 μL
Total	24.5 μL

 24.5 μL 3' Enrichment PCR Mix


49 Add 0.5 μL of i7 index primer (5 μM)





 0.5 μL i7 index primer (5 μM)

50 In a thermocycler with heated lid, incubate as follows:


Step	Temperature	Time	Cycles
Gap-fill	72 °C	3 min	1x
Initial Denaturation	95 °C	30 sec	
Denaturation	95 °C	10 sec	13x
Annealing	55 °C	30 sec	
Elongation	72 °C	1 min	
Final Elongation	72 °C	5 min	1x
Store	8 °C	∞	

51 Mix Index PCR with 50 μL **Clean-up Beads** (1/1 ratio)

 50 μL Clean-up Beads

- 52  00:05:00 binding of DNA onto the beads
- 53 Place on magnet until clear and discard supernatant
- 54 Wash twice with 150 μ L 80% ethanol (while on magnet) and discard supernatant
 150 μ L 80% ethanol (freshly prepared)
- 55  00:05:00 air dry beads
- 56 Elute cDNA in 20 μ L UltraPure Water & transfer to new tube
 20 μ L UltraPure Water

Size selection

- 57 Load complete library onto an 2% Agarose E-Gel EX and run for 10 minutes.
 00:10:00
- 58 As soon as the Gel run has finished open the Gel framing using the Gel opening tool
- 59 Excise the Library from 300bp to 900bp using a clean scalpel
- 60 Gel purify the slice using the Qiagen MinElute Kit following manufacturer's guidelines:
- Add 450 μ L Buffer QG
 - Dissolve the gel slice in QG for 10 min @ 42 °C
 - Add 150 μ L Isopropanol to the sample and mix by inverting
 - Transfer sample to spin column and centrifuge at 16 000 x g for 1 min
 - Discard flow through and add 500 μ L Buffer QG
 - Centrifuge at 16 000 x g for 1 min and discard flow through
 - Add 700 μ L Buffer PE
 - Centrifuge at 16 000 x g for 1 min and discard flow through
 - Centrifuge again at 16 000 x g for 1 min to remove residual ethanol
 - Transfer column to a new 1.5 mL microcentrifuge tube
 - Add 20 μ L H₂O to column and incubate for 1 min
 - Centrifuge at 16 000 x g for 1 min to elute and discard the spin column



Note

The Monarch DNA Gel Extraction Kit (NEB T1020L) can also be used.

Library Quantification

- 61 Quantify and quality control the library using the Agilent 2100 Bioanalyzer with High Sensitivity DNA Analysis Kits.


Load both the library and a 1/10 dilution on two different lanes of the chip.

Note

If cDNA was not quality controlled previously, it is strongly suggested to run cDNA on the same chip.

Expected result

Successful libraries will typically exceed 3-5 ng/ μ l concentration.



62 If pooling several libraries, combine equal molar amounts.

Sequencing

63 Sequence your library on any compatible Illumina sequencer. Dilute libraries to recommended molarity according to Illumina's recommendations (eg. 2 nM). Select the following paired-end read-length settings:



	Read	Cycles	Content
	Read 1	16	Cell barcode & UMI
	Index 1	8	i7 Index
	Index 2	0	
	Read 2	50	cDNA fragment

Primary data processing using zUMIs

- 64 Download and install zUMIs including all dependencies.

Software

zUMIs

NAME

Linux

OS

<https://github.com/sdparekh/zUMIs>

SOURCE LINK

- 65 Copy the sequencing data from the sequencer and run bcl2fastq without demultiplexing.

Command

```
bcl2fastq --use-bases-mask Y16,I8,Y50 --create-fastq-for-index-reads
```

- 66 Run zUMIs by specifying your parameters in the associated yaml file or using the **YAML config Rshiny application**. Below is an example yaml file for a typical gmcSCRB-seq run.

Command

Example zUMIs yamI

```
project: Example_project_name
sequence_files:
  file1:
    name: /data/R2.fq.gz
    base_definition: BC(1-8)
  file2:
    name: /data/R1.fq.gz
    base_definition:
      - BC(1-6)
      - UMI(7-16)
  file3:
    name: /data/R3.fq.gz
    base_definition: cDNA(1-50)
reference:
  STAR_index: /data/STAR5idx
  GTF_file: /data/Species.gtf
  additional_STAR_params: ''
  additional_files: /data/ERCC92.fa
out_dir: /data/zUMIs
num_threads: 15
mem_limit: 0
filter_cutoffs:
  BC_filter:
    num_bases: 1
    phred: 20
  UMI_filter:
    num_bases: 1
    phred: 20
barcodes:
  barcode_num: ~
  barcode_file: /data/bc.txt
  automatic: no
  BarcodeBinning: 0
  nReadsperCell: 100
counting_opts:
  introns: yes
  downsampling: '0'
  strand: 0
  Ham_Dist: 0
  velocvto: no
```



```
primaryHit: yes
twoPass: yes
make_stats: yes
which_Stage: Filtering
Rscript_exec: Rscript
STAR_exec: STAR
pigz_exec: pigz
samtools_exec: samtools
zUMIs_directory: /data/zUMIs2/zUMIs
read_layout: SE
```