

Apr 01, 2022

🌐 Golden Gate Assembly

DOI

<https://dx.doi.org/10.17504/protocols.io.x54v9yr3mg3e/v1>

Isaac Núñez¹, Tamara Matute¹, Fernan Federici¹


¹SynBioUC

Laboratorio de Tecnolog...



Isaac Núñez

Laboratorio de Tecnologías Libres, iBio



Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.x54v9yr3mg3e/v1>

Protocol Citation: Isaac Núñez, Tamara Matute, Fernan Federici 2022. Golden Gate Assembly. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.x54v9yr3mg3e/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 01, 2022

Last Modified: April 01, 2022

Protocol Integer ID: 60170

Keywords: DNA Assembly, Golden Gate, DNA fabrication, Assembly, Restriction Enzymes, TypeIIIS, Synthetic Biology, assembly of genetic sequence, assembly position, cyclic assembly, such as uloop assembly, assembly level, position of assembly, ordered assembly position, dna piece, golden gate technique, assembly, substrates of the next level reaction, final assembly product, uloop assembly, genetic sequence, cleavage of the component, typeiis restriction enzyme, restriction enzyme, enzyme recognition site, standardized dna component, competent bacterial cell, asymmetric dna sequence, enzyme, dna, synthetic syntaxes of these overhang, ligase enzyme, assembled vector, catalysis, golden gate, defining synthetic syntax, next level reaction, substrate

Abstract

The Golden Gate technique allows the assembly of genetic sequences from libraries of standardized basic components, which are cleaved from their donor vectors and concatenated in the acceptor vector in a defined order. This reaction is subsequently transformed into competent bacterial cells that are grown overnight to reveal positive colonies carrying the correctly assembled vector.

The cleavage of the components is done by a TypeIIS restriction enzyme, which recognizes asymmetric DNA sequences and cleaves outside of this site, leaving overhangs in a directional way (Figure 1A). This property allows the creation of ordered assembly positions by defining synthetic syntaxes of these overhangs (Figure 1B).




Figure 1: Golden Gate Assembly. **A:** Cleavage of TypeIIS restriction enzyme with recognition site indicated in red and overhangs indicated in gray. **B:** Ordered assembly of DNA pieces into an acceptor vector.

Within the different methodologies and libraries based on Golden Gate, some methods, such as uLoop Assembly (Pollak et al, 2020), allow the cyclic assembly of increasingly large and complex vectors in which the products of a reaction are the substrates of the next level reactions (**Figure 2**). By this means it is possible to build libraries of easily reusable components to perform assemblies in a combinatorial way of odd and even assembly levels.




Figure 2: Loop Assembly of complex vectors by cycling assembly of combinatorial components. Odd acceptor vectors are indicated in hot colors and Even vectors in cold colors.

The standardized DNA components have to be flanked by the enzyme recognition sites and the proper syntax according to the desired position of assembly. These DNA pieces can be obtained from previously created libraries or custom made.

The reaction operates cycling between 37°C and 16°C. At the first temperature, the restriction enzyme cuts the DNA pieces. These pieces are concatenated by their overhangs sequence homology and sealed by the catalysis of a ligase enzyme at the second temperature. Once a piece is correctly assembled it cannot be cutted again, and cycle by cycle, the pieces are concatenated in the final assembly product.



Materials

Reagents

Molecular grade H₂O

T4 DNA Ligase [400 U/μL] (NEB)

BsaI-HF®v2 [20 U/μL] (NEB)

SapI [10 U/μL] (NEB)

T4 DNA Ligase Buffer 10X (NEB) -- It comes with T4 DNA Ligase

CutSmart Buffer 10X (NEB) -- It comes with SapI

Purified BSA 20 mg/mL (NEB)

X-Gal 20 mg/mL in DMSO

Note

T4 ligase buffer has to be aliquoted upon arrival to avoid degradation by thaw and dethaw.  10 μL aliquots are recommended.

Materials

0.2 μL tubes

Pipette tips

Equipment & Tools

Thermocycler

P2 Pipette

P10 Pipette

Ice Bucket

Preparation of the DNA components

14h 20m

- 1 Prepare plasmids stock solutions of the desired components to be used.

14h

Note

This step is done by overnight growing of the strains in LB supplemented with proper antibiotics and performing the purification by any standard miniprep protocol. We use

 Wizard® Plus SV Minipreps DNA Purification System
Promega Catalog #A1460

- 2 Measure the concentration of purified plasmids.

5m

Note

Typically concentrations range from [M] 50-800 ng/μL depending on the plasmid.

- 3 Perform dilutions of the plasmids to working concentrations.

15m

It is:

- [M] 15 fmol/μL for donor vectors
- [M] 7.5 fmol/μL for receiver vectors

Molar concentration can be computed from mass concentration accord the next equation:

$$X \left[\frac{\text{fmol}}{\mu\text{L}} \right] = \frac{\text{Concentration} \left[\frac{\text{ng}}{\mu\text{L}} \right] * 10^6 \left[\frac{\text{fg}}{\text{ng}} \right]}{650 \left[\frac{\text{fg/fmol}}{\text{bp}} \right] * \text{length} [\text{bp}]}$$

Then, dilutions should be made in molecular grade H₂O using the standard equation:

$$v_x = \frac{c_f \left[\frac{\text{fmol}}{\mu\text{L}} \right] * v_f [\mu\text{L}]}{\left[X \frac{\text{fmol}}{\mu\text{L}} \right]}$$

**Note**

c_f : 7.5 or 15 fmol/ μ L

v_f : Final dilution volume. Typically 50 μ L to create a stock of the component

v_x : Volumen of the component at X concentration to be taken

Golden Gate Reaction

3h 21m

- 4 Prepare the **DNA mix**. In a pcr tube, mix  1 μ L of each component to be assembled (donor vectors) and  1 μ L of the acceptor vector.

5m

Note

The total volume of the DNA mix should be  5 μ L considering four components and one acceptor, however the number of components may be different. If this is the case, the volume difference must be corrected by adjusting the volume of water in the 2X Reaction Master Mix (next step).

- 5 Prepare the **2X Reaction Master Mix** (Bsal or SapI)  On ice according to the next table:

5m

Note

To avoid too small volume pipetting and variations, we recommend to do a master mix adding [1.1 * each component] and then divide into each tube.

Note

Bsal mix is for *odd level* assemblies and SapI mix for *even level* assemblies.

STEP CASE**2X Reaction Master Mix Bsal**

5 steps



Component	Volume (μL)
H ₂ O (HPLC grade)	3,00
T4 DNA Ligase Buffer 10X (NEB)	1,00
BSA [1mg/mL]	0,50
T4 DNA Ligase [400 U/μL] (NEB)	0,25
BsaI-HF®v2 [20 U/μL] (NEB)	0,25

*volumes per reaction

- 6 Mix the  5 μL of the DNA mix (step 4) with  5 μL of the 2X Reaction Master Mix (step 5). 1m

Note

Spin the tube to ensure the whole reaction is at the bottom of the tube and free of bubbles.

- 7 Bring the tubes to a thermocycler and run it according the next program: 3h 10m

Cycling parameters	
3 minutes at 37°C	Repeat 25X
4 minutes at 16°C	
5 minutes at 50°C	
10 minutes at 80°C	

Note


You can short-store the reactions at  4 °C or  -20 °C for long storage.

Transformation and colony selection 14h 33m

- 8 Directly transform  5 μL of the reaction into  50 μL chemo-competent cells. 2h 30m

- 9 Plate the transformant cells into plates supplemented with the proper antibiotic plus X-Gal **100** 40 µg/mL and growth overnight at **37** 37 °C 12h
- 10 Identify the colonies carrying the proper assembled vectors based on blue - white screening of the grown colonies (blue colonies carry intact acceptor vector and should not be selected), plus any other particular screening criteria related to your assembly (i.e. fluorescence expression, colony PCR, sequencing, etc). 3m

Expected result



Golden gate transformation results. LacZ(+) colonies become "blue" by reacting with X-Gal and white colonies should carry properly assembled vectors.