

Sep 23, 2019

Golden Gate Assembly (Protocol for NEB® Golden Gate Assembly Mix)

DOI

dx.doi.org/10.17504/protocols.io.7kfhktn

Alba Balletbó¹

¹Wageningen University

iGEM Wageningen 2019



Alba Balletbó

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.7kfhktn>

Protocol Citation: Alba Balletbó 2019. Golden Gate Assembly (Protocol for NEB® Golden Gate Assembly Mix). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.7kfhktn>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: September 23, 2019

Last Modified: September 23, 2019

Protocol Integer ID: 28007

Keywords: golden gate assembly, type iis recognition site distal, recognition sequence from the assembly, such that the type iis rease, type iis rease, ability of type iis restriction, type iis restriction, dna, assembly

Abstract

Golden Gate Assembly and its derivative methods exploit the ability of Type IIS restriction endonucleases (REases) to cleave DNA outside of the recognition sequence. The inserts and cloning vectors are designed to place the Type IIS recognition site distal to the cleavage site, such that the Type IIS REase can remove the recognition sequence from the assembly.

Guidelines

For complex (>10 fragment) assemblies, high efficiencies are achievable with increased ligase and BsaI-HFv2 levels (1000 units T4 DNA Ligase, 30 units BsaI-HFv2), as listed in this protocol. For assemblies involving 10 fragments and less, the standard amounts (500 units T4 DNA Ligase, 15 units BsaI-HFv2) are sufficient. Note the reaction volume of 25 µl is used to allow sufficient volume for precloned insert additions, if needed.

Materials

MATERIALS

☒ NEB 10-beta Competent E.coli (High Efficiency) - 6×0.2 ml **New England Biolabs Catalog #C3019I**

☒ NEB Golden Gate Assembly Mix - 15 rxns **New England Biolabs Catalog #E1600S**

☒ T4 DNA Ligase **New England Biolabs Catalog #M0202**

☒ NEB 10-beta/Stable Outgrowth Medium - 4×25 ml **New England Biolabs Catalog #B9035S**

☒ BsaI-HFv2 **New England Biolabs Catalog # R3733L**

pGGA Destination Plasmid*

LB Agar plates with chloramphenicol

* Included in the NEB Golden Gate Assembly Mix (**NEB #E1600**)

Troubleshooting

Safety warnings

- ! Wear laboratory coat, gloves and goggles. Always check the safety warnings indicated by the Mix supplier.



Assembly Reactions

1 Set up 25 μ l assembly reactions as follows:

		Assembly Reaction	Negative Control (If desired)
	pGG A Destination Plasmid*, 75 ng/ μ l	1 μ l (75 ng)	1 μ l (75 ng)
	24 precloned inserts cloned into pMin iT 2.0, 100 ng/ μ l each plasmid	0.75 μ l (75 ng) each, (18 μ l total)	-
	T4 DNA Ligase Buffer (10X)	2.5 μ l	2.5 μ l
	T4 DNA Ligase, 2000 U/ μ l	0.5 μ l (1000 units)	0.5 μ l (1000 units)



Bsal - HFv 2, 20 U/ μ l	1.5 μ l (30 units)	1.5 μ l (30 units)
Nucl ease -free H ₂ O	1.5 μ l	19.5 μ l


*or user provided

- 2 Mix gently by pipetting up and down 4 times.
- 3 Briefly microcentrifuge (1 second) to bring material to the bottom of tube.
- 4 Transfer to thermocycler and program as follows: (5 min 37°C → 5 min 16°C) x 30 cycles followed by 5 minutes 60°C. If reactions are done overnight, add a 4°C terminal hold to the protocol, but repeat the final 5 minutes 60°C step the next day before the transformations.


Transformation

- 5 For each assembly, thaw a 50 μ L tube of NEB 10-beta competent *E. coli* cells on ice for 5–10 minutes.
- 6 ▮ Add 2 μ L of the assembly reaction; gently mix by flicking the tube 4–5 times.
- 7 ▮ Incubate on ice for 30 minutes.
- 8 ▮ Heat shock at 42°C for 30 seconds.
- 9 ▮ Place back on ice for 5 minutes.



- 10  Add 950 μL of room temperature NEB 10-beta/Stable Outgrowth Medium. Incubate at 37°C for 60 minutes, shaking vigorously (250 rpm) or using a rotation device.

Plating

- 11 Warm LB agar plates containing chloramphenicol (for pGGA) at 37°C for 15 minutes.
- 12  Mix the cells thoroughly by flicking the tube and inverting, then spread 100 μL outgrowth onto each plate.
- 13 Incubate the plates overnight at 37°C, or 24 hours at 30°C, or 48 hours at 25°C.