



Jan 16, 2026

Golgi-Cox stain and neuronal morphological analysis

 [PLOS One](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.q26g7n643lwz/v1>

Mario Alberto Bautista-Carro¹

¹Cinvestav



Mario Alberto Bautista-Carro

Cinvestav

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q26g7n643lwz/v1>

External link: <https://doi.org/10.1371/journal.pone.0339028>

Protocol Citation: Mario Alberto Bautista-Carro 2026. Golgi-Cox stain and neuronal morphological analysis. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.q26g7n643lwz/v1>

Manuscript citation:

Bautista-Carro MA, Sánchez-Teoyotl P, Juárez-Serrano D, Iannitti T, Díaz A, Flores G, Morales-Medina JC (2026) Olfactory bulbectomy induces neurobiological alterations in the prefrontal cortex and hyperlocomotion in male rats. PLOS One 21(1). doi: [10.1371/journal.pone.0339028](https://doi.org/10.1371/journal.pone.0339028)



License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 08, 2025

Last Modified: January 16, 2026

Protocol Integer ID: 229355

Keywords: Golgi-Cox method, Neuronal morphology, Neuron reconstruction, Dendritic spines, Morphometric analysis, Sholl analysis, neuronal morphological analysis the golgi, cox stain on rat brain, detailed study of neuronal morphology, neuronal morphological analysis, neuronal morphology, neuronal organization, golgi, cox stain, performing golgi, help of optical microscopy, optical microscopy, rat brain, dendritic arborization, spine density

Abstract

The Golgi-Cox technique allows for detailed study of neuronal morphology. With the help of optical microscopy and Sholl analysis, dendritic arborization, spine density, and neuronal organization can be analyzed. This allows us to identify structural changes in physiological or pathological conditions. Below we describe protocol for performing Golgi-Cox stain on rat brains and Sholl analysis.

Guidelines

The entire tissue preparation process must be carried out in a darkroom. Be sure to work in a laboratory fume hood when necessary.

Materials

Perfusion pump, Ketamine, Xylazine, Sterile saline solution, 0.9% saline solution, Syringes, Falcon tubes, Golgi-Cox solution, Vibratome, Gelatinized slides, Ammonium hydroxide, Kodak film fixer, Alcohol, Xylene, Light microscope equipped with camera lucida, Permount mounting medium, Cover slips.

Before start

Ensure you have all the necessary laboratory materials and equipment.



Perfusion

- 1 Anesthetize the rats by intraperitoneal injection of a ketamine/xylazine cocktail (0.75 ml ketamine + 0.25 ml xylazine + 5 ml sterile saline), administering a dose of 0.125 ml per 20 g of body weight.
- 2 Perfuse the rats intracardially with 0.9% saline solution and quickly remove the brain.

Tissue preparation

4w 5d

- 3 Place the brains in Golgi-Cox solution for 30 days.
- 4 Incubate 30% sucrose solution for 3 days.
- 5 Using a vibratome (Leica, VT1000S microsystem, USA), obtain 200 μ m coronal sections of the brain region of interest and mount them on 2% gelatinized slides. Keep the slides in a humid chamber.

4w 2d



3d

Staining procedure (in darkroom)

1h 45m 25s

- 6 Rinse with distilled water 5 times.
- 7 Treat sections with ammonium hydroxide for 30 minutes.
- 8 Rinse with distilled water 10 times.
- 9 Incubate for 30 minutes in Kodak Film Fixer.
- 10 Rinse with distilled water 10 times and dehydrate in increasing concentrations of alcohol: 75% 5 minutes, 90% 5 minutes, and 100% 10 minutes twice.

5s

30m



10s

30m

30m 10s



11 Incubate with xylene for 15 minutes.

15m

12 Cover the slides with Permount mounting medium and coverslips.

Neuronal morphological analysis

13 Using an optical microscope (DM 2000 microscope, Leica Microsystems, USA) equipped with a camera lucida, reconstruct ten pyramidal neurons (5 in each hemisphere) from each brain in a two-dimensional plane.



14 Quantify dendritic tracing by Sholl analysis. Place a transparent grid with concentric rings spaced 10 μm apart over the dendritic drawing (centering the soma) and count the number of intersections. Calculate the total dendritic length by multiplying the total number of intersections in each ring by 10 μm . Calculate the dendritic arborization as the total number of dendritic branches present in each order of the dendritic axis.



15 To obtain the density of spines, count and draw the spines of ten distal and/or last-order dendrites with an approximate length of 10 μm (5 for each hemisphere) at 1000x magnification.



Protocol references

Flores, G., Morales-Medina, J. C., Rodríguez-Sosa, L., 26 Calderón-Rosete, G. (2015). Characterization of Cytoarchitecture of dendrites and fiber neurons using the golgi-cox method: An overview. *Horizons in neuroscience research*, 137-146.

Espinoza, I., Gómez-Villalobos, M. J., Aguilar-Hernández, L., Flores, G., 26 Morales-Medina, J. C. (2025). Cerebrolysin treatment improved short-term memory deficits while simultaneously increasing hippocampal spine density in hypertensive female rats. *Behavioural brain research*, 481, 115436.

Sholl D. A. (1953). Dendritic organization in the neurons of the visual and motor cortices of the cat. *Journal of anatomy*, 87(4), 387-406.