

May 24, 2025

Good practices in hPSC handling

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk287wl5r/v1>



Valeria Fernandez Vallone¹, Katarzyna Ludwik¹, Lyn Healy², Nathalie Lefort³, Tamer Onder⁴, Lisa Pavinato^{5,6}, Fatma Visal OKUR⁷, Harald Stachelscheid¹

¹Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany;

²The Francis Crick Institute;

³Université de Paris, Imagine Institute, iPSC Core Facility, INSERM UMR U1163, F-75015 Paris, France.;

⁴Koc University;

⁵Institute of Oncology Research (IOR), Bellinzona Institutes of Science (BIOS+), Bellinzona, Switzerland;

⁶Faculty of Biomedical Sciences, Università della Svizzera Italiana, Lugano, Switzerland.;

⁷Hacettepe University, Center for Stem Cell Research and Development (PEDISTEM) and Hacettepe University Faculty of Medicine, Department of Pediatrics

Valeria Fernandez Vallone: Contributed equally;

Katarzyna Ludwik: Contributed equally;

CorEuStem



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Inst...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk287wl5r/v1>

Collection Citation: Valeria Fernandez Vallone, Katarzyna Ludwik, Lyn Healy, Nathalie Lefort, Tamer Onder, Lisa Pavinato, Fatma Visal OKUR, Harald Stachelscheid 2025. Good practices in hPSC handling. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.j8nlk287wl5r/v1>

License: This is an open access collection distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this collection and it's working

Created: April 30, 2025

Last Modified: May 24, 2025

Collection Integer ID: 198125

Keywords: reproducibility, hPSC, pluripotent stem cells, hiPSC, hESC, standards, culture, E8 medium, passaging, cryopreservation, maintenance, survival factors, ROCK inhibitor, CEPT, CloneR2, single cells, reference images, matrix, manual picking, colony, mTeSR, E8 preparation, thawing, healthy hpsc morphology under various culture condition, standardization of stem cell culture technique, human pluripotent stem cell, stem cell culture technique, pluripotent stem cell, illustrating healthy hpsc morphology, stem cell, cell quality assessment, critical visual benchmarks for cell quality assessment, enhanced cell viability, extensive experience from established hpsc core facility, established hpsc core facility, cell viability, hpsc, cell viability during critical manipulation, foundational resource for laboratory, including rock inhibitor, cell

Funders Acknowledgements:

COST CorEuStem

Grant ID: CA20140



Abstract

This protocol collection provides a comprehensive and practical guide for the maintenance, passaging, freezing, and thawing of human pluripotent stem cells (hPSCs). Designed to support reproducibility and transferability, especially among new users, the protocols are informed by extensive experience from established hPSC core facilities.

In addition to detailed step-by-step procedures, the collection includes high-quality reference images illustrating healthy hPSC morphology under various culture conditions, as well as visual examples of spontaneous differentiation—serving as critical visual benchmarks for cell quality assessment. The protocols also cover the in-house preparation of Essential 8 (E8) medium and the use of key survival-supporting compounds, including ROCK inhibitors (ROCKi), CEPT cocktail, and CloneR2, ensuring enhanced cell viability during critical manipulations. By consolidating expert-validated practices and making them accessible through a user-friendly format, this collection aims to serve as a foundational resource for laboratories working with hPSCs and to facilitate the standardization of stem cell culture techniques across different research environments.

Guidelines

- Always work aseptically in a Class II biosafety cabinet.
- Ensure that you keep the safety cabinet clean and tidy.
- For incubation, place culture vessels on trays to contain any accidental spillage of media and reduce the risk of contamination.

To avoid cross contamination between cell lines:

- Only one cell line at a time should be handled in a biosafety cabinet when feeding or passaging the cells.
- Bottles or aliquots of media should be prepared separately for each cell line.
- Aliquots should be labelled with the name of the line and the date prepared and stored and handled appropriately.

This protocol collection describes the culture of hPSC under normoxia conditions (5%CO₂) due to availability in most cell culture laboratories. However, to reduce genomic instability during culture we recommend to incubate hPSC under **hypoxia conditions**: 5%CO₂ and 5% O₂.

Safety warnings

- ! The estimated time for an individual protocol does not account for the time required to complete supporting protocols.

Before start

Before starting this protocol, ensure that you have all of the materials required to complete the protocol with ease and efficiency.

Files

 SEARCH

Protocol

NAME



Maintenance of hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →

Protocol

NAME



Reference pictures of hPSC cultured in defined conditions

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →

Protocol

NAME



E8 media production

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →

Protocol

NAME



Survival factors for hPSC growth

VERSION 1



CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

[OPEN](#) →

Protocol

NAME



Coating of tissue Culture Vessels for hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

[OPEN](#) →

Protocol

NAME



Thawing of cryopreserved hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

[OPEN](#) →

Protocol

NAME



Non-enzymatic passaging of hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

[OPEN](#) →

Protocol

NAME



Single cell passaging of hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →

Protocol

NAME



Manual picking and passaging of hPSC colonies

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →

Protocol

NAME



Freezing of hPSC

VERSION 1

CREATED BY



Harald Stachelscheid

Core Unit pluripotent Stem Cells and Organoids - Berlin Institute of Health @ Charite, Berlin, Germany

OPEN →