


Dec 04, 2019

Version 1

# gRNA design and cloning with BbsI into Loop plasmid L1\_lacZgRNA-Ck2/3 V.1

DOI

[dx.doi.org/10.17504/protocols.io.94ah8se](https://doi.org/10.17504/protocols.io.94ah8se)



Eftychis Frangedakis<sup>1</sup>, Marta Marta Tomaselli<sup>2</sup>, Susana Sauret-Gueto<sup>2</sup>

<sup>1</sup>University of Cambridge; <sup>2</sup>Plant Sciences, University of Cambridge, OpenPlant

OpenPlant Project



**Susana Sauret-Gueto**

Plant Sciences, University of Cambridge, OpenPlant

## Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.94ah8se>

**Protocol Citation:** Eftychis Frangedakis, Marta Marta Tomaselli, Susana Sauret-Gueto 2019. gRNA design and cloning with BbsI into Loop plasmid L1\_lacZgRNA-Ck2/3. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.94ah8se>



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** December 04, 2019

**Last Modified:** December 04, 2019

**Protocol Integer ID:** 30562

**Keywords:** l1-grna transcription unit, cas9 transcription unit via l2 sapi loop assembly, grna target sequence, guide rna target sequence, plasmid l1-laczgrna, cas9 transcription unit, antibiotic resistance transcription unit, dual grna, grna design, grna, rna, cloning, l2 sapi loop assembly, bbsi into loop, mpef1 $\alpha$ , cas9


## Abstract

This protocol explains how to design and clone the guide RNA target sequence into a L1 plasmid ready to accept the gRNA by cloning with BbsI. L1 plasmids are L1\_lacZgRNA-Ck2 and L1\_lacZgRNA-Ck3.

If one gRNA target sequence is cloned into the Ck2 plasmid and another one into Ck3 one, the two L1\_gRNA transcription units can be combined with an antibiotic resistance transcription unit and a MpEF1 $\alpha$ :Cas9 transcription unit via **L2 SapI Loop assembly**. This allows for dual gRNA editing.

## Troubleshooting

1



Design of oligos for gRNA BbsI mediated cloning into L1\_lacZgRNA-Ck2 or L1\_lacZgRNA-Ck3 vectors. The L1\_lacZgRNA-Ck2/3 BbsI digested vectors have GAGC and TTTA overhangs. Therefore, oligos for gRNA should be designed such that the forward strand has a 5' overhang of CTCG and the reverse strand has a 5' overhang of gt-AAAT (addition of "gt" nucleotides is necessary to reconstitute the full sequence of the gRNA scaffold in pink). Light brown arrows: BbsI recognition site. Light brown dashed lines: BbsI cleavage site. LacZ: lacZα cassette for blue-white screening of colonies.

## 2 gRNA oligo design

Order two oligos that contain the forward and reverse guide sequence plus the overhangs necessary for ligation (highlighted with bold) into L1\_lacZgRNA-Ck2 or L1\_lacZgRNA-Ck3:

**oligo F:** 5'- CTCG-NNNNNNNNNNNNNNNNNNNNNN-gt 3'

**oligo R:** 5'- **TAAAc**-NNNNNNNNNNNNNNNNNNNNNNNN-3'

Note: Standard de-salted oligos are ok

### 3 Oligo annealing.

Mix oligos with water as follow:

oligo F (100μM) 1μl



oligo R (100 $\mu$ M) 1 $\mu$ l

water 8 $\mu$ l

Total volume 10 $\mu$ l

Anneal in a thermocycler using the following parameters: 37°C for 30 min, 95°C for 5 min and then ramp down to 25°C at 5°C per min. After annealing the gRNA can be directly cloned into L1\_lacZgRNA-Ck2 or L1\_lacZgRNA-Ck3 plasmid without the need of any further processing (step 4).

#### 4 **Cloning into backbone vector**

In a 0.2 mL tube set up the following reaction:

Component	Volume ( $\mu$ L)
Sterile water	6
10x T4 ligase buffer (NEB)	1
1 mg/mL bovine serum albumin (NEB)	0.5
T4 DNA ligase (5 U/ $\mu$ L) (Thermo Fisher)	0.25
SapI (Lgul) (5 U/ $\mu$ L) (Thermo Fisher)	0.25
L1_lacZgRNA-Ck2 or L1_lacZgRNA-Ck3 (25-50 ng)	1
annealed oligo	1
Final volume	10

- Place samples on the thermocycler and incubated using the following program:

Assembly: 15 cycles: 3 minutes at 37°C and 4 minutes at 16°C

Termination: 5 minutes at 50°C and 10 minutes at 80°C



- Transform chemically competent using 1  $\mu\text{L}$  of reaction and plate on LB agar plates with 100  $\mu\text{g/mL}$  spec and X-gal 40  $\mu\text{g/mL}$  for blue-white screening. Incubate at 37°C for 16 h.
- Confirm with sequencing