


Mar 20, 2019

Version 1

Growth Incremental analysis of thins sections from marine molluscs V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.y7efzje>



Niklas Hausmann¹, Harry K. Robson¹, Chris Hunt²

¹University of York, BioArCh; ²Liverpool John Moores University, Natural Sciences and Psychology,



Niklas Hausmann

University of York

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.y7efzje>

Protocol Citation: Niklas Hausmann, Harry K. Robson, Chris Hunt 2019. Growth Incremental analysis of thins sections from marine molluscs. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.y7efzje>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 14, 2019

Last Modified: March 20, 2019

Protocol Integer ID: 21446

Keywords: *Ostrea edulis*, Incremental growth-line analysis, Seasonality, ecological baseline, laser induced breakdown spectroscopy, geochemical data, structural context for the geochemical data, marine ecologies through time, marine ecology, growth incremental analysis of thins section, elemental map, geochemical pattern, marine mollusc, elemental maps with the result, sustainable marine resource, isotopic analysis, marine mollusc shell, overall robustness of palaeoecological interpretation, growth pattern, palaeoecological interpretation, annual growth pattern, annual growth increment, prehistoric seasonality practice, repeating growth pattern

Abstract

Annual growth patterns in marine mollusc shells are valuable indicators of the condition of marine ecologies through time. In archaeological contexts, the mollusc's time of death (i.e. the last season of growth) is an indicator of human exploitation patterns throughout the year, enabling the reconstruction of when and how often gathering occurred as well as when sites were occupied. Both pieces of information, growth rate and season of death, are vital for understanding exploitation pressure in the past, and building baselines for modern environmental policies that secure sustainable marine resources. Previously, these parameters have been determined by incremental growth-line or isotopic analyses, which are time consuming and often expensive techniques, thus restricting sample size and the overall robustness of palaeoecological interpretations.

Here, we apply Laser Induced Breakdown Spectroscopy (LIBS) to produce elemental maps (Mg/Ca) with the potential to trace and display growth patterns quickly, and at a reduced cost. We further compare the elemental maps with the results of incremental growth-line analysis to provide a structural context for the geochemical data, and demonstrate the utility of an integrated methodological approach.

Our study was undertaken on 12 European oysters (*Ostrea edulis*) from the Late Mesolithic shell midden at Conors Island, Co. Sligo in the Republic of Ireland. Our LIBS analysis enabled us to accurately and quickly determine repeating growth patterns, which were in agreement with the annual growth increments visible through the microscopic analysis. Based on this comparative dataset, including structural and geochemical patterns, the Late Mesolithic site of Conors Island had been occupied throughout the year. Moreover, our analyses highlight the applicability of LIBS to determine prehistoric seasonality practices as well as biological age and growth at an improved rate and reduced cost than was previously achievable.

Materials

STEP MATERIALS

-  Ethanol
-  Release Agent **Buehler Catalog #20-8185-016**
-  Epoxy Resin **Buehler Catalog #No 20-8140-12B**
-  Epoxy Hardener **Buehler Catalog #20-8142-016**
-  Polycrystalline suspension **Buehler Catalog #40-6016**
-  Adhesive **Catalog #n/a**

Protocol materials

-  Ethanol
-  Release Agent **Buehler Catalog #20-8185-016**
-  Epoxy Resin **Buehler Catalog #No 20-8140-12B**
-  Epoxy Hardener **Buehler Catalog #20-8142-016**
-  Polycrystalline suspension **Buehler Catalog #40-6016**
-  Adhesive **Catalog #n/a**

Producing thin sections

1

Equipment

ISOMET 1000 Precision Saw

NAME

Saw

TYPE

Buehler

BRAND

11-2180

SKU

Using a Buehler Diamond Wafering Blade (Series 15LC Diamond No. 11-4276)

SPECIFICATIONS



Section oysters following Milner (2001) and using a precision saw with a diamond blade. First section is at the hinge and in direction of growth, the second (optional) section is done to remove unnecessarily long portions of the shell section.

Milner, N., 2001. At the Cutting Edge: Using Thin Sectioning to Determine Season of Death of the European Oyster, *Ostrea edulis*. Journal of archaeological science, 28(8), pp.861–873.

2 Clean cylindric plastic moulds with

 Ethanol

leave to dry for  00:05:00 . Then coat the inside with

 Release Agent **Buehler Catalog #20-8185-016**

and dry again for  00:30:00

3

 Epoxy Resin **Buehler Catalog #No 20-8140-12B**

Place sectioned shells face down in 1- and 2-inch plastic moulds depending on the size of the cut hinge.

Then add resin (2 parts by volume of Epoxy Resin mixed with 1 part by volume Epoxy Hardener.

 Epoxy Hardener **Buehler Catalog #20-8142-016**

Leave to cure in an environment of constant humidity and temperature overnight (or

 12:00:00)

4

Equipment

Motorpol 2000

NAME

Polisher

TYPE

Buehler

BRAND

49-7268

SKU



Remove resin block from plastic moulds once hardened, and grind hinge-surface using successively finer metallographic grit papers (P600, P1200 and P2500 grades respectively) on a Buehler Motorpol 2000.

 Polycrystalline suspension **Buehler Catalog #40-6016**

Clean the surfaces of the resin blocks with water and polish further using a polishing cloth and polycrystalline suspension of 3µm.

Leave to dry overnight (or  12:00:00)

5

 Adhesive **Catalog #n/a**

Attach the polished facets of the resin blocks to a glass slide using Loctite 322 Adhesive, exposed to an ultraviolet light source and leave to harden for  24:00:00 .

6

Equipment

ISOMET 1000 Precision Saw

NAME

Saw

TYPE

Buehler

BRAND

11-2180

SKU

Using a Buehler Diamond Wafering Blade (Series 15LC Diamond No. 11-4276)

SPECIFICATIONS



Once the sample had bonded to the slide, cut the the bodies of the resin blocks from the slide using a Buehler ISOMET 1000 Precision Saw and leave a slice of the embedded hinge of approximately 200 to 500 µm in thickness on the glass slide.

7

Repeat Step 4, in order to lightly ground the shell and the resin to produce the final thin section of ~50 µm.