

Sep 10, 2021 Version 2

Growth of mixed *E. coli* colonies V.2

DOI

dx.doi.org/10.17504/protocols.io.bx3epqje

Wolfram Moebius¹

¹University of Exeter



Wolfram Moebius

University of Exeter

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.bx3epqje

External link: <https://doi.org/10.1128/mBio.01542-21>

Protocol Citation: Wolfram Moebius 2021. Growth of mixed *E. coli* colonies . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.bx3epqje> Version created by **Wolfram Moebius**

Manuscript citation:

Aranda-Díaz A, Rodrigues C, Grote A, Sun J, Schreck C, Hallatschek O, Souslov A, Möbius W, Huang KC, Bacterial Filamentation Drives Colony Chirality. *mBio* 12(6). doi: [10.1128/mBio.01542-21](https://doi.org/10.1128/mBio.01542-21)

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working.

Created: September 08, 2021

Last Modified: September 10, 2021





Protocol Integer ID: 53062

Abstract



Growth of mixed *E. coli* colonies on agar plates and between an agar plate and an agar pad



Make plates

- 1 Autoclave Lennox LB (10 g/l tryptone / casein digest peptone, 5 g/l NaCl, 5 g/l yeast extract) with 1.5 % agar (w/v).
- 2 **Optional:** Add antibiotics at desired concentration after medium has cooled down sufficiently.
- 3 Pipette  10 mL of LB and agar in Petri dishes with diameter of  6 cm .
- 4 Leave plates to dry  Overnight at  Room temperature . Optionally enclose to avoid nonuniform drying.




Store and prepare plates

- 5 If not used the day after pouring plates, store at  4 °C . Enclose in plastic container or bag to avoid further drying of plates.
- 6 Before usage of plates that have been refrigerated, warm up plates at  37 °C .

Prepare bacterial cultures

- 7 Grow overnight culture (colony picked or directly from glycerol frozen stocks) in Lennox LB.

Inoculate colonies and incubate

- 8 Mix  900 µL Lennox LB and  50 µL of each of the two overnight cultures in an Eppendorf tube.
Vortex.
- 9 Inoculate  1 µL in centre of plate and let dry.



- 10 **Optional:** Cover colonies with agar pad. Cut about x agar pad inside plate. Use spatula to lift pad. Place vertically next to colony and let fall upside down onto colony. Use spatula to remove bubbles by pressing on top.
- 11 Place plates into container, together with sufficiently wet paper towels. Incubate for 7-8 days in a dark environment.