

Dec 20, 2018

GUIDE-seq simplified library preparation protocol (CRISPR/Cas9 off-target cleavage detection)

DOI

dx.doi.org/10.17504/protocols.io.wikfccw

Nagendra Palani¹

¹University of Minnesota Genomics Center



Nagendra Palani

University of Minnesota

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.wikfccw

Protocol Citation: Nagendra Palani 2018. GUIDE-seq simplified library preparation protocol (CRISPR/Cas9 off-target cleavage detection). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.wikfccw>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We developed this protocol at the University of Minnesota Genomics Center (UMGC) and it is working. UMGC will prepare GUIDE-seq NGS libraries and sequence them as a service.

Created: December 20, 2018

Last Modified: December 20, 2018

Protocol Integer ID: 18732

Keywords: guideseq, crispr, cas9, offtarget, cleavage, NGS, library

Abstract

GUIDE-seq is an experimental method to detect off-target cleavages caused during CRISPR/Cas editing.

See publication <https://doi.org/10.1038/nbt.3117> for details about the method, the protocol used, and the analysis performed.

We identified several areas of improvement in the published protocol for GUIDE-seq NGS library preparation and developed a simplified protocol that is easier to perform & less expensive. The new protocol is also load-and-go compatible with standard sequencing runs on all Illumina instruments (using the kit SBS reagents & default instrument settings) compared to the NBT protocol that requires the libraries to be run on a MiSeq with some tweaks to the software configuration.

Notable changes:

1. Moved the UMI (Unique Molecular Identifier) from the Y-adapter to be inline with the insert. UMI is now sequenced as part of the sequencing read (R2 on a paired-end run). This eliminates the need to export index reads from the MiSeq.
2. Removed the sample index from the Y-adapter and changed the PCR scheme. Now, only a single adapter is required for ligation to any number of samples, PCR1 is for enrichment of the dsOligo - genome junction, PCR2 is to add sample-specific indexes to the enriched amplicons.

See oligo designs for further information.

3. Replaced sonication based shearing followed by end-repair with a one-pot, one-step enzymatic fragmentation + end-repair kit from New England Biolabs.

Guidelines

Time from genomic DNA sample to sequencing-ready library takes 2 days.

Protocol can be performed with 1.5 ml tubes or 96-well plates.

Materials

MATERIALS

✕ TE buffer

✕ NEBNext Ultra II Q5 Master Mix **New England Biolabs Catalog #E7649**

✕ Thermal cycler

✕ NEBNext Ultra II FS DNA Library Prep with Sample Purification Beads - 24 rxns **New England Biolabs Catalog #E6177S**

Oligonucleotides required:

TA_Adaptor_Top - /5Phos/CTCACCGCTCTTGTAGS NNNNNNNN CTGTCTCTTATACACATCTCCGAG*C

TA_Adaptor_Bottom - CTACAAGAGCGGTGAGT

dsODN_Enrich_Plus - TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG NNN
*GTTTAATTGAGTTGTCATATGTTAATAACG*G*

dsODN_Enrich_Minus - TCGTCGGCAGCGTCAGATGTGTATAAGAGACAG NNN
*CCGTTATTAACATATGACAACTCAATTAAA*C*

Replace the italicized sequences with sequences targeting the dsODN that you used in your experiments.

dsODN_Enrich_Adaptor - GTCTCGTGGGCTCGGAGATGTGTATAAGAGACA*G

Nextera Indexing Oligo design

Nextera_R1 : AATGATACGGCGACCACCGAGATCTACAC [i5] TCGTCGGCAGCGT*C

Nextera_R2 : CAAGCAGAAGACGGCATACGAGAT [i7] GTCTCGTGGGCTCG*G

You will need to replace the [i5] & [i7] with unique index barcodes. You can use the Hamming barcodes from <https://doi.org/10.1371/journal.pone.0036852>. Order the appropriate number of indexing oligos for your samples based on if you want to do combinatorial indexing (cheaper) or unique dual indexing (better with the patterned flowcells).

Note: When selecting index barcodes, make sure

- 1) there is no collision with other barcodes that will be run on the same flowcell as your samples
- 2) don't start with or have consecutive G nucleotides if the samples will be run on instruments that use two-color (NextSeq, NovaSeq) or one-color (iSeq) chemistry.



Safety warnings

! NOTE:

The NGS library generated by this protocol is NOT plug-and-play compatible with the guideseq analysis package (<https://github.com/aryeelab/guideseq>).



We have made slight modifications to the guideseq scripts to process these libraries. The output will be the same figure (of off-target sites & their frequencies) as you get with the default package. I will add a comment to this section when the modified package is made available on Github.

Before start

Make sure that the genomic DNA is clean ($260/230$ & $260/280 > 1.8$) and 400 ng of DNA is in a volume less than 26 μ l. Measure DNA concentration by fluorometry (Picogreen or Qubit).

Adaptor formation

1

- Use 0.1x TE to prepare 100 μ M stock of all the oligos.
- In a 0.2 ml PCR tube, add the following reagents, vortex well, and spin down. Place on the thermal cycler at  95 °C for  00:03:00 (heated lid).

 15 μ L Oligo TA_Adaptor_Top (100 uM)

 15 μ L Oligo TA_Adaptor_Bottom (100 uM)

 70 μ L dH₂O

 100 μ L Total




- After incubation, terminate the incubation and let the heat block cool down to room temperature. Remove the tube and place it on the bench for 1 hr.
- The oligos will have annealed and formed the adaptor.

DNA fragmentation, end-repair, ligation

2

This section of the protocol is from the manual for NEB Ultra II FS kit - Catalog # E7805.

(Refer to <https://www.neb.com/-/media/catalog/datacards-or-manuals/manuale7805.pdf> : Section for Input > 100 ng for detailed protocol)

- Aliquot  400 ng of genomic DNA into a well of a 96-well plate and use 1x TE to bring the volume to  26 μ L .
- Add the NEBNext Ultra II FS buffer + enzyme mix to the DNA (follow protocol from the linked manual)
- Set up the following program on the thermal cycler  75 °C heated lid



🌡️ 37 °C for ⌚ 00:10:00

🌡️ 65 °C for ⌚ 00:30:00

🌡️ 4 °C hold.

Note

This should fragment the library to have a distribution between 300 - 700 bp.

- Add the NEBNext Ultra II ligation mix + enhancer to the well. Add 🧪 2.5 µL of the adaptor to the mixture. Mix well and incubate at 🌡️ 20 °C for ⌚ 00:15:00 heated lid off . Post incubation, add 🧪 28.5 µL of 0.1x TE and bring volume to 🧪 100 µL .

- Perform bead clean-up of the ligation mix using the sample purification beads that came with the kit.

Use 🧪 20 µL of beads for 1st bead addition

🧪 10 µL of beads for 2nd bead addition (to select for an insert size distribution of 350 - 600 bp).

- Elute in 🧪 17 µL of dH₂O.

Enrichment PCR

3

Make 10 uM dilutions of the following primers.

dsODN_Enrich_Plus

dsODN_Enrich_Minus

dsODN_Enrich_Adaptor

For each sample set up the following PCRs in a 96-well plate.

PCR 1.1 - Plus_Strand_Enrichment



🧴 7.5 µL Adaptor ligated DNA

🧴 12.5 µL Ultra II Q5 master mix

🧴 2.5 µL dsODN_Enrich_Adaptor

🧴 2.5 µL dsODN_Enrich_Plus

PCR 1.2 - Minus_Strand_Enrichment

🧴 7.5 µL Adaptor ligated DNA

🧴 12.5 µL Ultra II Q5 master mix

🧴 2.5 µL dsODN_Enrich_Adaptor

🧴 2.5 µL dsODN_Enrich_Minus

PCR program:

🌡️ 98 °C ⌚ 00:00:30 Initial denaturation

🌡️ 98 °C ⌚ 00:00:10 Denaturation

🌡️ 70 °C Gradient -1 °C / cycle ⌚ 00:00:30 Annealing

🌡️ 72 °C ⌚ 00:01:00 Extension

7 cycles

🌡️ 98 °C ⌚ 00:00:10 Denaturation

🌡️ 63 °C ⌚ 00:00:30 Annealing

🌡️ 72 °C ⌚ 00:01:00 Extension

13 cycles

🌡️ 4 °C Hold

After end of PCR, use 🧴 1 µL of the PCR end-product for Picogreen / Qubit.

Use 🧴 2 µL to run the sample on a Tapestation (D5000 high-sensitivity) to check for library size.

Each sample library should be a unimodal distribution between 400 - 1300 bp. All libraries should have an average size distribution within 15% across all samples.



Indexing PCR

4

Make 10 uM dilutions of the **Nextera_R1** & **Nextera_R2** primers.

For each sample, mix

🧴 5 ng of **PCR 1.1** +

🧴 5 ng of **PCR 1.2**

Set up the following PCR:

🧴 10 µL Ultra II Q5 master mix

🧴 1 µL Nextera_R1

🧴 1 µL Nextera_R2

🧴 10 ng Template

🧴 0 µL dH₂O As required

🧴 20 µL Total

PCR Program:

🌡️ 98 °C ⌚ 00:00:30 Initial Denaturation

🌡️ 98 °C ⌚ 00:00:10 Denaturation

🌡️ 65 °C ⌚ 00:00:30 Annealing

🌡️ 72 °C ⌚ 00:01:00 Extension

10 cycles


🌡️ 4 °C Hold

After end of PCR, use 🧴 1 µL of the PCR end-product for Picogreen / Qubit.



Library pooling, QC, Sequencing

5

- Pool samples by equal mass (assuming the library size distribution for the samples are similar).
- Do a 1.2x SPRI bead clean-up. Elute in  30 µL 0.1x TE
- Check concentration of the pool by Picogreen/Qubit and quantify library using Kapa Illumina Library quantification kit. Dilute the pool to appropriate molarity required for the sequencing instrument.
- Load the pool on the sequencer and run with 15% PhiX spike-in for base diversity (works on MiSeq, should work with other Illumina 4-color instruments. Might need more PhiX for NextSeq.)