

Jun 18, 2018

Version 1

Gut_All Cells_Human V.1

DOI

<https://dx.doi.org/10.17504/protocols.io.q3xdypn>

Kylie James¹

¹Wellcome Sanger Institute



Kylie James

Wellcome Sanger Institute

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.q3xdypn>

Protocol Citation: Kylie James 2018. Gut_All Cells_Human. protocols.io <https://dx.doi.org/10.17504/protocols.io.q3xdypn>

Manuscript citation:

Jenke et al, 2013, PLOS one

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

I use this protocol for dissociation of non-immune and immune cells from human colon samples and it works well.



Created: June 18, 2018

Last Modified: June 18, 2018

Protocol Integer ID: 13143

Keywords: gut-all cells-human this protocol, gut-all cells-human, human adult colon sample, good results for scrnaseq, scrnaseq, cell

Abstract

This protocol has been tested by me on human adult colon samples and given good results for scRNAseq. I suggest that you do an enrichment step afterwards, because the dissociation can be quite variable results in terms of the fraction of immune vs non-immune cells that are retrieved.

Materials

MATERIALS

☒ D-PBS (Without Ca++ and Mg++) 500 mL **STEMCELL Technologies Inc. Catalog #37350**

☒ Liberase DH **Roche Catalog #5401054001**

☒ Hyaluronidase **Merck Millipore (EMD Millipore) Catalog #385931-25KU**

☒ HBSS Media **Merck MilliporeSigma (Sigma-Aldrich) Catalog #41122100-007**

☒ DNase II **Merck MilliporeSigma (Sigma-Aldrich) Catalog #D8764-30KU**



- 1 Receive samples. Take about 4mm² sample for OCT

- 2 Wash samples twice with HBSS medium. In petri dish, mince sample into 0.5cm pieces.

- 3 Incubate the samples in 2ml HBSS medium containing 1.07 Wunsch units/ml (2.14wu) Liberase DH (Roche Applied Science, Germany) and 70 U/ml hyaluronidase (Calbiochem/Merck, Germany) for 45 min at 37C on a shaking platform at 600 rpm.

- 4 Homogenise the sample every 15 min the solution by pipetting up and down with a 1ml pipette for 1 min.

- 5 Pass cells through a 40mM sieve. Spin them down at 300 g at 4C for 10 min.

- 6 Resuspended in 300ml PBS buffer containing 100 IU/ml DNase II (Sigma Aldrich, Missouri, USA).

- 7 Proceed with bead enrichment as per manufacturers protocols. Typically this will be a Miltenyi CD45+ enrichment.

- 8 Count cells and dilute to the desired concentration with PBS.

- 9 Load cells (aim for 5000/run) on to the 10X.

- 10