



Oct 11, 2022

H&E Staining for 10X Genomics Visium Imaging

DOI

dx.doi.org/10.17504/protocols.io.4r3l2owqqv1y/v1

Stephen Fisher¹, Marielena Grijalva¹, Rong Guo¹, Sarah A Johnston¹, Hieu Nguyen¹, John Renz², Jean G Rosario¹, Steven Rudich², Brian Gregory¹, Junhyong Kim¹, Kate O'Neill¹

¹University of Pennsylvania; ²Gift of Life Donor Program



Stephen Fisher

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.4r3l2owqqv1y/v1>

Protocol Citation: Stephen Fisher, Marielena Grijalva, Rong Guo, Sarah A Johnston, Hieu Nguyen, John Renz, Jean G Rosario, Steven Rudich, Brian Gregory, Junhyong Kim, Kate O'Neill 2022. H&E Staining for 10X Genomics Visium Imaging. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.4r3l2owqqv1y/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 25, 2022

Last Modified: October 11, 2022

Protocol Integer ID: 67567

Keywords: 10x genomics visium imaging, 10x genomics visium imaging this protocol, 10x genomics visium slide, imaging for visium spatial protocol, 10x genomic, 10x genomics protocol, visium slide, visium spatial protocol, genomic, imaging, fiducial frame of the visium slide, ensuring target tissue, target tissue

Funders Acknowledgements:

NIH

Grant ID: U54HD104392

Abstract

This protocol describes H&E staining of 10X Genomics Visium slides prior to imaging and is adapted from 10X Genomics protocol documented in "[Methanol Fixation, H&E Staining & Imaging for Visium Spatial Protocols](#)". H&E staining and imaging is essential to ensuring target tissue/region/cells are mounted within the fiducial frame of the Visium slide.



Materials

Eosin Mix: 1mL

- 100µL Eosin Y solution (Sigma-Aldrich, HT110216-500mL)
- 900µL Tris-Acetic Acid Buffer (see below)
- *Vortex to mix*
- *Prepare fresh for each use*

⊗ Eosin Y solution aqueous **Merck Catalog #HT110216-500ML**

Tris-Acetic Acid Buffer: 200mL

- Dissolve 11g Tris base in 100mL nuclease-free water (Fisher, BP152-500)
- Adjust pH to 6.0 using 100% Acetic Acid (Fisher, A38-212)
- Bring volume to 200mL with nuclease-free water
- *Filter through 0.2µm Corning 250mL Vacuum System*
- *Store at room temperature for up to 12 months*

⊗ Tris Base **Fisher Scientific Catalog #BP152**

⊗ Acetic acid Glacial **Fisher Scientific Catalog #A38-212**

⊗ Nuclease-free Water

Protocol:

- Tissue slice from OCT embedded tissue block, mounted on a 10X Visium slide
- ⊗ Methanol suitable for HPLC ≥99.9% **Merck MilliporeSigma (Sigma-Aldrich) Catalog #34860**
- ⊗ Mayer's Hematoxylin **Dako Catalog #S3309**
- ⊗ Shandon® Bluing Reagent **Thermo Fisher Catalog #6769001**
- ⊗ Eosin Y solution aqueous **Merck Catalog #HT110216-500ML**
- Tris-Acetic Acid Buffer (see above)
- Eosin Mix (see above)
- ⊗ Nuclease-free water or water filtered using a Milli-Q filtering system **Ambion Catalog #AM9932**
- 10X Visium Thermocycler Adaptor (part of 10X Genomics Visium Accessory Kit; 1000194)

⊗ Visium Accessory Kit **10x Genomics Catalog #1000194**

Equipment

Thermal cycler

NAME

T100 PCR thermal cycler

TYPE

BIO-RAD

BRAND

<https://www.bio-rad.com/pt-br/product/t100-thermal>

SKU

<https://www.bio-rad.com/pt-br/product/t100-thermal-cycler?ID=LGTWGIE8Z>^{LINK}



Troubleshooting

Before start

Prepare Eosin mix fresh for each use. The tris-acetic acid buffer can be stored for up to 12 months.



Tissue Fixation

- 1 Prechill methanol (40mL/slide, dispensed in a 50-mL centrifuge tube) to -20 °C .
- 2 Place a Thermocycler Adaptor on a thermal cycler set at 37 °C and equilibrate for 00:05:00 . Heating the Thermocycler lid is not required. 5m
- 3 Remove slide from -80°C and place on dry ice in a sealed container.
- 4 Place slide on the Thermocycler Adaptor with the active surface facing up and incubate 00:01:00 at 37 °C . 1m
- 4.1 DO NOT close the thermocycler lid. Maintain thermal cycler at 37 °C .
- 5 Remove slide from Thermocycler Adaptor and if necessary, wipe excess liquid from the back of the slide, without touching the tissue sections.
- 6 Completely immerse the slide in the prechilled -20 °C methanol.
- 6.1 Secure the tube cap to prevent methanol loss.
- 7 Incubate upright for 00:30:00 at -20 °C . 30m

Tissue H&E Staining

19m



8 Dispense the following volumes of Milli-Q water:

- a. 500mL in Beaker 1
- b. 800mL in Beaker 2
- c. 800mL in Beaker 3
- d. 800mL in Beaker 4

NOTE: Dispensed volume in each beaker can be used for two slides.

9 Prepare Eosin Mix.

9.1 DO NOT add pure eosin to tissue sections.




10 Remove slide from methanol and wipe excess liquid from the back of the slide, without touching the tissue sections.

11 Place on a flat, clean, nonabsorbent work surface.

Note: Some residual droplets may remain.

12 Add 500µL isopropanol to uniformly cover all tissue sections on the slide.

13 Incubate  00:01:00 at room temperature.

1m


Tip: When incubating the slide with reagents, ensure that the slide is not in contact with any absorbent surface, like laboratory wipes, which may absorb the reagents.

14 Discard reagent by draining and/or holding the slide at an angle with the bottom edge in contact with a laboratory wipe.

15 Wipe excess liquid from the back of the slide, without touching the tissue sections.

16 Place on a flat, clean, nonabsorbent work surface.


Note: Some droplets may remain.

17 Air dry the slide for  00:04:00 .

4m

18 Add 1mL Hematoxylin to uniformly cover all tissue sections on the slide.



19 Incubate  00:05:00 at room temperature.

5m

20 Discard reagent by draining and/or holding the slide at an angle with the bottom edge in contact with a laboratory wipe.

21 Immerse the slide 5x in the water in Beaker 1.

22 Immerse the slide 15x in the water in Beaker 2.


23 Immerse the slide 15x in the water in Beaker 3.

24 Wipe excess liquid from the back of the slide without touching the tissue section.

25 Place on a flat, clean, nonabsorbent work surface.

Note: Some droplets may remain.

26 Add 1mL Bluing Buffer to uniformly cover all tissue sections.

27 Incubate  00:02:00 at room temperature.

2m

28 Discard reagent by draining and/or holding the slide at an angle with the bottom edge in contact with a laboratory wipe

29 Immerse the slide 5x in the water in Beaker 3.


30 Wipe excess liquid from the back of the slide without touching the tissue section.

31 Place on a flat, clean, nonabsorbent work surface.



Note: Some droplets may remain.

32 Add 1mL Eosin Mix to uniformly cover all tissue sections.

33 Incubate  00:02:00 at room temperature.



2m

34 Discard reagent by draining and/or holding the slide at an angle with the bottom edge in contact with a laboratory wipe.

35 Immerse the slide 15x in the water in Beaker 4.

36 Wipe the back of the slide with a laboratory wipe.

37 Place on a flat, clean, nonabsorbent work surface and air dry until tissue is opaque.

38 Incubate slide on the Thermocycler Adaptor with the thermal cycler lid open for  00:05:00 at  37 °C .

5m

39 Proceed to tissue imaging.

Note: Ensure that the entirety of the tissue slice is in the same focal plane before imaging, to reduce the risk of stitching-induced image artifacts hindering downstream analyses.