

Dec 03, 2024

HA Immunoprecipitation

 [The Journal of Cell Biology](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4wbo8vo5/v1>

Amanda Bentley-DeSousa^{1,2,3,4,5}, Shawn Ferguson^{1,2,3,4,5,6}

¹Department of Cell Biology, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

²Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA;

³Program in Cellular Neuroscience, Neurodegeneration and Repair;

⁴Wu Tsai Institute Yale University School of Medicine, New Haven, Connecticut 06510, USA;

⁵Aligning Science Across Parkinson's (ASAP) Collaborative Research Network, Chevy Chase, MD, 20815, USA;

⁶Kavli Institute for Neuroscience, Yale University School of Medicine, New Haven, Connecticut 06510, USA



Amanda Bentley-DeSousa

Yale University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4wbo8vo5/v1>

Protocol Citation: Amanda Bentley-DeSousa, Shawn Ferguson 2024. HA Immunoprecipitation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4wbo8vo5/v1>

Manuscript citation:

<https://doi.org/10.1101/2023.10.31.564602>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 02, 2024

Last Modified: December 03, 2024

Protocol Integer ID: 113417

Keywords: Immunoprecipitation, HA, ha immunoprecipitation this protocol details step, ha tag immunoprecipitation, ha immunoprecipitation, protocol details step, protocol detail, protocol

Funders Acknowledgements:

Aligning Science Across Parkinson's Disease

Grant ID: ASAP-000580

Abstract

This protocol details steps taken to perform HA tag immunoprecipitations.

Materials

DMEM (Thermo Fisher Scientific, 11965-092)

FBS (Thermo Fisher Scientific, 16140-071)

Penicillin/Streptomycin (Thermo Fisher Scientific, 15140122)

CellStripper (Corning, 356230)

PBS [1.1 mM KH₂PO₄ (J.T. Baker, 3246-01), 155.2 mM NaCl (Sigma-Aldrich, 3624-05), and 3 mM Na₂HPO₄ (J.T. Baker, 3828-05)]

Lysis Buffer [50 mM Tris Base (American Bio, AB02000-05000), 150 mM NaCl (Sigma-Aldrich, 3624-05), 1% Triton X-100 (Sigma-Aldrich, X100), 1 mM EDTA (Sigma-Aldrich, 03690) supplemented with cOmplete mini EDTA-free protease inhibitor (Roche, 11836170001) and PhosSTOP (Roche, 4906837001)

Coomassie Plus Protein Assay Reagent (ThermoFisher Scientific, 23236)

Laemmli Buffer [80 mM Tris-HCl pH 6.8 (American Bio, AB020000-05000 / HCl, Sigma-Aldrich, 3624-06), 25.3% Glycerol (American Bio, AB00751), 2.67% SDS (American Bioanalytical, AB01920-00500), and Bromphenol Blue (Sigma-Aldrich, B5525) supplemented with 6.187% fresh β-mercaptoethanol (Sigma-Aldrich, M3148)]

anti-HA beads (Thermo Fisher Scientific, 88837)

TBST [10 mM Tris Base (American Bio AB02000-05000), 150 mM NaCl (Sigma-Aldrich, 3624-05), 0.1% Tween 20 (Sigma-Aldrich, P7949)]



Day 0

- 1 Split a ~90% confluent 15-cm dish of RAW 264.7 cells into 2 15-cm dishes.

Day 1

1h 19m

- 2 Add the indicated compound at the provided concentration and time point in fresh DMEM media containing FBS and P/S.
- 2.1 Proceed to next step if not treating cells.
- 3 Wash each dish 2X with cold PBS on ice.
- 4 Remove the PBS, scrape each well with  700 μL of ice-cold lysis buffer, and transfer the lysates to Eppendorf tubes.
- 4.1 Lysis Buffer: 50 mM Tris Base, 150 mM NaCl, 1% Triton X-100, and 1 mM EDTA supplemented with cOmplete mini EDTA-free protease inhibitor and PhosSTOP.
- 5 Incubate on  00:05:00 . 5m
- 6 Centrifuge  14000 rpm, 4°C, 00:08:00 lysates. 8m
- 7 Measure the protein concentration of the supernatant using Coomassie Plus Protein Assay Reagent as per the manufacturer's instructions.
- 8 To generate the Whole Cell Lysate: mix an aliquot of the supernatant/lysate 1:1 with Laemmli Buffer and heat at  95 °C for  00:03:00 . 3m
- 8.1 Laemmli mBuffer: 80 mM Tris-HCl pH 6.8, 25.4% Glycerol, 2.67% SDS, and Bromophenol Blue supplemented with 6.187% fresh B-mercaptoethanol.



- 9 Create a mastermix of anti-HA beads ( 10 μL of beads per sample) and wash 3X with ice-cold lysis buffer.
- 10 Immunoprecipitate lysates using  10 μL anti-HA beads by incubating rotating end-over-end for  01:00:00 at  4 $^{\circ}\text{C}$. 1h
- 10.1 Ensure the same amount of protein was used for each sample and where needed, supplement samples with lysis buffer to ensure the same volume is used.
- 11 Wash beads 3X with 0.1% TBST.
- 11.1 TBST: 10 mM Tris Base, 150 mM NaCl, 0.1% Tween 20.
- 12 Wash beads 1X with mqH₂O and transfer to new Eppendorf tubes.
- 13 Elute samples by incubating beads with  40 μL Laemmli buffer and warming at  95 $^{\circ}\text{C}$ for  00:03:00 . 3m
- 14 Transfer the eluant to a new Eppendorf tube and supplement with 6.187% B-mercaptoethanol.
- 15 Subject samples to electrophoresis and immunoblotting as in Protein Lysate and Immunoblotting.