

May 24, 2023

## Halo assay to assess mitophagy

DOI

[dx.doi.org/10.17504/protocols.io.dm6gpjzm8gzp/v1](https://dx.doi.org/10.17504/protocols.io.dm6gpjzm8gzp/v1)

nguyen.tha<sup>1</sup>

<sup>1</sup>Laboratory of Michael Lazarou, Walter and Eliza Hall Institute of Medical Research, Parkville, Victoria, Australia

ASAP workspace



nguyen.tha

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gpjzm8gzp/v1>

**Protocol Citation:** nguyen.tha 2023. Halo assay to assess mitophagy. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.dm6gpjzm8gzp/v1>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

We use this protocol and it's working

**Created:** December 05, 2022

**Last Modified:** June 01, 2024

**Protocol Integer ID:** 73537

**Keywords:** Halo assay, mitophagy, ASAPCRN, using halo assay, halo assay, mitophagy, mitophagy this protocol, assay, halo, mizushima lab

**Funders Acknowledgements:**

**Aligning Science Across Parkinson's**

Grant ID: ASAP-000350

## Abstract

This protocol describes how to assess mitophagy using Halo assay developed by Mizushima lab (DOI: [10.7554/eLife.78923](https://doi.org/10.7554/eLife.78923)).

## Attachments



[585-1231.docx](#)

20KB

## Materials

### Buffers and reagents:

#### *Growth media:*

A	B
DMEM with 10% FBS	
Glucose (Sigma, G8769)	4.5 g/l
1x GlutaMAX™ (ThermoFisher, 35050061)	
1x MEM NEAA (ThermoFisher, 11140-050)	
HEPES (1688449)	25 mM

- Antimycin-A (Sigma, A8674; made up in 100% Ethanol to 20 mg/ml),
- Oligomycin (Calbiochem, 495455; made up in DMSO to 10 mg/ml)
- qVD (MedChemExpress, HY-12305; made up in DMSO to 10 mM)
- TMR-conjugated Halo ligand (Promega, G8251)
- *Lysis buffer*: (diluted from 4x LDS (NP007; ThermoFisher)

A
1x LDS
0.1 M DTT

can be aliquoted and stored at -20 or -80°C.

- 4-12% Bis-Tris NuPAGE gels (ThermoFisher)
- NuPAGE™ Antioxidant (NP0005, ThermoFisher; use 0.5 ml/ 200 ml of gel running buffer)
- 20x NuPAGE™ MOPS SDS running buffer (NP001, ThermoFisher)
- 20x NuPAGE transfer buffer (NP00061, ThermoFisher)

#### *PVDF destain:*

A
40% Methanol
7% Acetic Acid



- 1x PBS
- 1x PBS/0.1% Tween20 (PBS/Tween)
- *Blocking buffer*: 5% skim milk in PBS/Tween (make fresh)
- VCP (Cell Signaling, 2649), HALO (Promega, G9211)

⊗ 45% D-( )-Glucose **Merck MilliporeSigma (Sigma-Aldrich) Catalog #G8769**

⊗ GlutaMAX<sup>®</sup> Supplement **Thermo Fisher Catalog #35050061**

⊗ MEM Non-Essential Amino Acids Solution (100X) **Thermo Fisher Scientific Catalog #11140050**

⊗ HEPES Buffer 1M Solution Cell Culture Grade MP Biomedicals **Fisher Scientific Catalog #ICN1688449**

⊗ Antimycin A from Streptomyces sp. **Merck MilliporeSigma (Sigma-Aldrich) Catalog #A8674**

⊗ Oligomycin **Merck MilliporeSigma (Sigma-Aldrich) Catalog #495455**

⊗ Q-VD-OPh **MedChemExpress Catalog #HY-12305**

⊗ HaloTag(R) TMR Ligand, 30ul **Promega Catalog #G8251**

⊗ NuPAGE<sup>™</sup> LDS Sample Buffer **Thermo Fisher Scientific Catalog #NP0007**

⊗ NuPAGE Antioxidant **Thermo Fisher Scientific Catalog #NP0005**

⊗ NuPAGE<sup>™</sup> MOPS SDS Running Buffer (20X) **Thermo Fisher Scientific Catalog #NP0001**

⊗ Tris-Glycine Transfer Buffer **Invitrogen - Thermo Fisher Catalog #NP00061**


⊗ VCP (7F3) Rabbit mAb **Cell Signaling Technology Catalog #2649**

⊗ Anti-HaloTag(R) Monoclonal Antibody **Promega Catalog #G9211**

## Troubleshooting



## Procedures

4h 57m 10s

- 1 Generating cells expressing mitochondrially targeted Halo-GFP using pSu9-Halo-mGFP from Mizushima lab (Addgene #184905; DOI: [10.7554/eLife.78923](https://doi.org/10.7554/eLife.78923)).
- 2 Seed HeLa cells the day before the treatment day in 6 well plates.
  - Each well contained  2 mL of growth media;
  - Seed 350,000 cells for penta KO expressing BFP-Parkin and GFP-OPTN or -NDP52 and 380,000 cells for other knockout lines such as ATG13 KO/penta KO expressing GFP-NDP52;
  - Adjust the number of cells of other cell lines, so that the next day they are all in similar confluency with penta KO expressing BFP-Parkin and GFP-OPTN or -NDP52.

### Note

- Remember to include a set of samples as your non-mitophagy-induced controls

- 3 The next day, make sure the seeded cells are spreading out (not concentrated in the middle of the well because this can affect the results).
- 4 Aspirate off the old media and treat each well with  1 mL growth media containing  50 nanomolar (nM) TMR-conjugated Halo ligand.

### Note

- Aliquot the ligand into small aliquots and avoid freeze-and-thaw cycles.
- The ligand is sensitive to light so put the leftover away immediately after use and keep the media in the dark if you need to do a time course treatment.

- 5 Incubate in a normal TC incubator for  00:20:00 .

20m



- 6 Aspirate off the media and wash thoroughly.







- 6.1 Aspirate off the media and wash thoroughly with 1x PBS. (1/2)



6.2 Aspirate off the media and wash thoroughly with 1x PBS. (2/2)


6.3 Harvest the non-mitophagy-induced samples immediately by scraping (see step 8).

7 For mitophagy-induced samples, treat each well with  2 mL of growth media containing  4 micromolar ( $\mu\text{M}$ ) Antimycin A,  10 micromolar ( $\mu\text{M}$ ) Oligomycin and  10 micromolar ( $\mu\text{M}$ ) QVD for desired period.




#### Note

Make sure all drugs are vortexed well, mix the media well after adding each drug.

8 After the treatment, harvest the cells  On ice by scraping.

#### Note



If you scan for fluorescence signal, from this step on consider keeping the samples away from light as much as possible.

8.1 Pre-chill eppies and 1x PBS  On ice .

#### Note

I normally put all the plates that need harvesting into a fridge and harvest one by one

 On ice .

8.2 Aspirate the media thoroughly from the wells, wash the wells with  1 mL of cold 1x PBS\*, aspirate off the PBS and add  1 mL of cold 1x PBS.



**Note**


\*Make sure swirl around after adding the PBS to wash the cells properly.

8.3 After that, use a plastic cell scraper to scrape all the cells off the wells (I use one scraper for each well. You can wash and reuse them again). Transfer the cells-containing PBS to eppies.

8.4 Centrifuge the eppies at  3000 x g, 4°C, 00:02:00 . Aspirate off PBS.



2m



8.5 Quickly centrifuge for  00:00:10 to spin down the residual PBS. Aspirate off all the PBS.

10s





9 Lyse the cell pellets in lysis buffer and boil the samples at  99 °C with shaking for  00:07:00 .

7m

**Note**

I use the plastic clips to make sure that the lids won't pop during heating.



10 Let the samples cool down and spin at max speed (  Room temperature ) for  00:01:00 .

1m





11 Estimate the protein concentration by nanodrop.

**Note**

Make sure the concentrations do not exceed  6 µL . If they do, dilute with lysis buffer and reheat them for a couple of minutes at  99 °C with shaking.





12 Aliquot  25 µg of each sample into a new eppie and add 1x LDS to make up to  15 µL .





13 Set up the gel tank with MOPs buffer.

- The inside chamber should be filled with 1x MOPs supplemented with antioxidants.
- The outside chamber doesn't need antioxidants.
- Wash each well with a glass syringe.





14 Load markers and samples into the wells and run at 100V for  00:10:00 and 190V for  00:55:00 .

1h 5m

15 Gels were then subjected to wet transfer onto PVDF membrane using cold NuPAGE transfer buffer containing 20% Methanol for  01:00:00 at  Room temperature .




1h

16 After transfer,

16.1 Incubate PVDF membrane with PVDF destain buffer on a shaker at  Room temperature for  00:02:00 .

2m



16.2 wash with PBS/Tween:  
wash with PBS/Tween for  00:05:00 . (1/3)  
wash with PBS/Tween for  00:05:00 . (2/3)  
wash with PBS/Tween for  00:05:00 . (3/3)

15m



16.3 Block with blocking buffer for  00:15:00 .

15m

17 Remove blocking buffer,


17.1 Rinse:  
Rinse with PBS/Tween. (1/2)  
Rinse with PBS/Tween. (2/2)

17.2 Wash:  
Wash with PBS/Tween for  00:05:00 . (1/2)

10m







Wash with PBS/Tween for  00:05:00 . (2/2)



17.3 Wash with 1x PBS for  00:05:00 .

5m



18 Cut the PVDF membrane and put appropriate parts into different antibodies (in this case, it's VCP (1/1000) and HALO (1/1000) antibodies made up in 3% BSA in PBS/Tween) to incubate on a  4 °C shaker  Overnight .

5m






#### Note

To make sure we don't lose antibodies, I wet the tubs with PBS/Tween before putting in the antibodies.


19 The next day, recycle the antibodies back to their tubes,

19.1 Wash:

- wash the blots with PBS/Tween for  00:05:00 . (1/3)
- wash the blots with PBS/Tween for  00:05:00 . (2/3)
- wash the blots with PBS/Tween for  00:05:00 . (3/3)

15m





19.2 Incubate with appropriate HRP-conjugated secondary antibodies (1/10000 for HALO and 1/5000 for VCP) made up in blocking buffer for  01:00:00 .

1h



20 Wash the blots and develop.



20.1 Wash the blots,  
wash the blots with PBS/Tween for  00:05:00 . (1/2)  
wash the blots with PBS/Tween for  00:05:00 . (2/2)

10m

20.2 Wash the blots with 1x PBS for  00:05:00 .

5m

20.3 Develop the blots with ECL prime.

