



Jun 28, 2023

Halo-LC3B processing assay to assess autophagy

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl4qexzvo5/v1>

Xuefeng Ren¹

¹Laboratory of James H. Hurley, University of California, Berkeley, CA



Xuefeng Ren

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl4qexzvo5/v1>

Protocol Citation: Xuefeng Ren 2023. Halo-LC3B processing assay to assess autophagy. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl4qexzvo5/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 12, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 83241

Keywords: Halo-LC3B processing assay, autophagy, ASAPCRN, autophagy this protocol details halo, autophagy, lc3b processing assay, lc3b, assay, protocol details halo

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-000350

Abstract

This protocol details Halo-LC3B processing assay to assess autophagy.

Attachments



[738-1854.pdf](#)

115KB

Guidelines

Reference: DOI: 10.7554/eLife.78923

Materials

Buffers and reagents:

Growth media: DMEM medium with GlutaMAX containing 10% FBS and 10% Pen-Strep.

1x PBS

Lysis buffer:

A
25 mM HEPES pH 7.5
200 mM NaCl
2 mM MgCl ₂
10% glycerol
1 mM TCEP
0.2% n-dodecyl-β-D-maltoside
pierce protease inhibitors
Benzonase

⊗ DMEM, high glucose, GlutaMAX[®] Supplement **Thermo Fisher Catalog #10566016**

⊗ Gibco[™] Fetal Bovine Serum value heat inactivated (formerly USDA-approved in North America or quali **Fisher Scientific Catalog #A5256801**

⊗ Penicillin-Streptomycin (10,000 U/mL) **Gibco - Thermo Fisher Scientific Catalog #15140122**

⊗ EBSS, calcium, magnesium, phenol red **Thermo Fisher Catalog #24010043**

⊗ Janelia Fluor[®] HaloTag[®] Ligands **Promega Catalog #GA1110**

⊗ Gibco[™] Trypsin-EDTA (0.05%) phenol red **Fisher Scientific Catalog #25-300-120**

⊗ Pierce Protease Inhibitor Tablets **Thermo Fisher Catalog #A32963**

⊗ n-Dodecyl-β-D-maltoside (DDM) **Gold Biotechnology Catalog # DDM5**



 Benzonase® Nuclease Purity > 90% **Merck Millipore (EMD Millipore) Catalog #70746-4**

 NuPAGE® 4-12% Bis-Tris Protein Gels, 1.0 mm, 12-well **Thermo Fisher Catalog #NP0322BOX**



Halo-LC3B processing assay to assess autophagy

1h 50m

- 1 Generating HeLa cells expressing HaloTag-LC3B using pMRX-IP-HaloTag7-LC3 from Mizushima lab (Addgene #184899; DOI: 10.7554/eLife.78923).
- 2 Seed HeLa cells at 100-150K cells/well in 12-well plate one day before.
- 3 Next day, incubate cells with complete DMEM medium and [IM] 50 nanomolar (nM) JF646 HaloTag ligand (Promega) for ⌚ 01:00:00 , and then wash twice with 1xPBS. The non-starved samples can be harvested immediately by trypsinization (step 5).
- 4 To induce autophagy by starvation, treat cells with EBSS buffer (Gibco) for desired period.
- 5 After the treatment, harvest the cells by trypsinization.
- 5.1 Wash the wells with 1x PBS.
- 5.2 Incubate the wells with 🧪 0.5 mL trypsin at 🌡️ 37 °C for ⌚ 00:05:00 .
- 5.3 Add 🧪 0.5 mL complete medium into each well.
- 5.4 Transfer cells into pre-chilled Eppendorf tube, spin 🌀 2000 x g, 00:05:00 .
- 5.5 Aspirate off the liquid.
- 6 Resuspend the cell pellets in 🧪 30 µL of lysis buffer and incubate 🌡️ On ice for ⌚ 00:30:00 .

1h



5m



5m



30m





- 7 Centrifuge the cell lysate at  21000 x g, 00:10:00 . Transfer the cleared lysate into another tube, and measure protein concentration nanodrop spectrophotometer (Thermo Fisher).  10m 
- 8 For each sample, load  20 µg clarified lysates onto NuPAGE 4-12% Bis-Tris Gel (Thermo Fisher).
- 9 For in-gel fluorescence imaging, the gel was immediately visualized with ChemiDoc MP imaging system (Bio-Rad) after SDS PAGE. Band intensities are acquired by exciting samples at 546 nm (mCherry signal) and 647 nm (JF646 HaloTag ligand signal).