


Sep 12, 2023

Harvest of adherend cells from Cytodex® 1 beads

DOI

dx.doi.org/10.17504/protocols.io.6qpvr4w7pgmk/v1



David König¹, Kim Kremser¹, Sina Bartfeld^{1,2}, Shirin Kadler²

¹Medical Biotechnology, Technische Universität Berlin, Berlin, Germany;

²Si-M / "Der Simulierte Mensch" a science framework of Technische Universität Berlin and Charité - Universitätsmedizin Berlin, Berlin, Germany



David König

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.6qpvr4w7pgmk/v1>

Protocol Citation: David König, Kim Kremser, Sina Bartfeld, Shirin Kadler 2023. Harvest of adherend cells from Cytodex® 1 beads. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.6qpvr4w7pgmk/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: March 13, 2023

Last Modified: September 12, 2023

Protocol Integer ID: 78658

Keywords: micro carrier, Cytodex1, adherent cells, MSC, Chondrocytes, Fibroblasts, dextranase, cell culture, dextran, cell harvest, harvest of adherend cell, adherend cells from dextran, expanding adherend cell, suitable for cell therapy, adherend cell, cells with high harvesting yield, cell therapy, other enzymatic schemes like trypsin, adherend cells to quantety, varouse cell type, cell, dextran, other enzymatic scheme, trypsin

Funders Acknowledgements:

ZIM Funding Project - Federal Ministry For Economic Affairs And Climate Action

Grant ID: KK5145501AD0

Abstract

Expanding adherend cells to quanteties suitable for cell therapy is a subject of increasing relevance. Varouse cell types either autolog or allogeneic - could in future be used for medical purposes / in clinical application as a key component in therapy. This protocol describes how to harvest adherend cells from dextran based Cytodex® 1 beads after expansion. Since other enzymatic schemes like Trypsin-EDTA or TrypLE™ have shown low efficiency and yield - here we present an approach that is easy in implementation and delivers cells with high harvesting yield and high viability.

Guidelines

This protocol can surely be further improved but besides that is pretty much straight forward - crucial steps are the centrifugation and supernatant discarding where many cells can be lost if samples are not handled with due care.



Materials

DMEM - 4.5g/L glucose - stable glutamine - sodium pyrovate - Biowest [L0103-500]

10% FBS - Corning

1x P/S - Corning

PBS

10x TrypLE Express - Gibco

Dextranase Plus L, 1,6- α -D-Glucan-6-glucanohydrolase from *Chaetomium erraticum*

Novozymes - >100 KDU-A/G [Dextranase units]

Cytiva - Cytodex1 beads [17044803]

Clean bench

Incubator - 37°C - 5% CO₂

Centrifuge

Centrifugation tubes 15/50mL

Eppendorf tubes 1.5mL

Tube rotator

Pipetboy + pipets [5-50mL]

Pipets - [10-1000uL]

NucleoCounter NC-200 + Solution A100 and B

Neubauer Chamber

Vacuum pump

Spinner Flask - Pfeiffer 500mL

Magnetic stirrer - Technomara

adherent primary cells - MSC / Chondrocytes / HUVEC / Fibroblasts

Troubleshooting


Safety warnings

 There are no particular safety risks exceeding the risk during cell culture.






sampling

1h 26m

- 1 withdraw homogenized cell solution from culture vessel and transfer to centrifuge tube  10-50 mL
- 1.1 take a sample from each tube for counting to determine available cell number






Micro carrier solving

1h 26m

- 2 add  0.1 % volume of dextranase (v/v)
- 3 transfer samples to a tube rotator and place in an incubator  10 rpm, 37°C, 01:00:00 1h
- 4 take out tubes and centrifuge  600 x g, 00:08:00 , acc 9 - dec 5 8m
- 5 carefully discard supernatant

single cells

1h 26m

- 6 add a small volume  1-5 mL of 2x TrypLE™ and completely solve the pallet
- 6.1 add more  9-45 mL of 2x TrypLE™, transfer tubes to a Tube rotator and place in an incubator  10 rpm, 37°C, 00:10:00 10m
- 7 take out tubes and centrifuge  600 x g, 00:08:00 , acc 9 - dec 5 8m
- 8 discard supernatent and solve pallet in adequat amount of culture medium  1-5 mL



8.1 filtrate cell solution through a 40 μ m Filter to exclude micro carrier fragments

Yield, viability and follow up

9 count cells and determine viability

10 use cells for further analysis