

Oct 16, 2019

Harvesting OP50 for C. elegans liquid culture

DOI

dx.doi.org/10.17504/protocols.io.xcmfiu6

Vidur Sabharwal¹

¹Tata Institute of Fundamental Research



Vidur Sabharwal

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.xcmfiu6>

Protocol Citation: Vidur Sabharwal 2019. Harvesting OP50 for C. elegans liquid culture. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.xcmfiu6>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: January 22, 2019




Last Modified: October 16, 2019

Protocol Integer ID: 19565

Keywords: harvesting op50, liquid culture

Troubleshooting


Harvesting OP50

- 1 After 8 hr incubation, use 2 ml of E. coli OP50 LB culture to inoculate each of 3 × 250 ml LB. Place flasks in shaking incubator (180 rpm) O/N.  17:00:00

 37 °C


- 2 Weigh 3 × 50 ml centrifuge tubes and write weight of the empty tube on the tube.

3

Day 2. Using a serological pipette, aliquot 40 ml of O/N E. coli OP50 culture into the pre weighed centrifuge tubes. Centrifuge in swing bucket centrifuge at 3,220 x g, 10°C, 10 min.  00:10:00

 10 °C

4

Carefully decant the supernatant, repeat step 4 till all the culture is used up. Keep the tube inverted with lid on and leave to stand for several minutes. Using a pipette remove any excess supernatant that may have collected in the lid. Weigh the tube and calculate the weight of the pellet.  [go to step #4](#)

- 5 Calculate the volume needed to provide a suspension of 30 g/L and mark this volume on the tube. Date, label and place tubes at -20 °C for use within 1-3 months or at - 80 °C if storing for longer than 3 months.
- 6 Prepare bacterial suspension for growing nematodes; allow bacterial pellet to thaw out and add required volume of S complete medium to each tube, vortex gently to resuspend pellet, pool contents of different tubes to obtain the volume required. Work under sterile conditions.