



Jan 22, 2026

HCR RNA-FISH protocol for the whole-mount antenna of *Drosophila*



Forked from [HCR RNA-FISH protocol for the whole-mount brains of *Drosophila* and other insects](#)

DOI

<https://dx.doi.org/10.17504/protocols.io.dm6gp3r3pvzp/v1>

Iván Méndez¹

¹Northwestern University



Iván Méndez

Northwestern University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.dm6gp3r3pvzp/v1>

External link: <https://doi.org/10.1371/journal.pgen.1010802>

Protocol Citation: Iván Méndez 2026. HCR RNA-FISH protocol for the whole-mount antenna of *Drosophila*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.dm6gp3r3pvzp/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 30, 2023

Last Modified: January 22, 2026

Protocol Integer ID: 84261

Keywords: drosophila, hcr, in situ, antenna, mount antenna of drosophila, hcr rna, rna fluorescent in situ hybridization, mount samples of drosophila antenna, drosophila antenna, rna fluorescent, drosophila, hcr reagent, mount antenna, fish protocol, hcr, rna, hybridization chain reaction, using hybridization chain reaction, molecular instrument, situ hybridization, probe

Abstract

This is a protocol to perform RNA fluorescent in situ hybridization (RNA-FISH) using hybridization chain reaction (HCR) on whole-mount samples of *Drosophila* antenna. Probes and HCR reagents are purchased from **Molecular Instruments**. This protocol is loosely based on the **"generic sample in solution" protocol published by Molecular Instruments**. Modifications include higher amount of probe and hairpins to boost signal, and the description of fixation conditions.

Citation

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall
(2020). Non-canonical odor coding ensures unbreakable mosquito attraction to humans. bioRxiv.

<https://doi.org/10.1101/2020.11.07.368720>

LINK

Safety warnings

- ⚠ This protocol uses solutions of paraformaldehyde and formamide, which are highly toxic chemicals. Consult the SDS sheets of the reagents used in this protocol for proper handling instructions.



Collection and fixation of antenna

2h 45m

- 1 Take aliquots of fixation solution (4%PFA in PBS+3%TritonX) from the -20C and thaw on ice. Put a tube of 1X PBS in ice to cool it down.
- 2 Fill a dissection glass with EtOH and a second one with cold PBS.
- 3 Collect flies in the fly pad, transfer them to a Falcon tube and put them in ice (this maintains the flies anesthetized but reducing their exposure to CO₂).
- 4 Before dissecting each fly, dip it quickly in EtOH (this removes the wax covering the cuticle making the fly and antenna to sink in the solutions).
- 5 Dissect antenna in cold PBS using forceps. Transfer the antenna using a 20ul pipet with the tip cut, to the fixation solution (keep it in ice). Usually, I dissect as many antenna as I can in 30-40 min - around 24 to 30 antenna is a good number for one in situ run.
- 6 After dissecting enough antenna fix for **30 min at RT** in fixation solution.

30m

🌡 Room temperature Gently shaking ⌚ 00:30:00

Pre-heat an aliquot of Probe Hybridization Buffer @37C if proceeding with hybridization on the same day OR
move to Step 7

- 7 Wash antenna before long term storage. To avoid losing samples during this process, use a dissection scope to remove liquid during washes.

30m

PBT ⌚ 00:05:00

PBT ⌚ 00:05:00

PBT ⌚ 00:05:00

MeOH wash ⌚ 00:05:00

MeOH wash ⌚ 00:05:00

MeOH wash ⌚ 00:05:00

- 8 Label with name of species, tissue, and date before storing in MeOH at 🌡 -20 °C



In theory, tissues can be stored like this for months, we usually proceed to the HCR within a week

Day 1

1h 45m

- 9 For each sample/tube, transfer 550µl of Probe Hybridization Buffer (store at -20C) to a 1.5 mL tube and warm to 37 °C

- 10 For all next steps, washes and rinse use around 1ml unless a specific volume is mentioned. Rinse means adding the solution and removing it as soon as the antenna settle in the bottom

15m

Rehydrate samples:

MeOH 75% /PBT 00:05:00

MeOH 50% /PBT 00:05:00

MeOH 25% /PBT 00:05:00

PBT rinse

- 11 **Post-fix:** Incubate for 10 minutes in 4% paraformaldehyde in PBT

- 12 **Wash:**

PBT 00:05:00

PBT 00:05:00

PBT 00:05:00

15m

- 13 Pre-hybridize samples in 300 µL of Probe Hybridization Buffer for:

00:30:00 37 °C

CAUTION: probe hybridization buffer contains formamide, a hazardous material.

30m

- 14 In the meantime, prepare a probe solution by adding 8 µL of probe to 200µL of warm Probe Hybridization Buffer. Even when using multiple probes, maintain the final of buffer to 200µL.

- 15 Remove the pre-hybridization solution from the sample and add the probe solution from step 9

3m



16 Add probe mixture and Incubate samples



Overnight around 24h

37 °C

Day 2

1h 20m

17 Pre-heat an aliquot of Probe Wash Buffer (2mL per tube/sample)

2m

37 °C

18 Equilibrate an aliquot of Amplification Buffer (stored at 4C) to room temperature (350ul per sample/tube)

2m

Room temperature

19 Remove and **SAVE PROBE SOLUTION**, which can be reused at least 3 times. Save used probe solutions at -20 °C

20 Remove excess probes by washing 37 °C :

50m

400 ul pre-heated Probe Wash Buffer 00:10:00

400 ul pre-heated Probe Wash Buffer 00:10:00

400 ul pre-heated Probe Wash Buffer 00:10:00

Before the next wash, take hairpins from -20C and thaw on ice

400 ul pre-heated Probe Wash Buffer 00:10:00

400 ul pre-heated Probe Wash Buffer 00:10:00

21 In the meantime, prepare **separately** each hairpin (h1 and h2) of each amplifier. Add 5 µL / sample of the 3 µM stock of each hairpin to a **separate** PCR tube. For example, if amplifiers B1 and B2 are being used, prepare 4 PCR tubes that contain B1-h1, B1-h2, B2-h1, and B2-h2, respectively. Incubate the tubes in a thermocycler at 95 °C for 90 sec. Immediately take them out of the machine (while it is still at 95 °C), place them in a rack and incubate them at room temperature **in the dark** for at least 30 min.

30m



00:30:00 or longer

Room temperature



- 22 Wash samples Room temperature : 15m
- 500 µL of 5X SSCT 00:05:00
- 500 µL of 5X SSCT 00:05:00
- 500 µL of 5X SSCT 00:05:00
- 23 Pre-amplify each sample with 250 µL of Amplification Buffer 30m
- 00:30:00
- Room temperature
- 24 In the meantime, add 5 µL of each hairpin from step 21 to 100 µL of Amplification Buffer (keep the volume of Amplification Buffer at 100 µL even if multiple hairpin sets are being used)
- Room temperature
- 25 Remove 250 µL liquid from the samples and add the Amplification Buffer + hairpins from step 18 2m
- 26 Incubate in **the dark (place tubes inside a cooler and in a drawer)**
- Overnight
- Room temperature

Day 3

3h 15m

- 27 **SAVE HAIRPIN MIXTURE.** Can be reused multiple times. Save used hairpin mixtures at -20°C. 20m
- 28 During the next washes minimize exposure to light by keeping the tubes in a drawer or inside a cooler. 20m
- Wash:
- 500 µL of 5X SSCT
- 500 µL of 5X SSCT
- 500 µL of 5X SSCT
- 500 µL of 5X SSCT 00:15:00
- 500 µL of 5X SSCT 00:05:00
- 29 Rinse with 1X Nuclease-Free PBS to remove the detergent



2m

30 Remove the PBS and add VectaShield.

2m

31 Mount on slides using circular holders.

10m

32 Proceed with imaging



Citations

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall. Non-canonical odor coding ensures unbreakable mosquito attraction to humans
<https://doi.org/10.1101/2020.11.07.368720>