Nov 05, 2021

# HCR RNA-FISH protocol for the whole-mount brains of *Drosophila* and other insects

DOI

dx.doi.org/10.17504/protocols.io.bzh5p386

Amanda A. G. Ferreira<sup>1</sup>, Bogdan Sieriebriennikov<sup>2</sup>, Hunter Whitbeck<sup>1</sup>

<sup>1</sup>New York University; <sup>2</sup>NYU Grossman School of Medicine

Desplan Lab

Amanda A G Ferreira

New York University





#### DOI: dx.doi.org/10.17504/protocols.io.bzh5p386

**Protocol Citation:** Amanda A. G. Ferreira, Bogdan Sieriebriennikov, Hunter Whitbeck 2021. HCR RNA-FISH protocol for the wholemount brains of Drosophila and other insects. **protocols.io** <u>https://dx.doi.org/10.17504/protocols.io.bzh5p386</u>

License: This is an open access protocol distributed under the terms of the <u>Creative Commons Attribution License</u>, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working We use this protocol and it's working

Created: October 26, 2021

Last Modified: November 05, 2021

Protocol Integer ID: 54557

Keywords: brain, drosophila, ant, hcr, in situ, larva, rna fish

# Abstract

This is a protocol to perform RNA fluorescent in situ hybridization (RNA-FISH) using hybridization chain reaction (HCR) on whole-mount samples of the brains of the fly *Drosophila melanogaster* and other insects, e.g. the jumping ant *Harpegnathos saltator*. Probes and HCR reagents are purchased from <u>Molecular Instruments</u>. This protocol is loosely based on the <u>"generic sample in solution" protocol published by Molecular Instruments</u>. Our modifications include the description of fixation conditions, counterstaining by Hoechst, and altered washes. Additionally, we use larger concentrations of probes and hairpins following the protocol described by Younger, Herre et al. 2020. We have successfully employed this protocol to stain insect brains with up to 4 different probe sets simultaneously (hairpins conjugated with Alexa Fluor 488, 546, 496, and 647).

#### CITATION

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall (2020). Non-canonical odor coding ensures unbreakable mosquito attraction to humans. bioRxiv.

LINK

https://doi.org/10.1101/2020.11.07.368720

# Guidelines

## WORKING PRACTICES:

Prepare all buffers using nuclease-free water. Use filter tips and nuclease-free tubes. If using spot plates, pre-clean them first with household bleach diluted 1:10 in water and then with 70% ethanol. Wear gloves and adhere to other practices aimed at minimizing RNA degradation in the sample. Working in a clean bench is not required if other RNAse-free practice are followed.

#### **PROBE DESIGN:**

We select the target sequence or isoform of the gene of interest and let Molecular Instruments design the probes. For genes with multiple isoforms, either target the isoform that includes as many as possible constitutive exons and as few as possible alternatively spliced exons, or the isoform that has the highest RNA-seq coverage (assessed visually in IGV) if RNA-seq data are available. Aim for the highest number of probes in a set, ideally 40, although we have successfully performed experiments with probe sets containing <20 probes.

#### AMPLIFIER CHOICE:

We routinely perform multiplexed stainings with up to 4 different probe sets and a Hoechst counterstain. We use amplifiers conjugated with Alexa Fluor 488, 546, 594, and 647. We are able to detect clearly distinguishable signals with minimal bleed-through on our confocal microscope (Leica SP8). However, be aware that simultaneously using fluorophores with partially overlapping spectra (e.g. AF 546 and 594) requires setting narrower detection ranges, which reduces the amount of signal detected.

## Materials

## **REAGENTS TO PURCHASE:**

- X Nuclease-Free Water (not DEPC-Treated) Thermo Fisher Scientific Catalog #AM9937
- 20X PBS (Phosphate Buffered Saline) pH 7.4 growcells.com Catalog #MRGF-6396
- Schneider's Drosophila Medium Thermo Fisher Catalog #21720024
- X Paraformaldehyde, 16% (wt/vol) Electron Microscopy Sciences Catalog #15710

#### Safety information

Paraformaldehyde is toxic, consult the SDS sheet for proper handling instructions

#### Note

Avoid long-term storage of the paraformaldehyde solution after opening the ampoule

X Triton X-100 Sigma Aldrich Catalog #X100

X 10% Tween 20 BIO-RAD Catalog #1662404

20X SSC Quality Biological Catalog #351-003-131

X SlowFade<sup>™</sup> Gold Antifade Mountant Invitrogen - Thermo Fisher Catalog #S36936 - or any other antifade

#### mountant

🔀 Hoechst 33258, Pentahydrate (bis-Benzimide), 100 mg Thermo Fisher Catalog #H1398 - dissolve in DMSO to

5 mg / mL, aliquot and store at -20 °C

- X Methanol Fisher Scientific Catalog #A412-4
- X HCR Probe Hybridization Buffer Molecular Instruments
- X HCR Probe Wash Buffer Molecular Instruments

#### Safety information

Hybridization and Wash Buffers contain formamide, consult the SDS sheet for proper handling instructions

X HCR Amplification Buffer **Molecular Instruments** 

#### **BUFFERS TO PREPARE ("Fly" version of the protocol, see Steps for details):**

#### Note

Prepare fresh using nuclease-free water, store at 4 °C if required after the 1st day of the protocol.

- 0.3% PBSTX (1X PBS with 0.3% v/v Triton X-100)
- Fixation Buffer (4% paraformaldehyde prepared from the 16% solution by diluting it with 0.3% PBSTX, i.e. 1 volume 16% PFA + 3 volumes 0.3% PBSTX)
- 0.3% PBST (1X PBS with 0.3% v/v Tween-20)
- 5X SSCT (5X SSC with 0.1% v/v Tween-20)
- 1X PBS

#### **BUFFERS TO PREPARE ("Ant" version of the protocol, see Steps for details):**

Note

Prepare fresh using nuclease-free water, store at 4 °C if required after the 1st day of the protocol.

- Fixation Buffer (1X PBS with **0.03%** v/v Triton X-100 and 4% paraformaldehyde)
- 0.1% PBST (1X PBS with 0.1% v/v Tween-20)
- 5X SSCT (5X SSC with 0.1% v/v Tween-20)
- 1X PBS

## Safety warnings

This protocol uses solutions of paraformaldehyde and formamide, which are highly toxic chemicals. Consult the SDS sheets of the reagents used in this protocol for proper handling instructions.

Day	1	1h 45m	
1	Prepare all solutions	15m	
2	Pre-heat an aliquot of Probe Hybridization Buffer if proceeding with hybridization on the same day (see step case below)	2m	
3	Dissect brains in cold Schneider's Medium	10m	
	Note		
	Dissection can also be done in 1X Nuclease-Free PBS.		
4	Choose the version of the protocol you would like to follow (Fly is default)		
	STEP CASE		
Dros	sophila melanogaster 28 steps		
We follow this path when working with <i>Drosophila melanogaster</i> brains. This path is 2 days shorter than the other one. Fixation is shorter, there is no dehydration-rehydration, and final washes are in both 5X SSCT and Probe Wash Buffer.			
5	<ul> <li>Fix brains in 800 μL of Fixation Buffer</li> <li>O0:20:00</li> <li>Commensative</li> </ul>	20m	
	<b>G</b> 24 rpm Nutator <b>or G</b> 60 rpm Orbital shaker		
	Note		
	For fixation and all subsequent steps, samples can be placed either in Eppendorf tubes or in wells of a spot plate (e.g. Pyrex spot plate with 9 depressions, Catalog #CLS722085). Tubes are incubated on a nutator and plates are incubated on an orbital (horizontal) shaker.		

6	Rinse 3x with 500 µL of PBST Room temperature	3m
7	Wash 3× 15 min with 500 µL of PBST <ul> <li> <li> <li> <li> <li> <li> <li> <l< td=""><td>45m</td></l<></li></li></li></li></li></li></li></ul>	45m
	3 24 rpm Nutator or 3 60 rpm Orbital shaker	
8	Pre-hybridize samples by incubating them with 250 μL of warm Probe Hybridization Buffer from step 2 O0:10:00 can be extended to 30 min 37 °C	10m
9	In the meantime, prepare a 16 nM probe solution by adding 8 pmol of each probe mixtu (e.g. 8 $\mu L$ of 1 $\mu M$ stock) to the warm Probe Hybridization Buffer for the total volume of 500 $\mu L$	re
	Note	
	4 $\mu L$ of each probe in 250 $\mu L$ total volume should work as well.	
10	Remove the pre-hybridization solution from the sample and add the probe solution from step 9	n 3m
	Take extra care while removing the liquid. Hybridization Buffer is viscous and brains may float.	
11	Incubate samples with the probes Overnight We usually do ~24 h, but incubation can be extended to ~48 h. Minimun  37 °C	The recommended is 12 h

	🗘 24 rpm Nutator or 🗘 60 rpm Orbital shake	r	
Day	/ 2		1h 20m
12	Pre-heat an aliquot of Probe Wash Buffer 37 °C		2m
13	Equilibrate an aliquot of Amplification Buffer to ro Room temperature	oom temperature	2m
14	Remove excess probes by washing 5× 10 min with step 12	th pre-heated Probe Wash Buffer from	50m
15	In the meantime, prepare <b>separately</b> each hairpin / sample of the 3 µM stock of each hairpin to a seamplifiers B1 and B2 are being used, prepare 4 P h1, and B2-h2, respectively. Incubate the tubes in Immediately take them out of the machine (while and incubate them at room temperature <b>in the date</b> 00:30:00 or longer	n (h1 and h2) of each amplifier. Add 7 μL eparate PCR tube. For example, if CR tubes that contain B1-h1, B1-h2, B2- n a thermocycler at 95 °C for 90 sec. it is still at 95 °C), place them in a rack ark for at least 30 min.	30m
16	Wash samples 2× 5 min with 500 µL of 5X SSCT ⓒ 00:05:00 ⓒ Room temperature		10m
	5 24 rpm Nutator or 5 60 rpm Orbital shake	r	

17 Pre-amplify each sample with 250 μL of Amplification Buffer

10m

•				
	👏 00:10:00 can be ex	tended to 30 min		
	Room temperature			
	C 24 rpm Nutator or	🗘 60 rpm Orbital shaker		
18	In the meantime, add 6 (keep the volume of Am used)	μL of each hairpin from step 15 plification Buffer at 100 μL eve	to 100 μL of Amplification Buffer n if multiple hairpin sets are being	
	Room temperature			
19	Remove 250 μL liquid fr step 18	om the samples and add the A	mplification Buffer + hairpins from	2m
	Note			
	Take extra care while re may float.	emoving the liquid. Amplification	Buffer is viscous and the brain	
20	All the next steps are li	aht sensitive and must be do	ne in the dark!	٨
		5		
21	Incubate in <b>the dark</b>			6
	👏 Overnight			
	Room temperature			
	<b>\$</b> 24 rpm Nutator or	🗘 60 rpm Orbital shaker		
Day	3			3h 15m
22	Equilibrate an aliquot of	Probe Wash Buffer to room ter	nperature	20m
23	Add 300 μL of 5X SSCT hairpins. This will make as possible.	to samples, which at this poin the solution inside the tubes le	t still contain Amplification Buffer + ss viscous. Remove as much liquid	2m
24	Wash 1× 5 min with 500	μL of 5X SSCT		5m
	Room temperature			
	<b>(3</b> 24 rpm Nutator or	<b>\$</b> 60 rpm Orbital shaker		

25	Wash 1× 15 min with 500 μL of Probe Wash Buffer	15m
	Room temperature	
	<b>C</b> 24 rpm Nutator <b>or C</b> 60 rpm Orbital shaker	
26	If not using Hoechst, skip this step and the next and proceed to step 28.	2h
	If using Hoechst: incubate samples with 500 $\mu$ L of Probe Wash Buffer + Hoechst (1 $\mu$ L of the 5 mg/mL stock - final concentration 10 $\mu$ g/mL) for 2 h	*
	Soom temperature	
	<b>(</b> ) 24 rpm Nutator <b>or (</b> ) 60 rpm Orbital shaker	
27	Wash 1× 15 min with 500 μL of Probe Wash Buffer	15m
	Room temperature	1
	<b>\$</b> 24 rpm Nutator <b>or \$</b> 60 rpm Orbital shaker	
28	Wash 1× 10 min with 500 μL of 5X SSCT	10m
	Room temperature	
	<b>\$</b> 24 rpm Nutator <b>or \$</b> 60 rpm Orbital shaker	
29	Rinse with 1X Nuclease-Free PBS to remove the detergent	2m
30	Remove the PBS and add Slowfade (or another antifade mountant)	2m
31	If necessary, finish dissecting to remove extra tissue. Mount brains on slides.	10m
32	Proceed with imaging	<u>द</u>

# Citations

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall. Non-canonical odor coding ensures unbreakable mosquito attraction to humans <u>https://doi.org/10.1101/2020.11.07.368720</u>

Meg A. Younger, Margaret Herre, Alison R. Ehrlich, Zhongyan Gong, Zachary N. Gilbert, Saher Rahiel, Benjamin J. Matthews, Leslie B. Vosshall. Non-canonical odor coding ensures unbreakable mosquito attraction to humans <u>https://doi.org/10.1101/2020.11.07.368720</u>