

Jul 25, 2023

HEK293 cell culture for co-immunoprecipitation experiments

DOI

<https://dx.doi.org/10.17504/protocols.io.kxygx3zeog8j/v1>

Dan Dou¹, Erika Holzbaur¹

¹University of Pennsylvania



Dan Dou

University of Pennsylvania

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.kxygx3zeog8j/v1>

Protocol Citation: Dan Dou, Erika Holzbaur 2023. HEK293 cell culture for co-immunoprecipitation experiments. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.kxygx3zeog8j/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: July 20, 2023



Last Modified: May 31, 2024

Protocol Integer ID: 85308

Keywords: ASAPCRN, hek293 cell culture, hek293 cell culture for the purpose, cell culture, cell

Funders Acknowledgements:

ASAP

Grant ID: ASAP-000350

Abstract

Here, we describe HEK293 cell culture for the purpose of co-immunoprecipitation experiments

Before start

Cells should be routinely tested for mycoplasma contamination using MycoAlert detection kit (Lonza, LT07)



- 1 Culture HEK293 cells in DMEM (Corning) supplemented with 10% fetal bovine serum (HyClone)
- 2 Maintain cells at 37 degrees Celsius with 5% CO₂
- 3 In preparation for co-immunoprecipitation experiments, passage HEK cells with Trypsin onto three 10 cm tissue culture dishes per experimental condition, at a density of 1:10
- 4 24 hours later, transfect using FuGENE 6 (6-12 µg total plasmid DNA; Promega)