

Sep 21, 2023

HEK3/LINC01509 Library Preparation

DOI

dx.doi.org/10.17504/protocols.io.5qpvorjxzv4o/v1



Mathew Chu¹

¹Walter and Eliza Hall Institute of Medical Research Melbourne



Mathew Chu

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5qpvorjxzv4o/v1>

Protocol Citation: Mathew Chu 2023. HEK3/LINC01509 Library Preparation. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.5qpvorjxzv4o/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: In development

We are still developing and optimizing this protocol

Created: May 24, 2023

Last Modified: September 21, 2023

Protocol Integer ID: 82388

Keywords: hek3 recording site, linc01509 library preparation library preparation, using pcr, hek3, pcr, linc01509, recording

Funders Acknowledgements:

Shalin Naik

Grant ID: 10771

Abstract

Library preparation for HEK3 recording site (LINC01509) using PCR and digestion.

Materials

- NGS primer F ("CTGGCCTGGGTCAATCCTTG")
- NGS primer R ("CCCAGCCAAACTTGTCAACC")
- Q5 HiFi polymerase 2X master mix (NEB)
- BclI (NEB)
- LMW agarose
- SYBR Safe gel stain
- Generuler 1kb PLUS
- SPRI beads (Nucleomag)

Troubleshooting



HEK293T Cell Lysis

1 Make up **1:50** Thermolabile Proteinase K (NEB) in Viagen (cell).

2 Resuspend cells in  25 μL of Viagen/Proteinase K.



Quickly transfer cells into a PCR strip to avoid viscosity.

3 Perform lysis in a PCR machine:

1h 30m

 37 °C  01:00:00

 55 °C  00:30:00

 On ice



Store lysates at  -20 °C in the pre-PCR freezer.

PCR Amplification of Genomic DNA

4 Make up a  100 μL PCR reaction per sample (Q5) using  10 μL of cell lysate.



4.1  98 °C  00:00:30

30s

4.2  98 °C  00:00:05

35s

 68 °C  00:00:10

 72 °C  00:00:20

4.3  go to step #4.2 24 X

4.4  72 °C  00:02:00

2m

SPRI Bead Concentration



5 Purify amplicons by **1X SPRI** bead cleanup, eluting in  10 μL of water.



Digestion of Zero-Length Recordings

6 Digest zero-length recordings with  5 μL BclI in **50 μL** volume:

1h 20m



 37 °C  01:00:00

 65 °C  00:20:00

 On ice

Gel Separation and Extraction of Recordings

7 Perform gel electrophoresis in **4% LMP agarose** at **3 V/cm** for  01:00:00 .

1h

8 Extract the region from **300 bp to 500 bp** and elute in  12 μL of water.



Store gel products at  -20 °C in the post-PCR freezer.

Indexing PCR of Recombinant DNA

9 Using  1 μL of gel product, perform indexing PCR in  12 μL volume with **Nextera primers**.



9.1  98 °C  00:00:30

30s

9.2  98 °C  00:00:05

35s

 67 °C  00:00:10

 72 °C  00:00:20

9.3  go to step #9.2 9 X



9.4

72 °C

00:02:00

2m

SPRI Bead Purification

10 Purify amplicons by **1X SPRI** bead cleanup, eluting in 10 µL of water.



D1000 Tapescreen

11 Take 1 µL of each indexed amplicon to D1000 Tapescreen.



Expected result

Recordings should be above 158 + 67 (adaptors) + 69 (indexes) = **294 bp**.

12 Pool indexed amplicons according to desired relative molarities for sequencing.