

Jan 06, 2023

Hestrin-Schramm (HS) medium

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly7xyelx9/v1>



Andreas Sagen¹

¹University of Oslo



Andreas Sagen

University of Oslo

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly7xyelx9/v1>

Protocol Citation: Andreas Sagen 2023. Hestrin-Schramm (HS) medium. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.eq2ly7xyelx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 05, 2023

Last Modified: January 06, 2023

Protocol Integer ID: 74795

Keywords: Komagataeibacter, bacterial cellulose, cellulose, growth of komagataeibacter, komagataeibacter, yeast extract, schramm medium, medium hestrin, yeast, primary nitrogen source, bacteria, hestrin, schramm, common medium composition, usable energy in aerobic organism, nitrogen, magnesium sulphate, phosphate, medium composition, glucose, citrate, sodium phosphate, magnesium, intermediate compound, aerobic organism

Abstract

Hestrin-Schramm medium is a common medium composition to grow *Komagataeibacter*. The medium contains glucose as a source for energy and carbon, while yeast extract and peptone are the primary nitrogen sources. Both sodium phosphate and magnesium sulphate is added to benefit the bacteria, but it is not necessary to support growth of *Komagataeibacter*. Citrate is used as an intermediate compound in the tricarboxylic acid cycle and an essential element to produce usable energy in aerobic organisms.

Guidelines

Follow step by step, unless stated otherwise. Equipment needed should be standard to a microbiology lab.

Materials

Analytical scale, autoclave, bottle(s), weight vessel, LAF bench

Protocol materials

☒ Dextrose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021

☒ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

☒ Magnesium sulfate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M7506

☒ Citrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S4641

☒ Peptone Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905

☒ Disodium phosphate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9763

☒ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296



Safety warnings

! You can mix Dextrose from the beginning with the other compounds, and autoclave together. While this is more time efficient and easier, it is important to take into account the possibility of toxic byproducts produced by the Millard reaction when autoclaving, producing Acrylamide, a probable human carcinogen (IARC Group 2A).

Furthermore, when removing autoclaved components, be sure to take care as these can be very hot. If using antibiotics, use sufficient PPE to protect yourself, as some can be toxic to humans.

Before start

Prepare glassware by cleaning it, and ensure that scale is sufficiently calibrated



100 mL Dextrose solution

20m

1 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

1.1 Fill the bottle with  60 mL double-distilled water

1.2 Measure  10000 mg Dextrose

Powdered compounds:

 Dextrose Merck MilliporeSigma (Sigma-Aldrich) Catalog #G7021

1.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for  00:05:00

5m

1.4 Adjust pH while mixing to  5.7 using concentrated sodium hydroxide

1.5 Add distilled water to a total of  100 mL

1.6 Autoclave liquid at  121 °C for  00:15:00

15m

500 mL HS (broth) medium

20m

2 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

2.1 Fill the bottle with  300 mL double-distilled water

2.2 Measure  2500 mg Peptone,  2500 mg Yeast extract,  1350 mg Disodium phosphate,  750 mg Citrate and  600 mg Magnesium sulfate

Powdered compounds:



⊗ Peptone Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905

⊗ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

⊗ Disodium phosphate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9763

⊗ Citrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S4641

⊗ Magnesium sulfate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M7506

2.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

⌚ 00:05:00

5m

2.4 Adjust pH while mixing to pH 5.7 using concentrated citric acid

2.5 Autoclave liquid at $^{\circ}\text{C}$ 121 $^{\circ}\text{C}$ for ⌚ 00:15:00

15m

2.6 In a LAF bench, add flask 100 mL sterile Dextrose solution

2.7 Add sterile water to a total of flask 500 mL

Note

Cool to 50°C and supplement with antibiotics as appropriate

500 mL HS (agar) medium

20m

3 All compounds are measured using a high precision analytical scale from powdered compounds. Each compound is measured to within 1% of the target weight. All compounds are mixed in a Duran bottle

3.1 Fill the bottle with flask 300 mL double-distilled water

3.2 Measure flask 2500 mg Peptone, flask 2500 mg Yeast extract, flask 1350 mg Disodium phosphate, flask 750 mg Citrate, flask 600 mg Magnesium sulfate and flask 7500 mg agar



Powdered compounds:

⊗ Peptone Merck MilliporeSigma (Sigma-Aldrich) Catalog #P5905

⊗ Yeast Extract Merck MilliporeSigma (Sigma-Aldrich) Catalog #Y0875

⊗ Disodium phosphate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S9763

⊗ Citrate Merck MilliporeSigma (Sigma-Aldrich) Catalog #S4641

⊗ Magnesium sulfate Merck MilliporeSigma (Sigma-Aldrich) Catalog #M7506

⊗ Agar Merck MilliporeSigma (Sigma-Aldrich) Catalog #A1296

3.3 Add powdered solids into bottle, and use a magnetic mixer with a stir bar to mix for

⌚ 00:05:00

5m

3.4 Adjust pH while mixing to pH 5.7 using concentrated citric acid

3.5 Autoclave liquid at $^{\circ}\text{C}$ 121 $^{\circ}\text{C}$ for ⌚ 00:15:00

15m

3.6 In a LAF bench, add flask 100 mL sterile Dextrose solution

3.7 Add sterile water to a total of flask 500 mL

Note

Cool to 50°C and supplement with antibiotics as appropriate

Agar can be stored, then reheated to 50°C to be poured