

Aug 23, 2018

Hi-C library preparation for the *Lateolabrax maculatus* genome

 [GigaScience](#)

 In 4 collections

DOI

<https://dx.doi.org/10.17504/protocols.io.ss4eegw>

Chang Li¹

¹BGI-Qingdao, BGI-Shenzhen, Qingdao, 266555, China.

GigaScience Press

BGI



Chang Li

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ss4eegw>

External link: <https://doi.org/10.1093/gigascience/giy114>

Protocol Citation: Chang Li 2018. Hi-C library preparation for the *Lateolabrax maculatus* genome. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ss4eegw>

**Manuscript citation:**

hangwei Shao, Chang Li, Na Wang, Yating Qin, Wenteng Xu, Qun Liu, Qian Zhou, Yong Zhao, Xihong Li, Shanshan Liu, Xiaowu Chen, Shahid Mahboob, Xin Liu, Songlin Chen, Chromosome-level genome assembly of the spotted sea bass, *Lateolabrax maculatus*, *GigaScience*, Volume 7, Issue 11, November 2018, giy114, <https://doi.org/10.1093/gigascience/giy114>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: August 21, 2018

Last Modified: November 21, 2020

Protocol Integer ID: 14908

Keywords: lateolabrax maculatus genome this protocol, lateolabrax maculatus genome, lateolabrax maculatus, genome, process of hi

Abstract

This protocol is used to clarify the process of Hi-C library preparation for L.maculatus genome.



Materials

STEP MATERIALS

☒ PBS **Invitrogen - Thermo Fisher**

☒ 37% formaldehyde **Merck MilliporeSigma (Sigma-Aldrich)**

☒ glycine **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mM NaCl **Beyotime**

☒ 1× PMSF **Merck MilliporeSigma (Sigma-Aldrich)**

☒ enzyme (Mbo I) **New England Biolabs**

☒ 10 mM dCTP, 10 mM dGTP, 10 mM dTTP **Invitrogen - Thermo Fisher**

☒ 5U/ µl DNA Polymerase I, Large (Klenow) Fragment **New England Biolabs**

☒ biotinylated residue (0.4 mM biotin-14-dATP) **Invitrogen - Thermo Fisher**

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mg/ml BSA **New England Biolabs**

☒ T4 DNA ligase **New England Biolabs**

☒ 10 mg/ml proteinase K **New England Biolabs**

☒ 1% SDS **Ambion**

☒ MPS **Invitrogen - Thermo Fisher**

☒ 70% ethanol **Sinopharm Chemical Reagent Co.**

☒ Dynabeads MyOne Streptavidin T1 **Invitrogen - Thermo Fisher**

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ 10 mM ATP **New England Biolabs**

☒ 25 mM dNTP mix **Enzymatics**

☒ 10 U/µl NEB T4 PNK, 3 U/µl NEB T4 DNA polymerase I , 5 U/µl NEB DNA polymerase I, Large (Klenow) Fragment **New England Biolabs**

☒ 10X NEBuffer 2 **New England Biolabs**

☒ 10 mM dATP **Invitrogen - Thermo Fisher**

☒ 5 U/µl NEB Klenow exo minus **New England Biolabs**

☒ 10X T4 PNK Reaction Buffer **New England Biolabs**

☒ 100mM ATP **FERMENTAS Inc.**



⊠ 600 U/ul T4 DNA Ligase **New England Biolabs**

⊠ 50% PEG8000 **RIGAKU**

⊠ PBS **Invitrogen - Thermo Fisher**

⊠ 37% formaldehyde **Merck MilliporeSigma (Sigma-Aldrich)**

⊠ glycine **Merck MilliporeSigma (Sigma-Aldrich)**

⊠ 10 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**

⊠ 10 mM NaCl **Beyotime**

⊠ 1× PMSF **Merck MilliporeSigma (Sigma-Aldrich)**

⊠ enzyme (Mbo I) **New England Biolabs**

⊠ 10 mM dCTP, 10 mM dGTP, 10 mM dTTP **Invitrogen - Thermo Fisher**

⊠ 5U/ μl DNA Polymerase I, Large (Klenow) Fragment **New England Biolabs**

⊠ biotinylated residue (0.4 mM biotin-14-dATP) **Invitrogen - Thermo Fisher**

⊠ 10X NEB T4 DNA ligase buffer **New England Biolabs**

⊠ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich)**

⊠ 10 mg/ml BSA **New England Biolabs**

⊠ T4 DNA ligase **New England Biolabs**

⊠ 10 mg/ml proteinase K **New England Biolabs**

⊠ 1% SDS **Ambion**

⊠ MPS **Invitrogen - Thermo Fisher**

⊠ 70% ethanol **Sinopharm Chemical Reagent Co.**

⊠ Dynabeads MyOne Streptavidin T1 **Invitrogen - Thermo Fisher**

⊠ 10X NEB T4 DNA ligase buffer **New England Biolabs**

⊠ 10 mM ATP **New England Biolabs**

⊠ 25 mM dNTP mix **Enzymatics**

⊠ 10 U/μl NEB T4 PNK, 3 U/μl NEB T4 DNA polymerase I , 5 U/μl NEB DNA polymerase I, Large (Klenow) Fragment **New England Biolabs**

⊠ 10X NEBuffer 2 **New England Biolabs**

⊠ 10 mM dATP **Invitrogen - Thermo Fisher**

⊠ 5 U/μl NEB Klenow exo minus **New England Biolabs**

⊠ 10X T4 PNK Reaction Buffer **New England Biolabs**

⊠ 100mM ATP **FERMENTAS Inc.**



⌘ 600 U/ul T4 DNA Ligase **New England Biolabs**

⌘ 50% PEG8000 **RIGAKU**



Protocol materials

✕ 10 mM NaCl **Beyotime**

✕ PBS **Invitrogen - Thermo Fisher**

✕ 37% formaldehyde **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 10X NEBuffer 2 **New England Biolabs**

✕ 10 mM dATP **Invitrogen - Thermo Fisher**

✕ 5 U/μl NEB Klenow exo minus **New England Biolabs**

✕ 5U/ μl DNA Polymerase I, Large (Klenow) Fragment **New England Biolabs**

✕ 1% SDS **Ambion**

✕ MPS **Invitrogen - Thermo Fisher**

✕ 25 mM dNTP mix **Enzymatics**

✕ 600 U/ul T4 DNA Ligase **New England Biolabs**

✕ 10 mg/ml BSA **New England Biolabs**

✕ 1× PMSF **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 10 mg/ml proteinase K **New England Biolabs**

✕ 10X NEB T4 DNA ligase buffer **New England Biolabs**

✕ biotinylated residue (0.4 mM biotin-14-dATP) **Invitrogen - Thermo Fisher**

✕ 10 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 10 U/μl NEB T4 PNK, 3 U/μl NEB T4 DNA polymerase I , 5 U/μl NEB DNA polymerase I, Large (Klenow) Fragment **New England Biolabs**

✕ 100mM ATP **FERMENTAS Inc.**

✕ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich)**

✕ T4 DNA ligase **New England Biolabs**

✕ 50% PEG8000 **RIGAKU**

✕ 10 mM ATP **New England Biolabs**

✕ 70% ethanol **Sinopharm Chemical Reagent Co.**

✕ 10X NEB T4 DNA ligase buffer **New England Biolabs**

✕ enzyme (Mbo I) **New England Biolabs**

✕ glycine **Merck MilliporeSigma (Sigma-Aldrich)**

✕ 10 mM dCTP, 10 mM dGTP, 10 mM dTTP **Invitrogen - Thermo Fisher**

✕ Dynabeads MyOne Streptavidin T1 **Invitrogen - Thermo Fisher**



10X T4 PNK Reaction Buffer **New England Biolabs**



Sample preparation

- 1) Blood sample was centrifuged at 2500g for 5min at 4°C, collect precipitated cell;
2) Add 1ml PBS to resuspended blood cell.

☒ PBS **Invitrogen - Thermo Fisher**

Formaldehyde fixation

- 1) 37% formaldehyde was added with the final concentration of 1 % and the reaction was stopped with glycine after standing for 10 min at room temperature.
2) The formaldehyde fixed powder was then re-suspended in nuclei isolation buffer (10 mM Tris-HCl pH 8.0 , 10 mM NaCl, 1× PMSF).

☒ 37% formaldehyde **Merck MilliporeSigma (Sigma-Aldrich)**

☒ glycine **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mM Tris-HCl **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mM NaCl **Beyotime**

☒ 1× PMSF **Merck MilliporeSigma (Sigma-Aldrich)**

Enzyme digestion

- Add enzyme (Mbo I) to digest the DNA.

☒ enzyme (Mbo I) **New England Biolabs**

DNA fragment end reparation

- Add 10 mM dCTP, 10 mM dGTP, 10 mM dTTP, 5U/ μl DNA Polymerase I, Large (Klenow) Fragment using a biotinylated residue (0.4 mM biotin-14-dATP).

☒ 10 mM dCTP, 10 mM dGTP, 10 mM dTTP **Invitrogen - Thermo Fisher**

☒ 5U/ μl DNA Polymerase I, Large (Klenow) Fragment **New England Biolabs**

☒ biotinylated residue (0.4 mM biotin-14-dATP) **Invitrogen - Thermo Fisher**

In situ ligation

- Add 10X NEB T4 DNA ligase buffer , 10% Triton X-100, 10 mg/ml BSA , T4 DNA ligase.

☒ 10X NEB T4 DNA ligase buffer **New England Biolabs**

☒ 10% Triton X-100 **Merck MilliporeSigma (Sigma-Aldrich)**

☒ 10 mg/ml BSA **New England Biolabs**

 T4 DNA ligase **New England Biolabs**

Reverse-crosslinked

- 6 Add 10 mg/ml proteinase K and 1% SDS to the tube and incubate at 56°C for overnight and purified ,put the Reverse-crosslinked DNA liquid into three tube equally, add 1.5× volumes of AMPure XP mixture to each tube, vortex and spin down briefly, incubate for 10 min at room temperature, place on the MPS for 5 min at room temperature, discard supernatant, wash the beads twice with 1 ml of freshly made 70% ethanol, air-dry the beads completely and re-suspend the beads in 30 µl of ddH₂O.

 10 mg/ml proteinase K **New England Biolabs** 1% SDS **Ambion** MPS **Invitrogen - Thermo Fisher** 70% ethanol **Sinopharm Chemical Reagent Co.**

Biotin-containing fragments capture

- 7 Shear 20 mg of DNA and capturing the biotin-containing fragments on streptavidin-coated beads using Dynabeads MyOne Streptavidin T1.

 Dynabeads MyOne Streptavidin T1 **Invitrogen - Thermo Fisher**

DNA fragment end repairation

- 8 Add 10X NEB T4 DNA ligase buffer with 10 mM ATP, 25 mM dNTP mix, 10 U/µl NEB T4 PNK, 3 U/µl NEB T4 DNA polymerase I , 5 U/µl NEB DNA polymerase I, Large (Klenow) Fragment.

 10X NEB T4 DNA ligase buffer **New England Biolabs** 10 mM ATP **New England Biolabs** 25 mM dNTP mix **Enzymatics** 10 U/µl NEB T4 PNK, 3 U/µl NEB T4 DNA polymerase I , 5 U/µl NEB DNA polymerase I, Large (Klenow) Fragment **New England Biolabs**

DNA fragment adenylation

- 9 Add 10X NEBuffer 2, 10 mM dATP, 5 U/µl NEB Klenow exo minus.

 10X NEBuffer 2 **New England Biolabs** 10 mM dATP **Invitrogen - Thermo Fisher** 5 U/µl NEB Klenow exo minus **New England Biolabs**



Adaptor ligation

- 10 Add 10X T4 PNK Reaction Buffer, 100mM ATP, 600 U/ul T4 DNA Ligase, 50% PEG8000 , 50 uM Ad153 barcode oligo_2B mix (BGI, Shenzhen, China), and followed by PCR (95°C 3 min.; [98°C 20 sec., 60°C 15 sec., 72°C 15 sec.] (8 cycles); 72°C 10 min.

 10X T4 PNK Reaction Buffer **New England Biolabs**

 100mM ATP **FERMENTAS Inc.**

 600 U/ul T4 DNA Ligase **New England Biolabs**

 50% PEG8000 **RIGAKU**