

Aug 09, 2023

🌐 High molecular weight DNA extraction for marine macroalgal tissue

DOI

dx.doi.org/10.17504/protocols.io.14egn2dnpg5d/v1



Malia Moore¹, Taylor S. Steele¹

¹Scripps Institution of Oceanography



Malia Moore

Scripps Institution of Oceanography, Salk Institute for Biol...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn2dnpg5d/v1>

Protocol Citation: Malia Moore, Taylor S. Steele 2023. High molecular weight DNA extraction for marine macroalgal tissue. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.14egn2dnpg5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 08, 2023

Last Modified: August 09, 2023

Protocol Integer ID: 81551

Keywords: Lyophilizing algal tissue, DNA extraction, Lysis and first precipitation, Final precipitation, Column cleanup, oxford nanopore hmw dna extraction from arabidopsis leaf, marine macroalgal tissue, oxford nanopore hmw dna extraction, high molecular weight dna extraction, cell culture dna midi kit for column cleanup, dna extraction, marine macroalgae, cell culture dna midi kit, sulfated polysaccharide, dna product, polysaccharide, arabidopsis leaf, polyphenolic contamination, dna, polyphenolic

Abstract

This protocol details high molecular weight DNA extraction for marine macroalgal tissue. Marine macroalgae contain a variety of unique cell wall components including sulfated polysaccharides and polyphenolics. These components often co-elute with high molecular weight (HMW) DNA and lead to reduced library prep and sequencing outcomes. This protocol incorporates polyvinylpolypyrrolidone (PVPP) and β -mercaptoethanol (BME) to reduce polyphenolic contamination, and an early salting out step with potassium acetate (KOAc) to address polysaccharides. This protocol is largely adapted from an Oxford Nanopore HMW DNA extraction from Arabidopsis leaves, which incorporates the QIAGEN Blood and Cell Culture DNA Midi Kit for column cleanup. The DNA product often requires additional cleanup after elution, and we suggest the BluePippin 15kb size selection for all HMW applications.

Attachments



[711-1533.pdf](#)



55KB



Guidelines

Marine macroalgae contain a variety of unique cell wall components including sulfated polysaccharides and polyphenolics. These components often co-elute with high molecular weight (HMW) DNA and lead to reduced library prep and sequencing outcomes. This protocol incorporates polyvinylpolypyrrolidone (PVPP) and β -mercaptoethanol (BME) to reduce polyphenolic contamination, and an early salting out step with potassium acetate (KOAc) to address polysaccharides.¹ This protocol is largely adapted from an Oxford Nanopore HMW DNA extraction from *Arabidopsis* leaves, which incorporates the QIAGEN Blood and Cell Culture DNA Midi Kit for column cleanup.² The DNA product often requires additional cleanup after elution, and we suggest the BluePippin 15kb size selection for all HMW applications.


Additional tips:

- In the field or in lab, it is vital to scrape off all surface epiphytes and wash the sample in clean water before flash freezing to reduce contaminants common in the marine environment that confound genome assembly.
- Marine macroalgae are incredibly diverse in biochemical content, so individual seaweeds may require troubleshooting. Suggested alterations include varying input tissue type or quantity, increasing CTAB or BME percent, or adding a second chloroform separation.
- It may be necessary to carry out extractions of the same tissue in parallel to yield sufficient DNA, especially when large losses from BluePippin are expected. It is not suggested to combine multiple extractions onto the same column, as this may lead to overloading and a dirty sample. This protocol as written, paired with BluePippin, has produced sequencing-quality DNA for Nanopore from a red alga *Porteria hornemanii* and a brown alga *Macrocystis pyrifera*. For *P. hornemanii*, a single  20 mL extraction produced sufficient DNA for sequencing, but for *M. pyrifera*, three parallel extractions of  20 mL were necessary.












Materials

Equipment:

- Lyophilizer
- Stir plate
- Heat block or water bath
- Vortex
- Mortar and pestle
- Refrigerated centrifuge for spins up to 3,500 xg with  50 mL
- Suggested: Sage Science BluePippin

Consumables:

- Stock solution:  1 Molarity (M) Tris-HCl,  9.5
- Stock solution:  5 Molarity (M) sodium chloride (NaCl)
- Stock solution:  500 millimolar (mM) ethylenediaminetetraacetic acid (EDTA)
- Stock solution:  5 Molarity (M) potassium acetate (KOAc)
- Cetyltrimethylammonium bromide (CTAB)
- Polyethylene glycol (PEG) 8000
- β -mercaptoethanol (BME)
- Polyvinylpyrrolidone (PVPP)
- RNase A,  100 mg/mL (eg. QIAGEN Mat. #1007885)
- 100% isopropanol
- 95-100% ethanol
- Nuclease-free water
-  Blood & Cell Culture DNA Mini Kit (25) **Qiagen Catalog #13323**
- Tris-EDTA (TE) buffer
-  50 mL Falcon Tubes
-  DNA LoBind Tube 1.5ml **Eppendorf Catalog #022431021**
- Suggested: Sage Science High Pass Plus Cassette (BPLUS10 or BPLUS03) for BluePippin

Troubleshooting



Lyophilizing algal tissue

- 1 Flash-freeze algal tissue in liquid nitrogen (target \geq 5 g wet tissue).
- 2 Quickly transfer sample to lyophilization container and freeze dry for 36-48 hours.
- 3 Macerate the tissue with a clean spatula to increase surface area and put on the lyophilizer for another 24:00:00 . 1d
- 4 Remove and refrigerate with desiccant for immediate use, or store at -80 °C for longer periods.

Setting up the DNA extraction

- 5 Prepare desired volume of Carlson lysis buffer (100 millimolar (mM) Tris-HCl, 9.5 , 2% CTAB, 1.4 Molarity (M) NaCl, 1% PEG 8000, 20 millimolar (mM) EDTA) and mix Overnight on a magnetic stirrer. The stock solutions suggested under consumables will yield a homogenous buffer with no precipitate. 1d
- 6 Pre-heat a heat block or water bath to 65 °C and place in a fume hood.
- 7 For each extraction, transfer 20 mL of Carlson lysis buffer to a 50-ml Falcon tube.
- 8 In a fume hood, add 400 μ L BME (originally 50 μ L) and mix by vortexing. Pre-warm the solution to 65 °C for 00:30:00 before starting the extraction. 30m
- 9 Scoop 0.5 teaspoons lyophilized plant tissue into a clean mortar and add 50-100 mg powdered PVPP. Grind with pestle for \sim 00:00:30 , until tissue is powdered and combined, but not long enough to introduce significant moisture. Move immediately into DNA extraction. 30s

Lysis and first precipitation





10 Pour tissue into the warm buffer. Invert 5 times.

11 Add  40 μL of RNase A and vortex for  00:00:05 .

5s



12 **Optional:** If using a heat block with mixing, set the block (still at  65 $^{\circ}\text{C}$) to mixing at  300 rpm, 00:05:00 .



13 Incubate for  01:00:00 at  65 $^{\circ}\text{C}$.

1h



13.1 Invert 10 times every 15 minutes.

13.2 At 30 minutes, add another  40 μL of RNase A, inverting 10 times to combine.



14 Allow the tubes to cool down to  Room temperature for  00:10:00 .

10m

15 Add  20 mL chloroform and vortex for two pulses of  00:00:05 each.

5s





16 Centrifuge the tubes at  3500 x g, 4 $^{\circ}\text{C}$, 00:15:00 .

15m








17 In a fume hood, transfer the top layer of lysate from each tube to a new 50-ml Falcon tube, without disturbing the interphase.


Note


Tip: The lysate layer should be  14–18 mL of solution, but it is recommended to use widebore tips, transferring  1 mL at a time. Tips can also be widened by cutting standard P1000 tips.



- 18 Mix supernatant with 0.4X 1M 5 Molarity (M) potassium acetate (KOAc) at Room temperature, inverting at least 10 times to combine, then incubate On ice for 00:20:00. 20m  
- 19 Centrifuge the tubes at 3500 x g, 4°C, 00:45:00. 45m 
- 20 Remove and retain the supernatant.
- Note














Tip: This may best be done by pouring slowly and observing the polysaccharide-salt pellet, which may be mobile. Leave some liquid behind in favor of avoiding the pellet.
- 21 Add 0.7X volumes of isopropanol. Invert 10 times. Incubate at -80 °C for 00:15:00. 15m  
- Note

Do not extend this incubation.
- 22 Centrifuge the sample at 3500 x g, 4°C, 00:45:00. 45m 
- Note

Tip: If available, a fixed-angle centrifuge will create a pellet on the wall of the tube that has greater surface area for dissolution in step 24 (as compared to a conical pellet at the base of a falcon tube from a swinging bucket).
- 23 Discard the supernatant without disturbing the pellet. Use sterile wipes to absorb the liquid on the tube walls, being careful not to disturb the pellet.
- 24 To each pellet, add 10 mL G2 buffer, from the QIAGEN kit. Incubate at 50 °C for 30-60 minutes, or until the pellet is dissolved. Swirl the pellet to mix but do not try to pipette or vortex. 








Column cleanup

- 25 Equilibrate a QIAGEN Genomic-tip 100/G column with  4 mL of Buffer QBT.
- 26 Pour the DNA in G2 buffer through the equilibrated column and allow it to flow through with just gravity.
- 27 Once all the lysate has passed through, wash the column with  8 mL of Buffer QC. 
- 28 Repeat the wash with another  8 mL of Buffer QC. 
- 29 Place the column over a clean 50-mL Falcon tube, and elute the genomic DNA with  5 mL of Buffer QF, pre-warmed to  55 °C .
- 30 Allow the eluate to cool down to  Room temperature .
- 31 Add  3.5 mL of isopropanol to the eluted DNA and mix by inverting the tube 10 times. 
- 32 Incubate the tube at  -20 °C for at least 3 hours, or  Overnight .  15m



Final precipitation

1h 10m

- 33 Centrifuge the tube at  3500 x g, 4°C, 00:45:00 .  45m
- 34 Discard the supernatant without disturbing the pellet. 
- 35 Add  4 mL of ice-cold 70% ethanol to the pelleted DNA and invert the tube 10 times. 



36 Centrifuge at  3500 x g, 4°C, 00:10:00 .




10m

Note

Tip: If using a swinging bucket centrifuge the DNA will pellet at the base of the tube and be easy to locate and resuspend. If using a fixed angle, mark the side of the tube that faces outwards in order to locate the pellet for washes and elution.




37 Discard the supernatant without disturbing the pellet. Use sterile wipes to dry the tube walls, being careful not to disturb the pellet.


38 Resuspend the DNA in  100 µL of TE buffer and incubate at  Room temperature , typically  Overnight .

15m



39 Transfer the DNA into a nuclease-free 1.5-mL tube (DNA LoBind tube preferred) using a wide-bore tip, and store at  4 °C .

Note

Tip: Often, waiting a further  48:00:00 before quantifying on Nanodrop and Qubit will allow the DNA to further relax and yield the most accurate results

40 Carry samples forward to BluePippin size selection, if available. This gel separation will retain DNA fragments greater than 15 kb and discard any residual contamination still evident on a Nanodrop trace. For these benefits, expect 50-70% loss of DNA.

Protocol references

Citations

1. Chekan, J. R. et al. Scalable Biosynthesis of the Seaweed Neurochemical, Kainic Acid. Angew Chem Int Ed Engl 58, 8454–8457 (2019).
2. Nanopore, Arabidopsis Leaf gDNA. https://community.nanoporetech.com/extraction_method_groups/plant-leaf-gDNA