


Mar 28, 2025

🌐 High molecular weight DNA extraction from algae : SoCK (Sorbitol + CTAB + KAc)

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlkdbzwg5r/v1>



Benoît acherie¹, Karine Labadie²

¹Genoscope/CEA; ²Genoscope

Benoit

HPM/Genoscope



Benoît Vacherie

Genoscope/CEA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkdbzwg5r/v1>

Protocol Citation: Benoît acherie, Karine Labadie 2025. High molecular weight DNA extraction from algae : SoCK (Sorbitol + CTAB + KAc). **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkdbzwg5r/v1>



License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: March 25, 2025

Last Modified: March 28, 2025

Protocol Integer ID: 124949

Keywords: extraction, DNA, algae, HMW, Long read sequencing, high molecular weight, high molecular weight dna extraction protocol from algae, high molecular weight dna extraction from algae, algae with sorbitol, high molecular weight dna extraction protocol, high molecular weight dna extraction, dna extraction, algae, flash frozen algal fragment, types of algae, frozen algal fragment, conventional extraction with ctab, algal fragment, conventional extraction, extraction, dna, purification, sorbitol, purification with potassium acetate, molecular weight dna, weight dna, high molecular weight

Abstract

High Molecular Weight DNA extraction protocol from algae for Long Read sequencing.

Extraction is performed from flash frozen algal fragments stored at -80° C.

The protocol is adapted for an extraction of 1 g of algal fragment.

The protocol is based on washing the algae with Sorbitol, followed by a conventional extraction with CTAB and then a purification with Potassium acetate coupled with chloroform.

This protocol is specifically adapted to all types of algae (green, red or brown).

Guidelines

Use only wide bore tips.

Work in a chemical hood when using 2-mercaptoethanol.

Allow the DNA to resuspend for at least of 24 hours before proceeding with QC.

Materials

Reagents :

 Tris HCl Buffer 1M Solution, Sterile pH 8.0 **Bio Basic Inc. Catalog #SD8127.SIZE.450ml**

 EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**

 CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**

 Sodium chloride **P212121**

 PEG-8000

 Sorbitol **P212121**

 Polyvinylpyrrolidone (PVP-40) **Bio Basic Inc. Catalog #PB0436.SIZE.250g**

 1X TE buffer (10 mM Tris-HCl pH 8.0 1 mM EDTA)  70% Ethanol

 Sodium acetate **P212121**

 2-Propanol (IsoPropanol) **Bio Basic Inc. Catalog #PC8601.SIZE.4L**

 Chloroform Isoamyl alcohol 24:1 **Merck**

 Potassium Acetate

 Ethanol absolute

 RNaseA **Thermo Fisher Catalog #12091021**

 2-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**

 Liquid Nitrogen

Equipment :

 Eppendorf ThermoMixer C **pipette.com Catalog #2231000667**

 Mortar and Pestels

 Vortex Mixer

Chemical fume hood



Protocol materials

- ✕ Liquid Nitrogen
- ✕ 2-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**
- ✕ Ethanol absolute
- ✕ Potassium Acetate
- ✕ Chloroform Isoamyl alcohol 24:1 **Merck**
- ✕ RNaseA **Thermo Fisher Catalog #12091021**
- ✕ 2-Propanol (IsoPropanol) **Bio Basic Inc. Catalog #PC8601.SIZE.4L**
- ✕ Sodium acetate **P212121**
- ✕ 70% Ethanol
- ✕ 1X TE buffer (10 mM Tris-HCl pH 8.0 1 mM EDTA)
- ✕ Eppendorf ThermoMixer C **pipette.com Catalog #2231000667**
- ✕ Vortex Mixer
- ✕ Polyvinylpyrrolidone (PVP-40) **Bio Basic Inc. Catalog #PB0436.SIZE.250g**
- ✕ EDTA (0.5 M), pH 8.0 **Life Technologies Catalog #AM9260G**
- ✕ Tris HCl Buffer 1M Solution, Sterile pH 8.0 **Bio Basic Inc. Catalog #SD8127.SIZE.450ml**
- ✕ CTAB (Hexadecyltrimethylamm onium bromide) **BBI Biotech Catalog #CB0108-100g**
- ✕ Mortar and Pestels
- ✕ Sodium chloride **P212121**
- ✕ PEG-8000
- ✕ Sorbitol **P212121**



Preparation of reagents

- 1 **Sorbitol pre-wash buffer** : the buffer can be stored for 6 months at 4°C

	Reagents	Final Concentration
	Tris HCl ph8	100 mM
	EDTA	5 mM
	Sorbitol	700 mM
	PVP 40	1%
	H2O	Qsp 100 ml

- 2 **Lysis Buffer**

	Reagents	Final Concentration
	Tris HCl ph8	100 mM
	EDTA	20 mM
	CTAB	2 %
	NaCl	1.4 M
	PEG 8000	1 %
	H2O	Qsp 20 ml

Sorbitol Wash

39m

- 3 Add **20 ml of Sorbitol pre-wash Buffer** in a 50 ml tube.

 20 mL Sorbitol buffer

- 4 Put a **mortar in ice**.
Cool the mortar and pestle with **liquid nitrogen** until the bubbling stops.

10m

 Liquid Nitrogen



- 5 **Grind 1g of frozen sample** to a fine powder (approx. 2 min), without adding liquid nitrogen.



1 g



00:02:00

2m

- 6 **Transfer powder** to sorbitol pre-wash buffer

- 7 Add **20 µl of 2-Mercaptoethanol** and vortex 5 s



2-mercaptoethanol Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250

- 8 **Centrifuge** 5 minutes with low deceleration



5000 x g, Room temperature, 00:05:00 , Acc 9 / Dec 7

5m

- 9 Carrefully **remove the supernatant** without resuspending the pellet.

2m

- 10 **Perform a second sorbitol wash** as above (step 3 + 7 + 8 +9)

20m

DNA extraction and purification

5h 30m

- 11 Resuspend the pellet with **20 ml of Lysis buffer**.

- 12 Add **50 µl of 2-Mercaptoethanol** and **40 µl of RNaseA** (100mg/ml) and vortex 5 s



2-mercaptoethanol Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250



RNaseA Thermo Fisher Catalog #12091021

- 13 **Incubate 1h at 65°C** with low agitation (300 rpm)



01:00:00



65 °C

1h

- 14 **Centrifuge** 10 minutes with low deceleration



8500 x g, 4°C, 00:10:00 , Acc 9 / Dec 7

10m

- 15 Gently **transfer the supernatant** to a new 50 ml tube.



- 16 Add :
- 0.1 vol of **absolute ethanol**  Ethanol absolute
- 0.25 vol of **KAc 3M**  Potassium Acetate
- 1 vol of **chloroform isoamyl alcohol**  Chloroform Isoamyl alcohol 24:1 **Merck**
- 17 Vortex 5 sec and **incubate 20 min at -20°C** 20m
 00:20:00  -20 °C
- 18 **Centrifuge** 20 minutes with low deceleration 20m
 8500 x g, 4°C, 00:20:00 , Acc 9 / Dec 7
- 19 Gently **transfer the aqueous phase** to a new 50 ml tube.
- 20 Add :
- 0.8 vol of **isopropanol**
 2-Propanol (IsoPropanol) Bio Basic Inc. Catalog #PC8601.SIZE.4L
- 0.1 vol of **NaAc 3M**  Sodium acetate **P212121**
- 0.2% of **2-Mercaptoethanol**
 2-mercaptoethanol **Merck MilliporeSigma (Sigma-Aldrich) Catalog #M6250**
- 21 mix gently by inversion 10 times and **incubate 1h at -20°C** 1h
 01:00:00  -20 °C
- 22 **Centrifuge** 20 minutes with low deceleration 20m
 8500 x g, 4°C, 00:20:00 , Acc 9 / Dec 7
- 23 Carefully **discard the supernatant, gently resuspend the pellet** with 1 ml of cold 70% ethanol.
 70% Ethanol
- 24 **Transfer the resuspended DNA** to a 1.5 ml DNA LoBind tube.



25 **Centrifuge** 10 minutes with low deceleration

10m

 5000 x g, 4°C, 00:10:00 , Acc 9 / Dec 7

26 Carefully **discard the supernatant**.

Air dry the pellet at RT for about 10 minutes.  00:10:00

10m

27 **Resuspend the pellet** with 50-100 µl of 1X TE buffer.

 1X TE buffer (10 mM Tris-HCl pH 8.0 1 mM EDTA)

28 **Allow the DNA to resuspend** for 2 hours at 55°C or ON at RT.

 02:00:00  55 °C

2h

29 **Store the DNA at 4°C.**

Sample quality control

30 Quantify your sample with a Qubit DNA HS assay.

Check the purity of the sample with a Nanodrop (measurements of 260/280 and 260/230 absorbance ratios).






Estimate the molecular weight of the sample with a Tapestation and/or a Femto pulse

31 Depending on the DNA concentrations and DNA length profiles, deplete short DNA molecules using SRE size selection kits (SRE XS, SRE or SRE XL kits).

Results

32 **QC results** obtained on different species. (3 green algae and 3 red algae)



Phylum	Species	DNA quantity (µg)	Extraction Yield (µg DNA/g organism)	Femto pulse kb	Profile Femto pulse
Chlorophyta	Codium saccatum	4,5	2,3	60	
	Codium sp. 1	9,5	9,5	56	
	Chaetomorpha linum	13	11,3	50	
Rhodophyta	Gracilaria dura	2,6	2,6	61	
	Phymatolithon calcareum	4,8	4,6	20	
	Griffithsia corallinoides	1,8	1,6	74	