



Mar 21, 2018

# High Salt Nuclear Extract Preparation

 [PLOS Pathogens](#)

DOI

[dx.doi.org/10.17504/protocols.io.kh2ct8e](https://dx.doi.org/10.17504/protocols.io.kh2ct8e)

Xavier Contreras<sup>1</sup>, Rosemary Kiernan<sup>1</sup>

<sup>1</sup>Institut de Genetique Humaine, UMR9002, Montpellier, France



**Rosemary Kiernan**

Institut de Genetique Humaine, Montpellier, France

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.kh2ct8e](https://dx.doi.org/10.17504/protocols.io.kh2ct8e)

External link: <https://doi.org/10.1371/journal.ppat.1006950>

**Protocol Citation:** Xavier Contreras, Rosemary Kiernan 2018. High Salt Nuclear Extract Preparation. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.kh2ct8e>

**Manuscript citation:**

Contreras X, Salifou K, Sanchez G, Helsmoortel M, Beyne E, Bluy L, Pelletier S, Rousset E, Rouquier S, Kiernan R, Nuclear RNA surveillance complexes silence HIV-1 transcription. PLoS Pathogens 14(3). doi: [10.1371/journal.ppat.1006950](https://doi.org/10.1371/journal.ppat.1006950)

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** October 31, 2017

**Last Modified:** March 24, 2018

**Protocol Integer ID:** 8474

**Keywords:** nuclear extract, Dignam, transcription



## Abstract

Adapted from Dignam JD, Lebovitz RM, Roeder RG. 1983. **Accurate transcription initiation by RNA polymerase II in a soluble extract from isolated mammalian nuclei.** Nucleic Acids Res. 11(5):1475-89. PMID:6828386

## Guidelines

The protocol is suitable for HeLa S3 cells. Parameters should be adapted for other cell types.

## Before start

Prepare the following buffers and store at +4°C:

<b>PMSF 0.2 M:</b> 3.48 g /100 ml DMSO			
<b>TGME:</b> Glycerol Tris, pH 7.9 MgCl <sub>2</sub> EDTA water	<b>Stock</b> 100 % 0.5 M 1 M 0.5 M	<b>Final Conc</b> 50 % 50 mM 5 mM 0.1 mM	<b>Vol (for 1 l)</b> 500 ml 100 ml 5 ml 0.2 ml to 1 l
<b>Low Salt Buffer:</b> Tris, pH 7.3 Glycerol MgCl <sub>2</sub> EDTA pH 8 water	<b>Stock</b> 1 M 50 % 1 M 0.5 M	<b>Final Conc</b> 20 mM 12.5 % 1.5 mM 0.2 mM	<b>Vol (for 1 l)</b> 20 ml 250 ml 1.5 ml 0.4 ml to 1 l
<b>High Salt Buffer:</b> Tris, pH 7.3 Glycerol MgCl <sub>2</sub> EDTA pH 8 KCl water	<b>Stock</b> 1 M 50 % 1 M 0.5 M 3 M	<b>Final Conc</b> 20 mM 12.5 % 1.5 mM 0.2 mM 1.2 M	<b>Vol (for 1 l)</b> 20 ml 250 ml 1.5 ml 0.4 ml 400 ml to 1 l
<b>Hypotonic Buffer (HB):</b> Tris, pH 7.3 KCl MgCl <sub>2</sub> water	<b>Stock</b> 1 M 3 M 1 M	<b>Final Conc</b> 10 mM 10 mM 1.5 mM	<b>Vol (for 1 l)</b> 10 ml 3.34 ml 1.5 ml to 1 l
<b>10 x Buffer</b> Tris, pH 7.3 KCl MgCl <sub>2</sub> water	<b>Stock</b> 1 M 3 M 1 M	<b>Final Conc</b> 30 mM 140 mM 3 mM	<b>Vol (for 1 l)</b> 30 ml 46.6 ml 3 ml to 1 l

Cold room set up:

thaw 0.2 M PMSF

prepare hypotonic buffer (500 ml HB + 0.35 ml beta-ME + 0.5 ml 0.2 M PMSF)

prepare freezing solution (20 ml TGME + 0.4 ml DTT, in a 50 ml tube)

prepare tubes

## Preparation of Cells

- 1 (Note: keep tubes on ice between centrifugations)  
Pellet cells at 1500 x g for 8 mins in appropriate tube (eg 50 ml tubes)  
Remove supernatant and wash cells once in PBS  
Centrifuge at 1000 x g for 10 mins  
Remove all supernatant and record the packed cell volume (PCV) in each tube  
Fill tube with HB and carefully resuspend cells  
Spin at 1000 x g for 5 mins  
Verify that the swollen cell volume (SCV) is greater than PCV  
Carefully pour off supernatant (pellet may be loose)  
Add 2 x PCV volume of HB and carefully resuspend the cells  
Place on ice for 10 mins

## Cell homogenization

- 2 homogenize cells by douncing 15 times (for HeLa S3, to be determined for other cell types), check for lysis by microscopy  
spin homogenized cells at 2600 x g for 15 mins  
record nuclear pellet volume (NPV)  
remove supernatant (cytoplasm)

## Nuclear Extraction

- 3 calculate NPV/2 and add the amount of low salt buffer and high salt buffer to separate tubes  
to each buffer add 0.0007 x NPV/2 of beta-ME, 0.01 x NPV/2 of 0.2 M PMSF, mix well  
resuspend nuclear pellets in ½ low salt buffer and dounce 6 times  
add the high salt buffer while mixing continuously  
mix for 30 mins after the addition of all of the high salt buffer  
spin at 15000 x g for 30 min  
recover supernatant  
for storage, snap freeze in liquid N<sub>2</sub> and keep at -80°C