


Jan 18, 2020

# 🌐 High-Throughput Beta-glucuronidase (GUS) assay for *Phaeodactylum tricornutum*

DOI

[dx.doi.org/10.17504/protocols.io.bbexijfn](https://dx.doi.org/10.17504/protocols.io.bbexijfn)



Erin Garza<sup>1</sup>, Vincent A Bielinski<sup>2</sup>

<sup>1</sup>J. Craig Venter Institute; <sup>2</sup>J. Craig Venter Institute, Synthetic Biology & Bioenergy Group

Protist Research to Opti...

JCVI West Protocols



Erin Garza

J. Craig Venter Institute

OPEN  ACCESS



DOI: [dx.doi.org/10.17504/protocols.io.bbexijfn](https://dx.doi.org/10.17504/protocols.io.bbexijfn)

**Protocol Citation:** Erin Garza, Vincent A Bielinski 2020. High-Throughput Beta-glucuronidase (GUS) assay for *Phaeodactylum tricornutum*. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.bbexijfn>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** January 16, 2020

**Last Modified:** January 18, 2020

**Protocol Integer ID:** 31927

**Keywords:** Gus assay, diatoms, *Phaeodactylum tricornutum*, high-throughput



## Abstract

A high-throughput method for measuring  $\beta$ -glucuronidase (GUS) activity in the diatom *Phaeodactylum tricornutum*. This protocol has been optimized for 250  $\mu$ l volumes. For larger volumes see the following protocol [dx.doi.org/10.17504/protocols.io.hefb3bn](https://doi.org/10.17504/protocols.io.hefb3bn), which this protocol was based off of.

## Materials

### MATERIALS

⊗ MUG **Gold Biotechnology Catalog #MUG**

⊗ Sodium carbonate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #222321**

⊗ B-PER<sup>®</sup>; Bacterial Protein Extraction Reagent **Thermo Fisher Catalog #78243**

Flat bottom transparent and opaque 96-well plates

GUS extraction buffer- 50 mM NaPO<sub>4</sub>H<sub>2</sub> (pH 7), 0.1% Triton X-100, and 10  $\mu$ M  $\beta$ ME + 1 mM 4-Methylumbelliferyl

$\beta$ -D-Glucuronide (MUG)

GUS stop buffer- 0.2 M Na<sub>2</sub>CO<sub>3</sub>

Plate reader

Swing bucket centrifuge with plate adapter

## Before start

*Phaeodactylum tricornutum* cultures were initially grown in 5 ml L1 + antibiotics in a 50-ml conical at 18°C until the cell concentration reached at least  $1 \times 10^6$  cells ml<sup>-1</sup>.



## Centrifuge

- 1 Transfer 250 µl of each *P. tricornutum* culture to a 96-well plate and centrifuge at 3000 x g for 10 min. Discard supernatant.

## Lyse

- 2 To lyse the cells, add 150 µl bacterial protein extraction reagent (B-PER, ThermoFisher) to each well and mix by pipetting.

## Centrifuge

- 3 Centrifuge plate for 10 min at 3000 x g. Transfer supernatants to a new 96-well plate, being careful not to disturb the cell debris.

## Extract

- 4 Transfer 50 µl of each lysate to a new plate and add 125 µl GUS extraction buffer + 1 mM MUG to each well. Incubate the plate for 1 h at 37°C.

GUS extraction buffer= 50 mM NaPO<sub>4</sub>H<sub>2</sub> (pH 7), 0.1% Triton X-100, and 10 µM βME

## Stop Reaction

- 5 To stop the reaction, add 150 µl GUS stop buffer (0.2 M Na<sub>2</sub>CO<sub>3</sub>) to each well and mix by pipetting.
- 6 Transfer 200 µl quenched reaction to an opaque 96-well plate.

## Read fluorescence

- 7 Determine fluorescence using a plate reader. Settings: excitation- 360 nm; emission- 440 nm.
- 8 If fluorescence readings are too high to get a readout, dilute with additional stop buffer.



## Normalization

---

- 9 Use remaining cell lysates to perform a BCA assay (or an equivalent assay) to normalize the GUS activity to total cell protein for each culture.