

Jun 19, 2025

🌐 High-throughput quantitative analysis of isoprenol (IP) by using an Agilent RapidFire-QqQ system

🔗 Forked from [High-throughput quantitative analysis of isoprenyl acetate \(IPA\) by using an Agilent RapidFire-QqQ system](#)



DOI

dx.doi.org/10.17504/protocols.io.5jyl8d6w6g2w/v1

Yan Chen¹, David Carruthers¹, Blake Simmons¹, Taek Soon Lee¹, Christopher Petzold¹

¹Lawrence Berkeley National Laboratory

LBNL omics



Christopher J Petzold

Lawrence Berkeley National Laboratory

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.5jyl8d6w6g2w/v1>

Protocol Citation: Yan Chen, David Carruthers, Blake Simmons, Taek Soon Lee, Christopher Petzold 2025. High-throughput quantitative analysis of isoprenol (IP) by using an Agilent RapidFire-QqQ system. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.5jyl8d6w6g2w/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 07, 2025

Last Modified: June 19, 2025

Protocol Integer ID: 117828

Keywords: Metabolites, RapidFire, High-throughput, Quantitative, Isoprenol, throughput quantitative analysis of isoprenol, throughput measurement of isoprenol, isoprenol, analytical assay, such as microbial fermentation culture, microbial fermentation culture, agilent quantitative analysis software, throughput screening of thousands sample, qq system this protocol details step, agilent 6460 qq system, throughput quantitation of ip

Funders Acknowledgements:

Department of Defense (DOD) Tri-Service Biotechnology for a Resilient Supply Chain (T-BRSC)

Grant ID: AWD00007196

Department of Energy (DOE) Office of Science

Grant ID: DE-AC02-05CH11231

Disclaimer

This protocol is for research purposes only.

Abstract

This protocol details steps in high-throughput measurement of isoprenol (IP) with the Agilent RapidFire 400 coupled with Agilent 6460 QqQ system, and subsequent data processing using the Agilent Quantitative analysis software. This analytical assay provides accurate and high-throughput quantitation of IP from a broad variety of sample types, such as microbial fermentation cultures, chemical overlays, intracellular cell lysate, etc. The protocol provides flexible options for measurements of small number of samples or ultra high-throughput screening of thousands samples per day.

Materials

⊗ Acetonitrile LCMS quality **JT Baker Catalog #9829-02**

⊗ Water LC-MS grade B&J Brand **VWR International (Avantor) Catalog #BJLC365-2.5**

⊗ Cis-3-Hexenyl Acetate **Merck MilliporeSigma (Sigma-Aldrich) Catalog #W317101**

⊗ 3-METHYL-3-BUTEN-1-YL ACETATE **Merck MilliporeSigma (Sigma-Aldrich) Catalog #S881171**

- Bio-Rad hard shell 96-well PCR plate (Bio-Rad part number HSP9601) or other equivalent RapidFire compatible 96-Well plate.
- LC compatible non adhesive alumni seal (Agilent part number: 06644-001)
- PlateLoc Thermal Microplate Sealer or equivalent type
- BioMek FX or any equivalent liquid transfer automation platform

Equipment

RapidFire 400

NAME

High-throughput liquid chromatography

TYPE

Agilent

BRAND

G9532AA

SKU

<https://www.agilent.com/en/product/liquid-chromatography-mass-spectrometry-lc-ms/lc-ms-instruments/high-throughput-lc-ms/rapidfire-400>

LINK



Troubleshooting

Safety warnings

- ❗
 - Wear proper PPE (gloves, safety goggles, and lab coat), and prepare solvents in a chemical fume hood.
 - Store organic solvents in a flammable solvent storage cabinet when not in use.
 - Always prepare necessary control samples to validate sample preparation procedures.
 - Dispose of used solvents and samples in accordance with your institute's Environmental, Health, and Safety requirements.

Before start

Refill RapidFire pumps with sufficient volumes of the following solvents to complete the analysis:

- LCMS grade water for pump 1.
- LCMS grade water for pump 2.
- LCMS grade methanol for pump 3.

Verify MS status:

- Perform Check-tune of QqQ system to ensure it passes calibration criteria. Service the QqQ and perform autotune if check-tune fails.
- Perform a quality control (QC) run to ensure QC compound peak full width at half maximum (FWHM), shape, and intensity are within quality standard range.



Isoprenol (IP) standard preparation

- 1 Prepare internal standard (IS) stock of Cis-3-Hexenyl Acetate (CAS: 3681-71-8) in methanol at 10 millimolar (mM) .
- 2 Prepare general diluent solution of 50 micromolar (μM) IS in 5 % volume Methanol (MeOH). This general diluent is used for diluting samples and preparing calibration standard diluent.
- 3 Prepare sample matrix that contains everything that is present in the typical sample except IP.

Note

A typical sample suitable for this protocol could be but not limited to microbial fermentation broth, IP trapping overlay such as dodecane, or extracted intracellular lysate of a microbial culture.

- 4 Prepare the calibration standard diluent (diluent used for generating calibration standards at a serial of concentrations) by spiking the same volume ratio of sample matrix into general diluent solution used to prepare test samples.

Note

Proper determination of sample dilution factor is important to make sure IP measurement is in calibration linear range.

- 5 Prepare IP standard stock at 100 millimolar (mM) in MeOH. Prepare the highest concentration IP calibration standard at 5 millimolar (mM) (430 mg/L) by diluting IP standard stock in calibration standard diluent.

Note


IP standard calibration range can be adjusted to fit estimated sample IP amount and chosen sample dilution factor.

- 6 Perform serial dilution of the highest concentration standard in calibration standard diluent by a dilution factor of 2 .

- 7 Leave one well containing Calibration standard diluent as Blank. At least  50 μL of blank and each standard are needed.

Sample preparation

- 8 Calculate the volumes of general diluent solution and test samples that are needed to prepare  100 μL Per sample with predetermined sample dilution factors. The following steps use dilution factor of 2 as an example.
- 9 Transfer 50 μL of general diluent solution (5% Methanol with the internal standard) to wells of Bio-Rad 96-well PCR plate that were designated for samples. Leave wells designated to Blank, QC, Calibration standards, and controls empty.



Example analytical plate layout. Calibration standards could be arranged together with samples in one plate, or be prepared in a separate plate. For every plate that contains test samples, controls and QC are required, but their well locations could be random.

Note

Using automated liquid handler reduces human error and increases accuracy and throughput, especially when multiple dilution factors are used and sample locations are randomized.

- 10 Transfer the test samples (at least 300 uL per well) to the sample well locations of 96 deep-well plates, then centrifuge at 3,200 rcf at 4 °C for 30 minutes.

Note

The centrifugation time could be reduced when spinning at a higher centrifugation force. Alternatively, samples matrix could be removed through filtration process using commercially available filter plates, such as Agilent Captiva EMR-Lipid plate (Part Number: 5190-1001).

- 11 Transfer 50 uL of test samples from 96 deep-well sample plate, into the Bio-Rad 96-well PCR plate containing the general diluent solution (see Step 9).

Note

Using automated liquid handler reduces human error and increases accuracy and throughput, especially when multiple dilution factors are used and sample locations are randomized.

- 12 Finish preparing analysis plates by manually filling Blank, QC, Calibration standards and controls, then sealing the sample plate with LC compatible non adhesive alumni seal. Store the plates in -20 °C if they are not processed immediately.

Prepare RapidFire-MS system

- 13 Condition RapidFire pumps by purging pumps 1, 2, and 3 at 5 mL/min flow rate, and pump 4 at 10 RPM for 5 minutes.
- 14 Change the flow rates of pump 1, 2, and 3 to that of the running method, and change pump 4 speed to 3 RPM. Close the purge valves of pump 1, 2, and 3 afterwards.

	A	B
	RapidFire pump	Flow rate (mL/min)

	A	B
	Pump 1	1.25
	Pump 2	0.5
	Pump 3	0.5

RapidFire pump flow rates for IP and internal standard quantitation

- 15 Check the back pressure of RapidFire high-throughput mass spectrometry system. Monitor all pump's pressure until they are stabilized at expected range.

	A	B	C	D
	RapidFire state	Pump 1 pressure (Mpa)	Pump 2 pressure (Mpa)	Pump 3 pressure (Mpa)
	Aspirate	4.06 ~ 4.15	0.67 ~ 0.72	1.55 ~ 1.70
	Load	4.08 ~ 4.15	0.67 ~ 0.72	1.55 ~ 1.70
	Extra Wash	0.72 ~ 0.74	1.46 ~ 1.47	1.55 ~ 1.70
	Elute	0.72 ~ 0.74	0.72 ~ 0.79	2.20 ~ 2.23

RapidFire-QqQ system pumps back pressures at the method flow rates in each cycle state

Note

Pay close attention to the pressure of all the pumps. Make sure they are stable at the expected range before starting your run. Diagnose and service RapidFire if pump pressure exceeds normal range.

- 16 Load all sample plates in RapidFire stacker.



RapidFire plate stacker with loaded sample plates

RapidFire and QQQ method parameters

- 17 Create and save a RapidFire method using the following RapidFire Cycle parameters and pump flow rates:

	A	B	C
State		Time (ms)	Flow rate (mL/min)
Aspirate		600	
Load		3000	1.25
Extra Wash		0	0.5

	A	B	C
	Elute	5000	0.5
	Reequilibrate	500	1.25

RapidFire Cycle duration and flow rates for IP and internal standard quantitation

- 18 Install and condition RapidFire C18  4 µL Type C cartridge (Agilent Part Number: G9205A).

Note

Tracking sample counts on the working C18 cartridge and Install multiple ones if the number of total injection samples is greater than the maximum allowed injections per cartridge. Use the same cartridge type name if multiple C18 cartridges are installed.

- 19 Create and save a Agilent 6460 QqQ acquisition method that uses the following source parameters:

	A	B
	Gas Flow (L/min)	10
	Gas Temperature	300
	Vaporizer Temperature	300
	Nebulizer pressure	30
	Capillary Voltage	3500
	Corona Current (uA)	4

Source parameters setting for IP and internal standard quantitation

- 20 Update the Agilent 6460 QqQ acquisition method with the following MRM acquisition parameters:

	A	B	C	D	E	F
Compound name	Precursor (m/z)	Product (m/z)	Dwell time (ms)	CE	Polarity	
Cis-3-Hexenyl Acetate	143.1	83.2	50	5	Positive	
Isoprenol	69.2	41.2	50	12	Positive	

MRM acquisition parameters setting for IP and internal standard quantitation

Set up RapidFire-MS analysis


21 Create plate maps that contain the desired RapidFire injection sequences.

Additional Functions: Delay 30 s, Flush 60 s

Current Sequence:
 WASH1 WASH1 WASH2 WASH2 A1 A2 A3 A4 A5 A6 A7 A8 A9 A10 A11 A12 B1
 B2 B3 B4 B5 B6 B7 B8 B9 B10 B11 B12 C1 C2 C3 C4 C5 C6 C7 C8 C9 C10 C11 C12
 D1 D2 D3 D4 D5 D6 D7 D8 D9 D10 D11 D12 E1 E2 E3 E4 E5 E6 E7 E8 E9 E10
 E11 E12 F1 F2 F3 F4 F5 F6 F7 F8 F9 F10 F11 F12 G1 G2 G3 G4 G5 G6 G7 G8 G9
 G10 G11 G12 H1 H2 H3 H4 H5 H6 H7 H8 H9 H10 H11 H12 WASH1 WASH1
 WASH2 WASH2


Example platemap that contains RapidFire injection sequence of a full 96-well sample plate

- 22 Create and save a batch file when the application runs single plate analysis or plates with multiple different injection sequences.




Example Batch file that contains plate info(Barcode etc) and their respective injection sequence, RapidFire method, Cartridge choice, and MS method.

- 23 Load the saved batch file under RapidFire UI File drop down menu. Click the green arrow button and follow the prompt to start the plateset.



RapidFire UI plateset run control panel

- 24 Create and save a Run Template file when the application runs high-throughput screen assay that analyze numerous full sample plates or plates with same injection sequence.




Example Run Template file that contains RapidFire injection sequence, RapidFire method, Cartridge choice, and MS method.

- 25 Create sample path file according to analytical sample plate layout. Fill in applied barcode information, well position, Sample name (File name), Sample type, assigned Calibration standard Levels and sample dilution factor(s).

Barcode	Well Position	sample name	Sampe name_id	Sample type	Level	DI	Run Template
RF_000043	A1	Culture_Tube_24C_0hr_001	Culture_Tube_24C_0hr_001	Sample			2_Yan_FullPlate_IP
RF_000043	A2	Culture_Tube_24C_0hr_002	Culture_Tube_24C_0hr_002	Sample			2_Yan_FullPlate_IP
RF_000043	A3	Culture_Tube_24C_0hr_003	Culture_Tube_24C_0hr_003	Sample			2_Yan_FullPlate_IP
RF_000043	A4	Culture_Tube_24C_0hr_004	Culture_Tube_24C_0hr_004	Sample			2_Yan_FullPlate_IP
RF_000043	A5	Culture_Tube_24C_0hr_005	Culture_Tube_24C_0hr_005	Sample			2_Yan_FullPlate_IP
RF_000043	A6	Culture_Tube_24C_0hr_006	Culture_Tube_24C_0hr_006	Sample			2_Yan_FullPlate_IP
RF_000043	A7	Culture_Tube_24C_0hr_007	Culture_Tube_24C_0hr_007	Sample			2_Yan_FullPlate_IP
RF_000043	A8	Culture_Tube_24C_0hr_008	Culture_Tube_24C_0hr_008	Sample			2_Yan_FullPlate_IP
RF_000043	A9	Culture_Tube_24C_0hr_009	Culture_Tube_24C_0hr_009	Sample			2_Yan_FullPlate_IP
RF_000043	A10	Culture_Tube_24C_0hr_010	Culture_Tube_24C_0hr_010	Sample			2_Yan_FullPlate_IP
RF_000043	A11	Culture_Tube_24C_0hr_011	Culture_Tube_24C_0hr_011	Sample			2_Yan_FullPlate_IP
RF_000043	A12	Culture_Tube_24C_0hr_012	Culture_Tube_24C_0hr_012	Sample			2_Yan_FullPlate_IP
RF_000043	B1	Culture_Tube_24C_0hr_013	Culture_Tube_24C_0hr_013	Sample			2_Yan_FullPlate_IP
RF_000043	B2	Culture_Tube_24C_0hr_014	Culture_Tube_24C_0hr_014	Sample			2_Yan_FullPlate_IP
RF_000043	B3	Culture_Tube_24C_0hr_015	Culture_Tube_24C_0hr_015	Sample			2_Yan_FullPlate_IP
RF_000043	B4	Culture_Tube_24C_0hr_016	Culture_Tube_24C_0hr_016	Sample			2_Yan_FullPlate_IP
RF_000043	B5	Culture_Tube_24C_0hr_017	Culture_Tube_24C_0hr_017	Sample			2_Yan_FullPlate_IP
RF_000043	B6	Culture_Tube_24C_0hr_018	Culture_Tube_24C_0hr_018	Sample			2_Yan_FullPlate_IP
RF_000043	B7	Culture_Tube_24C_0hr_019	Culture_Tube_24C_0hr_019	Sample			2_Yan_FullPlate_IP
RF_000043	B8	Culture_Tube_24C_0hr_020	Culture_Tube_24C_0hr_020	Sample			2_Yan_FullPlate_IP
RF_000043	B9	Culture_Tube_24C_0hr_021	Culture_Tube_24C_0hr_021	Sample			2_Yan_FullPlate_IP
RF_000043	B10	Culture_Tube_24C_0hr_022	Culture_Tube_24C_0hr_022	Sample			2_Yan_FullPlate_IP
RF_000043	B11	Culture_Tube_24C_0hr_023	Culture_Tube_24C_0hr_023	Sample			2_Yan_FullPlate_IP
RF_000043	B12	Culture_Tube_24C_0hr_024	Culture_Tube_24C_0hr_024	Sample			2_Yan_FullPlate_IP
RF_000043	C1	Culture_Tube_24C_0hr_025	Culture_Tube_24C_0hr_025	Sample			2_Yan_FullPlate_IP
RF_000043	C2	Culture_Tube_24C_0hr_026	Culture_Tube_24C_0hr_026	Sample			2_Yan_FullPlate_IP
RF_000043	C3	Culture_Tube_24C_0hr_027	Culture_Tube_24C_0hr_027	Sample			2_Yan_FullPlate_IP
RF_000043	C4	Culture_Tube_24C_0hr_028	Culture_Tube_24C_0hr_028	Sample			2_Yan_FullPlate_IP
RF_000043	C5	Culture_Tube_24C_0hr_029	Culture_Tube_24C_0hr_029	Sample			2_Yan_FullPlate_IP
RF_000043	C6	Culture_Tube_24C_0hr_030	Culture_Tube_24C_0hr_030	Sample			2_Yan_FullPlate_IP
RF_000043	C7	Culture_Tube_24C_0hr_031	Culture_Tube_24C_0hr_031	Sample			2_Yan_FullPlate_IP
RF_000043	C8	Culture_Tube_24C_0hr_032	Culture_Tube_24C_0hr_032	Sample			2_Yan_FullPlate_IP


Example sample path file. Barcode column information needs to match applied barcode to corresponding plate. Well position, sample name, sample name_id and sample type are parsed from sample plate layout information. Calibration standard Levels are required for setting up quantitation method in Agilent QQQ quantitative analysis software. Dilution factor is a variable that could be modified accordingly. Run template file is generated and updated according to step 24.

26 Generate a Map File under RapidFire UI File drop down menu.



Example map file that parses information from sample path file to data file.

- Import sample data file based on the map file generated above. Click import sample data under RapidFire UI file drop down menu. Select sample path file generated in step 25 and Map file generated in step 26 to fill in Sample Data File section.



Example result after import sample data.

Note

Check the end of the import log to make sure that all samples are imported without errors, warnings and skipped samples.

28 Click the green arrow button and follow the prompt to start the plateset.

Data analysis

- 29 Convert RapidFire sequence data file(s) to individual injection data files.

Note

Sample data files could miss if 'Do not covert missed sips' option in the converting process is checked.

- 30 Create a new batch in Agilent MassHunter QQQ quantitative analysis software, and load all individual injection data files. If this is the first analysis of IP, create a new quantitative method from acquired IP standard data file.


Note

Other software such as the Agilent *RapidFire* Integrator could be also used for analyzing RapidFire MS data. Refer to each software user guide for details.

- 31 Validate and save the quantitative method after completing all tasks in the method setup sections. Apply the quantitative method and analyze the sample batch.

Note

The calibration standard levels and expected concentrations must match the levels defined in the sample batch table. Refer to MassHunter QQQ quantitative analysis software user manual for more details on setting up the quantitative method.



Example IP calibration curve that shows linear fitting ranging from 7.8 mg/L to 1000 mg/L.

- 32 Save sample batch and Export quantitative reports from batch table for further data processing and visualization.

Acknowledgements

This work was supported by the DOD Tri-Service Biotechnology for a Resilient Supply Chain (T-BRSC) program through Award AWD00007196 and by the Joint BioEnergy Institute (<http://www.jbei.org>), U.S. Department of Energy, Office of Science, Biological and Environmental Research Program under Award DE-AC02-05CH1123 to Lawrence Berkeley National Laboratory.

Competing Interests:

BAS has a financial interest in Illium Technologies, Caribou Biofuels, and Erg Bio. None of the other authors have an outside financial interest to disclose.