

Jan 23, 2025

Hippocampal slice preparation for electrophysiology

DOI

<https://dx.doi.org/10.17504/protocols.io.eq2ly6j9egx9/v1>

Quinn Pauli¹, Robert Bonin¹

¹Leslie Dan Faculty of Pharmacy, University of Toronto, Toronto, Ontario, Canada



Quinn Pauli

University of Toronto

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.eq2ly6j9egx9/v1>

Protocol Citation: Quinn Pauli, Robert Bonin 2025. Hippocampal slice preparation for electrophysiology . **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.eq2ly6j9egx9/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 05, 2024

Last Modified: January 23, 2025

Protocol Integer ID: 114347

Keywords: Electrophysiology, hippocampus, mouse, slice, brain dissection, hippocampus slices, maintenance of acute hippocampal slice, acute hippocampal slice, hippocampal slice, mice for electrophysiology, acute slice, slice preparation, electrophysiology, healthy slice, synaptic plasticity, day electrophysiology experiment, slicing, dissection, preparation

Funders Acknowledgements:

Natural Sciences and Engineering Research Council of Canada

Grant ID: RGPIN-2024-06215

Abstract

This protocol describes the preparation and maintenance of acute hippocampal slices from mice for electrophysiology. We describe the dissection, slicing and incubation methods to produce healthy slices. Slice preparation can be followed by same-day electrophysiology experiments (e.g. to assess synaptic plasticity). Acute slices can be further processed for protein detection (e.g. immunofluorescence, immunoblotting) as shown by our lab.

Materials

Solution Recipes:

	A	B	C	D	E
	Reagent	Mass (g) for 1 L	Mass (g) for 500 mL	MW	mM
	NaCl	7.2466	3.6233	58.44	124
	Glucose	1.8016	0.9008	180.156	10
	NaHCO ₃	2.1843	1.0922	84.007	26
	KCl	0.2236	0.1118	74.55	3
	NaH ₂ PO ₄ •H ₂ O	0.1932	0.0966	137.99	1.4
	MgSO ₄	1 mL	0.5 mL	1M Solution	1
	CaCl ₂	2 mL	1 mL	1M Solution	2
	Osmolarity 300-310 mOsm				

ACSF Recipe

	A	B	C	D	E
	Reagent	Mass (g) for 1 L	Mass (g) for 500 mL	MW	mM
	Sucrose	17.1145	8.5573	342.3	50
	NaCl	5.3765	2.6883	58.44	92
	Glucose	2.7024	1.3512	180.156	15
	KCl	0.3728	0.1864	74.55	5
	NaH ₂ PO ₄ •H ₂ O	0.1987	0.0994	137.99	1.4
	NaHCO ₃	2.1843	1.0922	84.007	26

	A	B	C	D	E
	CaCl ₂	0.5 mL	0.25 mL	1M Solution	0.5
	MgCl ₂ •6H ₂ O	1.4231	0.7116	203.30	7
	Kynurenic Acid	0.1892	0.0946	189.17	1
	Osmolarity 322-330 mOsm				

Sucrose Solution Recipe

Dissection tools:

- Fine scissors
- Angled fine scissors
- Surgical scissors
- Dumont Forceps
- Hemostat
- Dressing forceps
- Scalpel handle and blade (#21)
- Spatula
- Razor blades
- Glass petri dish

Other supplies:

- 25 or 50 mL syringe
- Tygon tubing
- 27g needle (for i.p. injection)
- 16g needle (for perfusion)
- Anesthetic (e.g. Avertin)
- Filter paper
- Tissue adhesive (e.g. 3M Vetbond)
- Ice
- 95% oxygen/5% carbon dioxide (carbogen tank and lines)
- Cell strainers, falcon tubes, super glue and a 250 mL beaker (for homemade brain slice keeper)
- Beakers (250-500 mL)

Equipment:

- Vibratome
- Brain slice keeper



Brain slice keeper

Troubleshooting

Safety warnings

- ⚠ If you choose to use this protocol, you do so at your own risk and must ensure that any local guidance is adhered to.

Ethics statement

Experiments involving animals must be conducted according to internationally-accepted standards and should always have prior approval from an Institutional Animal Care and Use Committee (IACUC) or equivalent ethics committee(s). Animal care and experimental procedures were reviewed and approved by the Animal Ethics and Compliance Program at the University of Toronto and conducted in accordance with the Canadian Council on Animal Care (CACC) guidelines.



Set up

55m

- 1 Prepare solutions (see Materials) 1 week to 1 day in advance (store at $4\text{ }^{\circ}\text{C}$) II
- 2 Bubble 100 mL sucrose solution and 200 mL ACSF for 20 minutes in 95% oxygen/5% carbon dioxide for $00:20:00$ 20m
- 3 Freeze sucrose solution at $-20\text{ }^{\circ}\text{C}$ until a slurry begins to form (~ $00:15:00$) 15m
- 4 Prepare fresh ACSF (day of) that will be used for recording 20m
- 5 Once prepared, continue bubbling all solutions with 95% oxygen/5% carbon dioxide throughout the duration of slice preparation
- 6 Add ACSF to a brain slice keeper (e.g. 250 mL beaker with cell strainers affixed to cut falcon tubes to ensure continuous solution flow; see Materials) and heat in a water bath to $28\text{ }^{\circ}\text{C}$. Ensure that the cell holders are submerged and the solution is bubbled with 95% oxygen/5% carbon dioxide for the duration of the procedure.
- 7 Prepare dissection area with dissection tools, an absorbent pad, and bench protectors. Add ice to vibratome well (around buffer tray)

Perfusion and brain dissection

12m

- 8 Weigh and anesthetize mouse with the appropriate dosage (i.p.) and place back into home cage. 5m
- 9 Assess depth of anesthesia using toe-pinch response (ensure response is gone before beginning procedure)
- 10 Add ice cold sucrose solution to 50mL syringe attached to tubing and needle for perfusion
- 11 Make a lateral incision through the abdomen at the xiphoid process using fine scissors



- 12 Pulling up on the xiphoid process, carefully cut the diaphragm and then cut along the base of the ribcage on both sides to expose the heart
- 13 Clamp the xiphoid process using a hemostat to expose the chest cavity
- 14 Securing the heart with dressing forceps, insert 18g needle attached to a 50mL syringe with tubing into left ventricle
- 15 Make a small incision in the right ventricle using fine scissors
- 16 Compress the syringe slowly to perfuse the mouse until the liver is visibly clear or approximately 30-40 mL of solution has been used 2m
- 17 Decapitate the mouse with large surgical scissors
- 18 Pull the skin back and make two lateral incisions at the base of the skull using fine scissors
- 19 Make a cut down the midline of the skull using angled fine scissors, taking care not to damage the brain
- 20 Gently peel away the two sections of the skull with Dumont forceps
- 21 Scoop the brain out of the skull using a spatula and immediately submerge in cold, bubbling ACSF for  00:05:00 5m



Slicing

1h 45m

- 22 Wet filter paper in ice cold ACSF and place in a glass petri dish
- 23 Remove the brain from solution and place on the wet filter paper. Use a scalpel to remove the cerebellum and olfactory bulb, then make two thin cuts along the sagittal plane of the brain bilaterally to provide a flat surface for slicing

- 24 Make a cut down the midline of the brain and glue the two hemispheres of the brain to the vibratome platform, with the hippocampus facing up and cortex towards the blade
- 25 Fill the buffer chamber with ice cold ACSF and continuously bubble the solution. Attach the razor blade to the vibratome and position it right above the tissue as the starting position
- 26 Slice $\pm 400 \mu\text{m}$ slices at 0.08 mm/s (on continuous slicing mode). Carefully remove the CA3 region of each slice using a scalpel. Remove any excess tissue around the dorsal hippocampal slices. Pause the vibratome when performing microdissections of slices.
- 27 Transfer each slice to $28 \text{ }^\circ\text{C}$ bubbling ACSF in the brain slice keeper to rest for a minimum of 01:30:00 before recording in a chamber maintained at $28 \text{ }^\circ\text{C}$
- 28 If recording from multiple slices, turn off the water bath after 03:00:00 of incubation and allow slices to gradually come to room temperature.

15m

1h 30m

