


Oct 24, 2020

## His-tag purification

DOI

[dx.doi.org/10.17504/protocols.io.bnw5mfg6](https://dx.doi.org/10.17504/protocols.io.bnw5mfg6)



Andreea S<sup>1</sup>

<sup>1</sup>University of Groningen

iGEM Groningen 2020



a.stan.6

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bnw5mfg6>

**Protocol Citation:** Andreea S 2020. His-tag purification . protocols.io <https://dx.doi.org/10.17504/protocols.io.bnw5mfg6>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



**Protocol status:** Other

The protocol is developed based on literature and has not been tested yet.

**Created:** October 24, 2020

**Last Modified:** October 24, 2020

**Protocol Integer ID:** 43709

**Keywords:** tag purification, purification, immobilised metal affinity chromatography, metal affinity chromatography, affinity of the protein, histidine residue, tag with high purity, protein, other protein, histidine, metal ion such as nickel, amino acids constituting protein, metal ion, chelating agent, end of the protein, number of histidine, genetic engineering, agent such as iminodiacetic acid, transition metal ion, tag, coordinate bond with metal ion, ion, iminodiacetic acid


## Abstract

His tag purification uses the technique of immobilised metal affinity chromatography. In this technique, transition metal ions are immobilized on a resin matrix using a chelating agent such as iminodiacetic acid. It has been studied that among amino acids constituting proteins, histidine is strongly involved in the coordinate bond with metal ions. Therefore, if a number of histidines are added to the end of the protein by genetic engineering, the affinity of the protein for the metal ion is remarkably increased and the basic idea is that purification can be easily carried out. When a protein having a His tag is brought into contact with a carrier on which a metal ion such as nickel is immobilized, the histidine residue chelates the metal ion and binds to the carrier. Since other proteins do not bind to the carrier, they can be washed off with a buffer. Thereafter, it is possible to recover the protein having the His tag with high purity.

## Troubleshooting







## Separate proteins from soil matrix

- 1 Collect soil samples of  5 g
- 2 Extract total proteins using NoviPure Soil Protein Kit or other commercially available kits for total protein soil extraction












## His-tag separation of NLP

1h

- 3 Wash the  $\text{Ni}^{2+}$ -sepharose column material with 12 CVs of MQ and 4 CVs of Column Wash Buffer (  10 millimolar (mM) imidazole, KPi  50 millimolar (mM)  7 , NaCl  200 millimolar (mM) ).




### Note

Use  $\pm 0.5$  ml of  $\text{Ni}^{2+}$ -sepharose column material per 10 mg of total protein.

- 4 Apply the sample, add imidazole (10mM final concentration) and the washed  $\text{Ni}^{2+}$ -sepharose column material. Nutate in at  4 °C for  01:00:00
- 5 Pour column, collect flow through to apply on SDS gel.
- 6 Wash column with 20 CVs of Wash Buffer (  50 millimolar (mM) imidazole, KPi  50 millimolar (mM)  7 , NaCl  200 millimolar (mM) ).
- 7 Elute protein with Elution Buffer (  500 millimolar (mM) imidazole, Kpi  50 millimolar (mM)  7 , NaCl  200 millimolar (mM) ) in  200  $\mu\text{L}$  fractions. Check elution fractions Absorbance by NanoDrop.
- 8 Run an SDS gel to check purification:

1h



- Soil suspension & Flow through: dilute 15x, apply  5  $\mu\text{L}$
- Wash: dilute 1.25x, apply  10  $\mu\text{L}$
- Elution fractions: dilute to  $\pm 0.2$  mg/ml, apply  5  $\mu\text{L}$