



Jul 18, 2025

## HIV-1 partial-*pol* nanopore sequencing protocol

DOI

<https://dx.doi.org/10.17504/protocols.io.261gekojdg47/v1>

Ganna Kovalenko<sup>1,2</sup>, Anna Yakovleva<sup>3</sup>, Tetyana Vasylyeva<sup>4</sup>

<sup>1</sup>University of California, Irvine; <sup>2</sup>ARTICNetwork; <sup>3</sup>University of Oxford, UK; <sup>4</sup>University of California Irvine



Ganna Kovalenko

University of California, Irvine, ARTICNetwork

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.261gekojdg47/v1>

**Protocol Citation:** Ganna Kovalenko, Anna Yakovleva, Tetyana Vasylyeva 2025. HIV-1 partial-*pol* nanopore sequencing protocol. protocols.io <https://dx.doi.org/10.17504/protocols.io.261gekojdg47/v1>

**Manuscript citation:**

Kovalenko, G., Liulchuk, M. G., Filippovych, M., Smyrnov, P., Strathdee, S. A., & Vasylyeva, T. I. (2026). VANTAGE: van-based real-time HIV sequencing for transmission mapping and drug resistance profiling in war-affected Ukraine. *AIDS (London, England)*, 40(1), 123–126. <https://doi.org/10.1097/QAD.0000000000004360>

Kovalenko, G., Yakovleva, A., Smyrnov, P., Redlinger, M., Tymets, O., Korobchuk, A., Kolodiazieva, A., Podolina, A., Cherniavska, S., Skaathun, B., Smith, L. R., Strathdee, S. A., Wertheim, J. O., Friedman, S. R., Bortz, E., Goodfellow, I., Meredith, L., & Vasylyeva, T. I. (2023). Phylodynamics and migration data help describe HIV transmission dynamics in internally displaced people who inject drugs in Ukraine. *PNAS nexus*, 2(3), pgad008. <https://doi.org/10.1093/pnasnexus/pgad008>

**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** July 16, 2025

**Last Modified:** July 18, 2025

**Protocol Integer ID:** 222678

**Keywords:** HIV-1, pol gene, drug resistance, nanopore sequencing, genomic surveillance, molecular epidemiology, native barcoding, MinION sequencing, HIV transmission , targeted sequencing, amplicon sequencing, nanopore library preparation, pol nanopore, reference matching to the hxb2 genome, hxb2 reference genome, protocol reverse transcription, sequencing protocol, hxb2 genome, oxford nanopore technology, transcriptase, nested pcr, reverse transcription, bioinformatics workflow suggestion, transcription, pcr



## Abstract

Reverse transcription and amplification of the HIV-1 partial-*pol* gene were performed using primers and PCR conditions adapted from [Zhou et al. \(2011\)](#), originally optimized for drug resistance surveillance in resource-limited settings. The resulting amplicons (~1.3 kb, covering the protease (PR) and reverse transcriptase (RT) regions) were used for nanopore library preparation.

*Optionally, a nested PCR can be performed using the second-round primers described in the same study, yielding a product of 1,084 nt. However, we found that a single round of RT-PCR was sufficient for successful amplification, reducing both time and cost.*

Library preparation followed a modified version of the [ARTIC LoCost protocol](#), adapted for partial-*pol* amplicons and compatible with the Native Barcoding Kit 96 V14 (SQK-NBD114.96, Oxford Nanopore Technologies).

### Bioinformatics Workflow Suggestion:

[RAMPART](#) can be run concurrently with ONT MinKNOW sequencing to provide real-time visualisation of read mapping, coverage, and reference matching to the HXB2 genome for each barcode. The RAMPART configuration protocol is available at [https://github.com/AnyaKovalenko/rampart-protocol-HIV-1\\_partial-pol](https://github.com/AnyaKovalenko/rampart-protocol-HIV-1_partial-pol). For reference-based consensus generation, align reads to the HXB2 reference genome (GenBank accession: [K03455.1](#)), or specifically to the *pol* region (coordinates: 2030–3400), using Minimap2, followed by polishing with Racon and Medaka.

*Please cite the associated publications when using or adapting this workflow.*

## Materials

| A                                                                                     | B                   | C             |
|---------------------------------------------------------------------------------------|---------------------|---------------|
| Reagent                                                                               | Supplier            | Catalog #     |
| QIAamp Viral RNA Mini Kit (50)                                                        | Qiagen              | 52904         |
| Nuclease free water                                                                   | Omega Bio-Tek       | PD092         |
| SuperScript III One-Step RT-PCR System with Platinum Taq High Fidelity DNA Polymerase | Invitrogen          | 12574035      |
| PCR Clean - DX                                                                        | Aline BioSciences   | C-1003-5      |
| QuantiFluor ONE dsDNA System                                                          | Promega             | E4870         |
| Thermo Scientific GeneRuler 1kb DNA Ladder, ready-to-use                              | Thermo Scientific   | 10809360      |
| All-in-one tablets with pre-measured agarose, TBE, and SeeGreen stain                 | GELATO              | RG-1500-20    |
| 10X TBE Buffer (Tris Borate EDTA) - 500 mL                                            | ufcbio              | BTBE-10X      |
| HIV pol primers                                                                       | IDT                 |               |
| NEBNext Ultra II End Repair/dA-tailing module                                         | New England BioLabs | E7546S        |
| Blunt/TA Ligase Master Mix                                                            | New England BioLabs | M0367S        |
| NEBNext Quick Ligation Module                                                         | New England BioLabs | E6056S        |
| Native Barcoding Kit 96 V14                                                           | ONT                 | SQK-NBD114.96 |
| Flow Cell Wash Kit                                                                    | ONT                 | EXP-WSH004    |
| Flow cells R10.4.1                                                                    | ONT                 | FLO-MIN114    |
| Ethanol (96–100%)                                                                     | Generic             |               |

## I. RNA Extraction

- 1 QIAamp Viral RNA Minikit (Qiagen), following the manufacturer's instruction

## II. One-Step RT-PCR

- 2 NOTE: Prepare primer:

|  | A | B           | C                              | D       | E            | F                | G              |
|--|---|-------------|--------------------------------|---------|--------------|------------------|----------------|
|  |   | Primer Name | Sequence 5'-3'                 | Type    | Product Size | Position in HXB2 | Purpose        |
|  | 1 | PRTM-F1a*   | TGAARGAITGYACTGARAGR CAGGCTAAT | forward | 1,300 nt     | 2057–2085        | One-Step RTPCR |
|  |   | PRTM-F1b*   | ACTGARAGRCAGGCTAATT TTTTAG     |         |              | 2068–2092        |                |
|  | 2 | RT-R1       | ATCCCTGCATAAATCTGAC TTGC       | reverse |              | 3370–3348        |                |

\*PRTM-F1 is a mixture of primers F1a and F1b at a ratio of 1:1 (w/w) and is used as the forward primer.

### Note

*The table includes primers for the One-Step RT-PCR only. If you also wish to perform the nested PCR step, please refer to [Zhou et al. \(2011\)](#). Once amplification is complete, you may proceed with the library preparation steps outlined in this protocol. We found that a single round of RT-PCR was sufficient for successful amplification, reducing both time and cost.*

- 3 Prepare 100 µM primer stocks for all primers: add NFW or TE to the primer based on its amount (in nanomoles) to achieve the desired concentration.  
Mix equal volume of the 100 µM primer stocks of PRTM-F1a and PRTM-F1b in one tube:  
ei., 40ul + 40ul =  
80ul (now this is your mixture of the forward primer 100 µM). Prepare 10 µM working solution for both forward and reverse primers.
- 4 Prepare RTPCR reaction mix for the required number of reactions plus 1-2 additional reaction for overage. A negative control (without template RNA) should be included in



every experiment:

| A                                              | B                   |
|------------------------------------------------|---------------------|
| Component                                      | Volume (ul) per rxn |
| 2X Reaction Mix                                | 25                  |
| SSIII RT/Platinum Taq High Fidelity Enzyme Mix | 1                   |
| 10 $\mu$ M forward PRTM-F1ab                   | 1                   |
| 10 $\mu$ M reverse RT-R1                       | 1                   |
| template RNA                                   | 20*                 |
| NFW                                            | 2                   |
| <b>Total</b>                                   | <b>50</b>           |

\*for high viral load samples 10 ul RNA can be used.

5 Start the following program on the PCR thermocycler:

| A                   | B           | C      | D      |
|---------------------|-------------|--------|--------|
| Step                | Temperature | Time   | Cycles |
| RT                  | 50°C        | 45 min | 1      |
| Enzyme inactivation | 94°C        | 2 min  | 1      |
| Denature            | 94°C        | 15 sec | 40     |
| Anneal              | 50°C        | 20 sec |        |
| Extend              | 72°C        | 2 min  |        |
| Final extension     | 72°C        | 10 min | 1      |
| Hold                | 4°C         | -      | -      |

#### Note

*PCR product can be confirmed by 1% agarose gel electrophoresis. Based on our experience, the bands may appear faint or barely visible for samples with low viral load, but this is typically sufficient for successful sequencing.*



### III. Amplicon Clean Up

11m 30s

- 6 Add  50  $\mu\text{L}$  of magnetic beads (1 $\times$  ratio) to each well/tube. Mix gently by pipetting (or by flicking). Briefly spin down and incubate for  00:05:00 at  Room temperature . 5m
- 7 Place on magnetic rack and incubate for  00:05:00 or until the beads have pelleted and the supernatant is completely clear. 5m
- 8 Carefully remove and discard the supernatant, being careful not to touch the bead pellet.
- 9 Keep on magnet, wash the pellet with  200  $\mu\text{L}$  of freshly prepared 80% ethanol without disturbing the pellet. After  00:00:30 remove the ethanol using a pipette and discard being careful not to touch the bead pellet. 30s
- 10 Repeat the ethanol wash step  [go to step #9](#)
- 11 Pulse centrifuge to collect all liquid at the bottom of the tube and carefully remove as much residual ethanol as possible using a P10 pipette.
- 12 With the tube lid open allow the pellet to dry for  00:01:00 or until the pellet loses its shine, but do not dry the pellet to the point of cracking. 1m
- 13 Remove the tube from the magnetic rack and resuspend pellet in  25  $\mu\text{L}$  NFW, mix gently by flicking and incubate for  00:05:00 at  Room temperature .
- 14 After the incubation place on magnet. Pellet the beads on the magnet until the eluate is clear and colourless (~3min) and elute  25  $\mu\text{L}$  to a clean 1.5 mL Eppendorf LoBind tube ensuring no beads are transferred into this tube.

### IV. Quantification

15

#### Note

use only thin-wall, clear, 0.5 ml PCR tubes (acceptable tubes include Qubit/Quantus assay tubes).  
The fluorometer should be calibrated before use. Please follow the calibration instructions specific to your instrument.

Take  1  $\mu\text{L}$  of sample for quantification with Quantus using QuantiFluor ONE dsDNA Dye. Record the final sample concentration results.

#### Note

PCR product concentrations for the ~1.3 kb HIV-1 pol amplicon may be low (e.g.,  $<1 \text{ ng}/\mu\text{L}$ ), especially when working with degraded RNA, such as that extracted from dried blood spots (DBS) or low viral load plasma samples. However, successful sequencing with high read depth (e.g.,  $>1000\times$ ) is often achievable.

## V. End Prep

1h

16 In a new PCR strip-tube/plate set up the following reaction for each sample:

| A                              | B                                       |
|--------------------------------|-----------------------------------------|
| Component                      | Volume (ul) per rxn                     |
| NFW                            | (optional, if high viral load samples)* |
| End Prep Reaction Buffer       | 1.2                                     |
| End Prep Enzyme Mix            | 0.5                                     |
| Amplicons (from previous step) | 8.3                                     |
| <b>TOTAL</b>                   | <b>10 ul</b>                            |

\*Nuclease-Free Water (NFW) is optional and used only to adjust the final reaction volume to 10  $\mu\text{L}$ . If the amplicon concentration is high, you may reduce the amplicon input volume and compensate with NFW (e.g., 3.3  $\mu\text{L}$  of amplicon + 5  $\mu\text{L}$  of NFW). However, if the amplicon concentration is low (e.g.,  $<5 \text{ ng}/\mu\text{L}$ ), it is recommended to use the maximum possible volume of amplicon (up to 8.3  $\mu\text{L}$ ) without NFW to ensure sufficient input for effective library preparation.





17 Mix gently by flicking and briefly spin down.

18 Set-up the following incubation using PCR thermocycler:

1h

🔥 20 °C for ⌚ 00:30:00

🔥 65 °C for ⌚ 00:30:00

🔥 4 °C Hold or place on ice

#### Note

*For best results EP products shouldn't be frozen before barcoding. For best results proceed directly to barcode ligation.*

## VI. Barcode Ligation

40m

19 Select a unique ONT Native barcode (Native Barcoding Kit 96 V14 (SQK-NBD114.96) ONT) for every sample to be run together on the same flow cell.

20 In a new PCR strip-tube set up the following reaction for each sample:

| A                          | B                                       |
|----------------------------|-----------------------------------------|
| Component                  | Volume (ul) per rxn                     |
| Blunt/TA Ligase Master Mix | 5                                       |
| Native Barcode (ONT)       | 1.25                                    |
| End-preped amplicons       | 3.75                                    |
| NFW                        | (optional, if high viral load samples)* |
| <b>TOTAL</b>               | <b>10</b>                               |

\*the same idea. See the previous step #16

NFW is optional and used only to adjust the final reaction volume to 10  $\mu$ L. If the amplicon concentration is high, you may reduce the end-preped amplicon input volume and compensate with NFW (e.g., 1.5  $\mu$ L of end-preped amplicon + 2.25  $\mu$ L of NFW).

21 Mix gently by flicking and briefly spin down.



22 Set-up the following incubation using PCR thermocycler:

40m

🌡️ 20 °C for ⌚ 00:30:00

🌡️ 70 °C for ⌚ 00:10:00

🌡️ 4 °C Hold or place on ice

#### Note

You can freeze at 🌡️ -20 °C (eg., overnight) the samples at this point if you need to pause the experiment.

## VII. Barcodes Pooling & Clean Up 1

24m

23 POOL barcoded amplicons in a clean 1.5 mL Eppendorf DNA LoBind tube:

- If processing <32 samples, pool all 🧪 10 µL from each barcoding reaction.
- If processing 33-47 samples, pool 🧪 8 µL from each barcoding reaction.
- If processing 48-80 samples, pool 🧪 5 µL from each barcoding reaction.
- If processing 81-96 samples, pool 🧪 3 µL from each barcoding reaction.

24 Vortex the SPRI magnetic beads to resuspend. Add 0.4× volume of beads to the pooled reaction and mix by flicking.

25 Incubate the reaction for ⌚ 00:05:00 at 🌡️ Room temperature .

5m

26 Place the tube into a magnetic rack for ⌚ 00:05:00 .

5m

27 Keep on magnet and remove all supernatants.



- 28 Wash the beads with  250  $\mu\text{L}$  SFB (ONT). Remove a tube from the magnetic rack and fully resuspend the pellet in SFB (gently by either flicking or pipetting). Then pulse centrifuge and return the tube to the magnet and allow the beads to pellet. Allow  00:05:00 contact time with the magnet. 5m
- 29 Repeat the SFB wash  go to step #27 and  go to step #28
- 30 Keep on magnet and remove all supernatants. Pulse centrifuge and remove any residual SFB. There is no need to air dry the beads after SFB washes.
- 31 Keep tube on the magnet and add  200  $\mu\text{L}$  of 80 % Ethanol to bathe the pellet. Carefully remove and discard ethanol, being careful not to touch the bead pellet.
- Note**
- Only perform 1x 80% ethanol wash*
- 32 Pulse centrifuge to collect all liquid at the bottom of the tube and return the tube on a magnet. Carefully remove as much residual ethanol as possible using a P10 pipette.
- 33 With the tube lid open incubate for  00:01:00 or until the pellet loses its shine (do not allow the pellet dries completely, it will crack and become difficult to resuspend). 1m
- 34 Remove the tube from the magnetic rack and add  32  $\mu\text{L}$  NFW water. Mix gently by flicking and incubate for  00:05:00 at  Room temperature . 5m
- 35 Place the tube into a magnetic rack for  00:03:00 or until the eluate is clear and colourless. 3m
- 36 Elute  32  $\mu\text{L}$  sample to a clean 1.5 mL Eppendorf DNA LoBind tube. Be careful not to take any beads! Use light!
- 37 Take  1  $\mu\text{L}$  for quantification using a fluorometer such as a Qubit or Quatus. Record the metric.

**Note**

You can store the pooled sample at 4 °C overnight and continue the library preparation the next day, if needed.

**VIII. Adapter Ligation**

20m

38 Set up the following Adapter Ligation reaction:

| A                         | B          |
|---------------------------|------------|
| Component                 | Volum (ul) |
| NEB Quick Ligation Buffer | 10         |
| Adapter mix ONT           | 5          |
| Quick T4 Ligase           | 5          |
| Pooled library            | 30         |
| <b>TOTAL</b>              | <b>50</b>  |

39 Mix gently by flicking and briefly spin down.

40 Incubate the reaction for 00:20:00 at Room temperature  
or if using PCR thermocycler: 20°C for 20 minutes and then 4°C Hold

20m

**IX. Final Clean Up**

22m

41

**Note**

1. **Do NOT use Ethanol!!!** for the final clean-up step as it would denature motor proteins on adaptors. SFB will remove excess adapter without damaging the adapter-protein complexes.
2. Vortex magnetic beads thoroughly before use to ensure they are well resuspended; the solution should be a homogenous brown colour.



Add  50  $\mu\text{L}$  magnetic beads (1 $\times$  ratio), and mix gently by flicking the tube. Pulse centrifuge to collect all liquid at the bottom of the tube.

42 Incubate for  00:05:00 at  Room temperature .

5m

43 Place the tube into a magnetic rack for  00:05:00

5m

44 Remove all supernatants and wash beads with  250  $\mu\text{L}$  SFB (ONT). Resuspend beads in SFB completely by gently pipette mixing/flicking. Pulse centrifuge to collect all liquid at the bottom of the tube.

45 Return on the magnet and allow the beads to pellet. Allow  00:05:00 contact time with the magnet.

5m

46 Repeat the SFB wash step  [go to step #44](#) and  [go to step #45](#)

47 After second SFB wash remove all supernatant. Then spin down and carefully remove any residual SFB. You do not need to allow to air dry with SFB washes.

48 Remove the tube from the magnet. Add  15  $\mu\text{L}$  of EB (ONT) and gently resuspend by flicking. Incubate  00:05:00 at  Room temperature .

5m

49 Place the tube into the magnetic rack for  00:02:00 or until the eluate is clear and colourless.

2m

50 Transfer  15  $\mu\text{L}$  of final library to a new 1.5mL Eppendorf DNA LoBind tube. Be careful not to take any beads!

51 Take  1  $\mu\text{L}$  for quantification using a fluorometer such as a Qubit or Quatus. Record the metric.

## X. MinION Sequencing

52 Prime and load the flow cell with the final library, following the current ONT protocol for the Native Barcoding Kit you are using (Native Barcoding Kit 96 V14).



#### Note

*For detailed instructions on priming and loading a MinION flow cell, see this [video](#)*

- 53 Start the sequencing run using MinKNOW.

#### Note

*If using live basecalling, make sure to enable "Barcode both ends" in the basecalling settings.*

*We also recommend using High Accuracy basecalling mode and enabling the "Trim barcodes" option in the barcoding settings.*

**HAPPY SEQUENCING!!!**