

Jan 09, 2026

HMW DNA extraction from red algae for Long Read Sequencing

DOI

<https://dx.doi.org/10.17504/protocols.io.e6nvwn577vmk/v1>



Benoît acherie¹, Karine Labadie¹

¹Genoscope/CEA

Benoît

HPM/Genoscope

1 more workspace



Benoît Vacherie

Genoscope/CEA

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.e6nvwn577vmk/v1>

Protocol Citation: Benoît acherie, Karine Labadie 2026. HMW DNA extraction from red algae for Long Read Sequencing . protocols.io <https://dx.doi.org/10.17504/protocols.io.e6nvwn577vmk/v1>



License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: January 07, 2026

Last Modified: January 09, 2026

Protocol Integer ID: 237186

Keywords: rhodophyta, algae, DNA extraction, HMW DNA extraction, sarkosyl, PVP, hmw dna extraction from red algae, hmw dna extraction, high molecular weight dna extraction, dna extraction, extraction of dna, red algae, read sequencing, long read sequencing, extraction, dna, high molecular weight, purification, mercaptoethanol,

Abstract

This protocol describes High Molecular Weight DNA extraction for Long Read Sequencing.

This protocol was adapted from Ramakrishnan et al. (2017; doi: 10.1007/s13205-017-0992-2), with some modifications.

This protocol is designed for the extraction of DNA from 200 mg of material that has been flash-frozen in liquid nitrogen and stored at -80°C .

The procedure involves cell lysis in a Sarkosyl-based buffer supplemented with PVP and β -mercaptoethanol, followed by purification using potassium acetate and chloroform:isoamyl alcohol.

Guidelines

Use only wide bore tips.

Work in a chemical hood when using β -mercaptoethanol and chloroform.

Allow the DNA to resuspend for at least of 24 hours before proceeding with QC.



Protocol materials

- ✕ Liquid Nitrogen
- ✕ RNase A (DNase free) 100mg/ml **Qiagen Catalog #19101**
- ✕ 70% ethanol
- ✕ DNA LoBind Tubes 1.5 ml **Eppendorf Catalog #022431021**
- ✕ Beta-mercaptoethanol **Catalog #97622-10X1ML**
- ✕ Isopropanol
- ✕ Sodium acetate, 3M **Fisher Scientific Catalog #J61928**
- ✕ Chloroform:Isoamyl Alcohol [24:1] **Merck MilliporeSigma (Sigma-Aldrich) Catalog #25666**
- ✕ Ethanol absolute
- ✕ Potassium Acetate
- ✕ Buffer TE 1x



Preparation of reagents

30m

- 1 **Extraction Buffer** : Prepare 10 mL of extraction buffer (per 200g of material) in a 50 mL tube

30m

	Reagents	Final concentration
	Tris-HCl ph8	100 mM
	EDTA	20 mM
	NaCl	1.5 M
	Sarkosyl	1%
	PVP40	2%
	H2O	QSP 10 ml

DNA Extraction

1h 26m

- 2 Add **10 ml of extraction buffer** to a 50 ml tube.
 10 mL Extraction buffer
- 3 Put a mortar in ice.
Cool the mortar and pestle with liquid nitrogen until the bubbling stops.
 Liquid Nitrogen
- 4 **Grind 200 mg of frozen sample** to a fine powder (approx. 2 min).
 200 mg Sample
- 5 **Transfer the powder** to the extraction buffer
- 6 Add **20 µl of Beta-mercaptoethanol** to the extraction buffer and vortex 5 sec.
 Beta-mercaptoethanol **Catalog #97622-10X1ML**
- 7 **Incubate at room temperature for 1 h** with agitation (300 rpm).

1m

5m

2m

1m

1m

1h



8 **Centrifuge** with low deceleration 15m
8500 x g, 4°C, 00:10:00 , Acc 9 / Dec 6

9 **Transfer the supernatant** to a new 50 ml tube. 1m

DNA Purification

2h 4m

10 Add to the supernatant : 1m
0.1 volume (1 ml) of **absolute ethanol** Ethanol absolute
0.25 volume (2.5 ml) of **3M potassium acetate** Potassium Acetate
1 volume (10 ml) of **chloroform isoamyl alcohol**

Chloroform:Isoamyl Alcohol [24:1] Merck MilliporeSigma (Sigma-Aldrich) Catalog #25666

Vortex 2× 5 sec

11 **Incubate for 20 min at -20°C** 20m

12 **Centrifuge** with low deceleration 25m
8500 x g, 4°C, 00:20:00 , Acc 9 / Dec 6

13 Gently **transfer the aqueous phase** to a new 50 ml tube. 1m

14 Add **2µl of RNase A** (100mg/ml) 1m
 RNase A (DNase free) 100mg/ml Qiagen Catalog #19101

15 **Incubate for 30 min at 37°C** 30m

16 Add 1 volume (10 ml) of **chloroform isoamyl alcohol** 1m
 Chloroform:Isoamyl Alcohol [24:1] Merck MilliporeSigma (Sigma-Aldrich) Catalog #25666
Vortex 2× 5 sec

17 **Incubate for 20 min at -20°C** 20m



18 **Centrifuge** with low deceleration

25m

8500 x g, 4°C, 00:20:00 , Acc 9 / Dec 6

19 Gently **transfer the aqueous phase** to a new 50 ml tube.

DNA précipitation

4h 2m

20 Beta-mercaptoethanol **Catalog #97622-10X1ML** Add to the supernatant :

1m

0.8 volume ($\approx$ 8 ml) of **isopropanol** Isopropanol

0.1 volume ($\approx$ 1 ml) of **3M sodium acetate**

Sodium acetate, 3M Fisher Scientific **Catalog #J61928**

0.2% ($\approx$ 2 μ l) of **β -mercaptoethanol**

Beta-mercaptoethanol **Catalog #97622-10X1ML**

21 **Homogenise by inversion and incubate for 1 hour at -20°C.**

1h

22 **Centrifuge** with low deceleration

25m

8500 x g, 4°C, 00:20:00 , Acc 9 / Dec 6

23 Carefully **discard the supernatant**, gently resuspend the pellet with **1 ml of cold 70% ethanol**.

5m

70% ethanol

24 **Transfer the resuspended DNA** to a 1.5 ml DNA LoBind tube.

5m

DNA LoBind Tubes 1.5 ml **Eppendorf Catalog #022431021**

25 **Centrifuge** with low deceleration

15m

5000 x g, 4°C, 00:10:00 , Acc 9 / Dec 6

26 Carefully **discard the supernatant**.
Air dry the pellet at RT for about 5-10 minutes.

10m

27 **Resuspend the pellet** in 50-100 μ l of 1X TE buffer.

1m

Buffer TE 1x



28 **Allow the DNA to resuspend** for 2 hours at 55°C or ON at RT.

2h

29 **Store the DNA at 4°C.**








Sample quality control

- 30 Quantify your sample with a **Qubit DNA HS assay**.
Check the purity of the sample with a **Nanodrop** (measurements of 260/280 and 260/230 absorbance ratios).
Estimate the molecular weight of the sample with a **Tapestation** and/or a **Femto pulse**.
- 31 Depending on the DNA concentrations and DNA length profiles, deplete short DNA molecules using **SRE size selection kits** (SRE XS, SRE or SRE XL kits).

Results

- 32 QC results obtained on different species.



Phylum	Taxon	DNA quantity µg	Extraction yield µg/g	Femto pulse
Rhodophyta	Gulsonia nodulosa	15	37,5	38 kb 
Rhodophyta	Dasya sp. 1	16	26,7	82 kb 
Rhodophyta	Lythophyllum incrustans	5,1	8,5	40 kb 
Ochrophyta	Sporochnus pedunculatus	4,3	15,0	14 kb 
Ochrophyta	Striaria attenuata	1,5	3,75	18 kb 
Chlorophyta	Cladophora laetevirens	21	52,5	41 kb 
Porifera	Cymbastela cantharella	62	310	73 kb 
Streptophyta	Prunus ramburei	5,6	14	128 kb 