


Oct 24, 2025

Version 2

HMW DNA extraction protocol for ferns v.2 V.2

DOI

<https://dx.doi.org/10.17504/protocols.io.bp2l69kb1lqe/v2>



Geferson Fernando Metz¹, Rafael Plá Matielo Lemos¹, Cristiane Barbosa D'Oliveira Matielo²,
Mariele Cristine Tesche Küster¹, Filipe Victoria¹

¹Núcleo de Estudos da Vegetação Antártica - NEVA (UNIPAMPA);

²Laboratório de Produtos Naturais e Fitoterápicos - LANAFi (UFRGS)



Rafael Plá Matielo Lemos

Universidade Federal do Pampa, Núcleo de Estudos da Vegetaçã...

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.bp2l69kb1lqe/v2>

Protocol Citation: Geferson Fernando Metz, Rafael Plá Matielo Lemos, Cristiane Barbosa D'Oliveira Matielo, Mariele Cristine Tesche Küster, Filipe Victoria 2025. HMW DNA extraction protocol for ferns v.2. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.bp2l69kb1lqe/v2> Version created by **Rafael Plá Matielo Lemos**

**Manuscript citation:**

Geferson Fernando Metz, Tiego Ferreira, Cristiane Barbosa D'Oliveira Matielo, Rafael Plá Matielo Lemos, Filipe de Carvalho Victoria. HMW DNA extraction protocol for ferns. protocols.io

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: October 23, 2025

Last Modified: October 24, 2025

Protocol Integer ID: 230592

Keywords: Ferns, Nanopore, HMW, long reads, Plant DNA, hmw dna extraction protocol for fern, disinfecting fern spore, fern spore, disinfection protocol for spore, hmw dna extraction protocol, dna extraction, fern, spore, read sequencing, germinated plant, disinfection protocol, large amounts of dna, sequencing, dna, dna in quantity, plant, extraction, growth free of contamination, extracting large amount, contamination, plants in the axenic state, quality hmw dna from recalcitrant plant tissue, quality hmw dna, preserving dna integrity, dna integrity, sequencing technology, antioxidants such as polyvinylpyrrolidone, many plant species due to high level, sodium dodecyl sulfate

Funders Acknowledgements:

Auxílio a Pesquisa GENO-ISLAND: Adaptações moleculares das plantas aos ambientes insulares

Grant ID: CNPQ 443237/2019-0

Coordenação de Aperfeiçoamento de Pessoal de Nível Superior - Brasil

Grant ID: CAPES 001

Fundação de Amparo à pesquisa do Estado do RS

Grant ID: 23/2551-0000849-1

Abstract

Obtaining high-molecular weight (HMW) DNA suitable for long-read sequencing remains challenging in many plant species due to high levels of secondary metabolites and the risk of mechanical shearing during extraction. To overcome these limitations, we developed an accessible, column-free protocol based on sodium dodecyl sulfate (SDS) lysis, optimized for preserving DNA integrity and purity. The method incorporates antioxidants such as polyvinylpyrrolidone, sodium metabisulfite, and β -mercaptoethanol to neutralize oxidative compounds, employs low NaCl concentration and potassium acetate precipitation to reduce co-purification of contaminants, and minimizes mechanical stress by avoiding vortexing and using wide-bore pipette tips. This optimized SDS-based approach provides a reproducible, cost-effective, and efficient solution for isolating high-quality HMW DNA from recalcitrant plant tissues, enabling reliable downstream applications in long-read sequencing technologies.

Image Attribution

Illustration by Crix D'Oliveira.

Guidelines

To reduce secondary metabolite production and optimize DNA extraction, we recommend placing the sample in a paper bag and storing it in the dark at 2–8 °C for approximately three days.

We also suggest using wide-bore pipette tips to minimize DNA shearing during handling.

For DNA extractions intended for PCR, the incubation times in Steps 18 and 22 may be reduced.

Materials

 Liquid nitrogen

 Chloroform: Isoamyl Alcohol (24:1)

 Isopropanol

 TE Buffer

 3M sodium acetate

 Ethanol 100%

 nuclease free water

 Sodium Hypochlorite Solution

 70% Ethanol

β-mercaptoethanol

PVP 40

TRIS 1 M pH 8.0

EDTA 0.5 M pH 8

Sodium Metabisulfite

NaCl 2.5M

Sodium dodecyl sulfate (SDS) 20%

2 ml LoBind tubes

1.5 ml LoBind tubes

Mortar & pestle

Water Bath

Centrifuge

Protocol materials

 70% Ethanol

 Sodium Hypochlorite Solution

 Liquid nitrogen

 Chloroform: Isoamyl Alcohol (24:1)

 Isopropanol

 TE Buffer

 3M sodium acetate

 Ethanol 100%

 nuclease free water

Safety warnings



Safety information

Work under fume hood when add β -mercaptoethanol and Chloroform.

Safety information

Be careful when handling liquid nitrogen.

Before start

- Prepare the SDS Lysis Buffer on the day of the experiment;

	A	B	C
	Reagent	Reagent Stock Concentration	FINAL
	PVP 40	100%	1%
	Sodium Metabisulfite	1%	1%
	NaCl	5 M	0,5 M
	TRIS HCL pH 8	1 M	100 mM
	EDTA pH 8	0.5 M	50 mM
	DDH ₂ O	-	~
	Sodium dodecyl sulfate (SDS)	20%	1,5%
	β-MERCAPTOETANOL	-	2% (v/v)

- When preparing the SDS Lysis Buffer, add SDS at last this will avoiding bubble formation.
- Preheat the water bath; Keep the SDS Lysis Buffer at 65°C until the tissue powder is added.
- Washing solution (EtOH 70%), fresh
35 ml Ethanol 100% + 15 ml H₂O
- Potassium Acetate 5M
Dissolve 4.9 gr KAc in 10 ml ddH₂O (4.9 gr KAc + ≈ 7.5 ml H₂O).
Adjust the pH with glacial acetic acid.
- Prepare TE buffer (10 mM Tris pH 8 and 1 mM EDTA pH 8)

Plant material sterilization

3d 0h 45m




Figure 1. Step by step plant spore and sporangia sterilization. Illustration by Crix D'Oliveira.

- 2 Sample and store the leaf tissue (fronds) in paper envelopes for three days to induce dehiscence.
- 3 Recover and store the spores together with sporangia in 1.5 mL microtubes until one-third of the tube was filled and then stored at -20°C until disinfection.

4 Add 1 mL $\otimes$ 70% Ethanol , homogenize by inversion for $00:00:20$ and briefly centrifuge them.

20s

5 Discard the supernatant and wash by inversion with 1 mL of autoclaved purified water, briefly centrifuge and discard the supernatant.

2m

6 Add 1 mL of an $\otimes$ Sodium Hypochlorite Solution (active chlorine $\sim 2\%$) at $10\% \text{ (v/v)}$, homogenize by inversion for $00:20:00$, centrifuge at 5.000 rpm $00:03:00$.

23m

- 7 Discard the supernatant and wash by inversion with 1 mL of autoclaved purified water, centrifuge at 5 5.000 rpm 00:03:00 . Repeat this step 4 times. 3m
- 8 Add 1 mL autoclaved purified water, homogenize and pipete into a Petri Dish with BCD medium (MgSO₄·7H₂O - 0.1 mM, KH₂PO₄ - 1.84 mM, KNO₃ - 1M, FeSO₄·7H₂O - 4.5mM). 1m

Extraction of high-molecular-weight DNA

2h 31m




Figure 2. Step by step extraction of high-molecular-weight DNA illustrated protocol. Illustration by Crix D'Oliveira.

- 10 Weigh between 50 mg and 100 mg of leaf tissue. Grind must be done with a crucible and pestle (previously exposed to Liquid nitrogen) until a very fine powder is obtained. 7m
- 11 Pre-heat the Lysis Buffer and aliquot 600 µL for each sample individually in 1.5 mL microtube. With the aid of a spoon or spatula, transfer the macerate to the microtube. 1m
- 12 Homogenize by inversion 10 times and incubate at 60 °C for 00:10:00 . 10m
- 13 Add 4 µL of RNase A (100 mg/mL); add 4 µL of proteinase K (> 40 U/mg);



- 14 Homogenize by inversion 10 times and incubate at 60 °C for 00:20:00 gently homogenizing by inversion every 00:05:00 . 25m
- 15 Add 600 µL of Chloroform: Isoamyl Alcohol (24:1) and homogenize by inversion for at least 00:03:00 until an off-white emulsion forms. 3m
- 16 Centrifuge at 14.000 rpm for 00:05:00 at room temperature. 5m
- 17 Carefully aspirate the upper phase of the tube and transfer to a new 1.5 mL microtube. 1m
- 18 Add 400 µL of Isopropanol -20 °C and incubate it for at least 01:00:00 at -20 °C . (Can be stored overnight) 1h
- 19 Centrifuge at 14.000 rpm for 00:05:00 at room temperature and discard the supernatant. 5m
- 20 Add 50 µL of TE Buffer . Add 5 µL of 3M sodium acetate . 1m
- 21 Add 120 µL of iced Ethanol 100% (-20 °C) and mix well by flicking the tube.
- 22 Store tubes at -20 °C for at least 00:20:00 . 20m
- 23 Centrifuge at 14.000 rpm for 00:15:00 at room temperature and carefully discard the supernatant. 15m
- 24 Add 1 mL of freshly prepared 70% Ethanol , mix gently by inversion, centrifuge briefly and carefully discard the supernatant. 1m



- 25 Dry the pellet for approximately  00:10:00 at Room Temperature on the bench with the tube upside down. 10m
- 26 Elute in  50 µL volume in  nuclease free water . (Elution volume can be adjusted for your needs) 1m

Quality Assessment

- 27 Quantify the DNA on a Qubit® fluorometer (dsDNA high sensitivity assay). DNA yield can be 500 - 1500 ng.
- 28 Check DNA integrity though electrophoresis in a 1% Agarose Gel.

DNA Size Selection

- 29 Remove short DNA fragments with Circulomics® Short-Read Eliminator Kit.

Protocol references

Metz GF, Ferreira TV, Ferreira RV, Matielo CBD, Lemos RPM, Victoria FC. The Complete Chloroplast Genome of Tree Fern *Cyathea delgadii* and Its Comparison to Other Cyatheales. *Biochem Genet.* 2025 Sep 11. doi: 10.1007/s10528-025-11248-3. Epub ahead of print. PMID: 40936070.

Mayjonade B, Gouzy J, Donnadiou C, Pouilly N, Marande W, Callot C, Langlade N, Muños S. Extraction of high-molecular-weight genomic DNA for long-read sequencing of single molecules. *Biotechniques.* 2016 Oct 1;61(4):203-205. doi: 10.2144/000114460. PMID: 27712583.