

Jan 14, 2025

 HNOE_15N_HSQC.nan

DOI

<https://dx.doi.org/10.17504/protocols.io.j8nlk88kwl5r/v1>



NAN KB¹, Alex Eletsy², John Glushka²

¹Network for Advanced NMR (NAN); ²University of Georgia



NAN-KB UGA

Network for Advanced NMR

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlk88kwl5r/v1>

Protocol Citation: NAN KB, Alex Eletsy, John Glushka 2025. HNOE_15N_HSQC.nan. protocols.io

<https://dx.doi.org/10.17504/protocols.io.j8nlk88kwl5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: May 29, 2024

Last Modified: January 14, 2025

Protocol Integer ID: 100901

Keywords: protein nmr, assignment, backbone, amides, ^{15}N , hsqc, hnoe hsqc experiment, hnoe hsqc experiment with sensitivity enhancement, hsqcnoef3gpsi3d, proteins with additional ^{13}C , saturation delay, topspin library pulseprogram, 2h labeling, protein, pulseprogram, d1 delay, hsqc experiment, hno, ^{15}N

Funders Acknowledgements:

NSF

Grant ID: 1946970

Disclaimer

This protocol is part of the knowledge base content of NAN: The Network for Advanced NMR (<https://usnan.nmrhub.org/>)

Specific filenames, paths, and parameter values apply to spectrometers in the NMR facility of the Complex Carbohydrate Research Center (CCRC) at the University of Georgia.

Abstract

This protocol describes running a Pseudo-3D ^{15}N $\{^1\text{H}\}$ HNOE HSQC experiment with sensitivity enhancement, gradient coherence selection, water flip-back, and optional ^{13}C decoupling.

Two datasets, reference and saturation, are acquired as a single pseudo-3D experiment, with TD1 fixed at 2. Both saturation delay and reference delay are determined internally from D1 delay.

^{15}N labeling is required; can also be applied to proteins with additional ^{13}C and ^2H labeling.

It uses the Topspin library pulseprogram 'hsqcnoef3gpsi3d'.

Attachments



[biotop.pdf](#)

2.5MB

Guidelines

Number of directly acquired points (2 TD) should be set to keep acquisition time $t_{2,max}$ (2 AQ) ~50–120 ms. For large proteins (~25 kDa) set closer to the lower range, higher values for small proteins.

^{15}N labeling is required; can also be applied to proteins with additional ^{13}C and ^2H labeling.

For samples with ^{13}C labeling use -DLABEL_CN flag to enable ^{13}C decoupling during ^{15}N evolution. ^{13}C channel offset O2P should be set ~110 ppm (middle of ^{13}C aliphatic and ^{13}CO shift range).

D1 delay serves both as saturation and recycle delay. For accurate results, set to at least 6 s. See paper by Renner et al. (<https://link.springer.com/article/10.1023/A:1015385910220>) for details

Before start

A sample must be inserted in the magnet either locally by the user after training, or by facility staff if running remotely.

This protocol requires a sample is locked, tuned/matched on ^1H , ^{13}C and ^{15}N channels, and shimmed. At a minimum, ^1H 90° pulse width and offset O1 should be calibrated and a 1D proton spectrum with water suppression has been collected according to the protocol **PRESAT_bio.nan**. Prior acquisition of a **2D ^{15}N HSQC** is also recommended, according to protocol **HSQC_15N.nan**.

General aspects of the use of the **bioTop** module are described in more detail in the protocols listed below.

It is recommended to calibrate ^1H carrier offset, ^1H H_2O selective flip-back pulse, as well as ^1H , ^{13}C , and ^{15}N 90° pulse widths using the "Optimization" tab of BioTop. Alternatively, ^1H 90° pulse width and offset can be calibrated using other methods, such as pulsecal or calibo1p1. Additional parameters, like ^{15}N and ^{13}C offsets and spectral width can be either optimized or manually entered in the "Optimization" tab of BioTop. Note that since BioTop optimizations are saved in the dataset folder, all experiments should be created under the same dataset name when using BioTop for acquisition setup.

Refer to protocols

1. Acquisition Setup Workflow, Solution NMR Structural Biology
2. PRESAT_bio.nan
3. HSQC_15N.nan
4. HNCO_3D.nan
5. Biotop-Calibration and Acquisition setup



Create HNOE_15N_HSQC experiment file

- 1 Start with existing dataset containing 1D proton or 15N-HSQC data for this sample.
 - 1.1 Open create dataset window with **edc** command or through menu.
Edit text in Title box.
 - 1.2 To select starting parameter set, check 'Read parameterset' box, and click Select.
For standard NAN parameter sets, change the Source directory at upper right corner of the window:
Source = /opt/NAN_SB/par
Click 'Select' to bring up list of parameter sets.
Select HNOE_15N_HSQC_xxx.par, where xxx=900,800 or 600.
Click OK at bottom of window to create the new EXPNO directory.
 - 1.3 If not done, tune Nitrogen and Carbon channels:
-on the Acquire menu, click Tune (or type atma on command line).

Calibrate pulses and adjust parameters

2

Note

Loading the HNOE_15N_HSQC_xxx.par, parameter set enters the default parameters into the experiment directory. While they represent a good starting point, they may not be fully optimal or accurate for your particular sample or spectrometer hardware. The probe- and solvent-specific parameters, specifically the ^1H 90° pulse length, and possibly the ^{13}C and ^{15}N 90° pulse lengths, along with other dependent pulse widths and powers may need to be updated.

There are two ways of automatically updating experimental parameters:

- 1) Use **getprosol** command, which typically only updates proton pulse widths and power levels. It is most useful for running routine experiments using the default parameters.
- 2) Or use **bioTop** module that organizes calibrated and defined parameters for a dataset.

2.1 Using getprosol:

Use the calibrated proton P1 value obtained from the proton experiment (protocol PRESAT_bio.nan) and note the standard power level attenuation in dB for P1 (PLW1); otherwise type calibo1p1 and wait till finished.

Then execute the getprosol command:

getprosol 1H [calibrated P1 value] [power level for P1]

e.g. getprosol 1H 9.9 -13.14.

Where for example, the calibrated P1=9.9 at power level -13.14 dB attenuation

This also loads 15N and 13C pulse widths and power levels from the 'edprosol' table, and are assumed to be sufficiently accurate.

Skip to step 2.3.**2.2 Use bioTop**

If you previously performed parameter calibrations using the "Optimizations" tab of the BioTop GUI, or entered parameters manually in the "Optimizations" tab, you can **type btprep** at the command line.

See protocol 'BioTop : calibration and acquisition parameters' and attached Bruker manual 'biotop.pdf' for details.

Inspect Parameters

- 2.3** Examine parameters by typing '**eda**', or select the 'Acqpars' tab and get the 'eda' window by clicking on the '**A**' icon. This view shows the three dimensions, F3 (1H) and F2 (15N) and F1 (2 points) in columns.

Parameters to check:

3 SW: 1H spectral width (~12-15 ppm, defined in BioTop)

2 SW: 15N amide spectral width (~25-40 ppm, defined in BioTop)

O1P: 1H offset defaults to 4.7ppm (calibrated with bioTop)

O2P: 13C offset (~110 ppm, for decoupling when **ZGOPTNS -DLABEL_CN**)

O3P: 15N amide offset (~115-120 ppm, defined in bioTop)

3 TD: Number of 1H time domain real points (~1024-2048, preferably 2^N , keep **3 AQ** at ~50-120 ms)

2 TD: Number of 15N time domain real points; keep 2 AQ at ~30-60 ms

1 TD: set to 2

- 2.4** Then examine the parameters in Acqpars 'ased' mode (click 'pulse' icon), or **type 'ased'**;

Most of the default parameters are suitable, however it's useful to compare values in the fields against those proposed by the original parameter file and the pulseprogram.

CNST4 - $^1J_{NH}$ value (≥ 93 Hz, used to calculate INEPT transfer delays)

CNST11 - multiplicity factor (NH: 4; NH/NH2: 8)

CNST12 - multiplicity factor (NH: 4; NH/NH2: 8)

D1: Recycle/Saturation delay. Set to >3 s (for accurate results should be >6s).

DS: 32-128 'dummy' scans that are not recorded ; allows system to reach steady state equilibration.

NS: multiple of 4; increase for increased signal to noise (S/N increases as $\sqrt{NS}$)

ZGOPTNS flags

-DLABEL_CN : enable ^{13}C decoupling during ^{15}N evolution for $^{13}\text{C},^{15}\text{N}$ -labeled samples (default)

channel f1:

P1 - ^1H 90° high power pulse (calibrated with calibo1p1 or BioTop)

SPdB1/SPW1 - power level [dB] for ^1H H_2O flip-back shaped pulse (calibrated with BioTop)

channel f2:

CNST21 - ^{13}CO offset (~173 ppm, calibrated with BioTop) for ^{13}C decoupling

CNST22 - ^{13}CA offset (~53 ppm, calibrated with BioTop) for ^{13}C decoupling

channel f3:

P21 - ^{15}N 90° high power pulse (calibrated with BioTop) at **PLW3**, high power

Acquire and Process Data

- 3 Type '**rga**' or click on 'Gain' in Topspin Acquire menu to execute automatic gain adjustment.
Type '**expt**' to calculate the expected run time
Type '**zg**' or click on 'Run' in Topspin Acquire menu to begin acquisition.
- 3.1 You can always check the first FID by typing 'efp' to execute an exponential multiplied Fourier transform. It will ask for a FID #, choose the default #1. You can evaluate the 1D spectrum for amide proton signal to noise and water suppression.
- 3.2 2D planes with NUS or uniformly sampled can be processed with Topspin: **xfb**

Protocol references

J.Cavanaugh, W.Fairbrother, A.Palmer, N.Skelton: Protein NMR Spectroscopy: Principles and Practice. Academic Press 2006 ; Hardback ISBN: 97801216449189, eBook ISBN: 9780080471037