



Dec 12, 2025

Holographic calibration

DOI

<https://dx.doi.org/10.17504/protocols.io.ewov11rqpv2/v1>

Ines Rodrigues-Vaz^{1,2,3}

¹Zuckerman Mind Brain Behavior Institute, Departments of Neuroscience and Neurology, Columbia University;

²Champalimaud Neuroscience Programme, Champalimaud Centre for the Unknown;

³Aligning Science Across Parkinson's (ASAP) Collaborative Research Network



Ines Rodrigues-Vaz

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.ewov11rqpv2/v1>

Protocol Citation: Ines Rodrigues-Vaz 2025. Holographic calibration. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.ewov11rqpv2/v1>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: December 02, 2025

Last Modified: December 12, 2025

Protocol Integer ID: 234021

Keywords: ASAPCRN, holographic calibration, bruker microscope, vivo ultima bruker microscope, photon, grin lens

Funders Acknowledgements:

Aligning Science Across Parkinson's

Grant ID: ASAP-020551

Abstract

This is the standard protocol for 2-photon holographic calibration with In vivo Ultima Bruker microscope through GRIN lens used for Rodrigues-Vaz and Athalye et al, 2025.

Materials

- Fluorescent slide
- Distilled water
- 2-photon laser scanning microscope with spatial light modulator module (Ultima In Vivo, Bruker)
- Femtosecond Ti:sapphire laser (Chameleon Vision-S, Coherent)
- Low repetition photostimulation laser (Monaco 1035-40, Coherent)

Safety warnings

 -Wear appropriate PPE as required by your institution.

Before start

-2-photon lasers (imaging and stimulation) and Pockels cell controller were turned on.



General calibration

- 1 Confirm imaging laser and photostimulation laser are aligned with target at objective output.
- 2 Place fluorescent sample under 16X water-immersion Nikon objective.
- 3 Focus on top of the sample.
- 4 Follow PrairieView software's calibration routine to ensure alignment between photostimulation and imaging lasers: burn spots 
- 5 Burn single spot and confirm calibration worked. If not, repeat step 4 as many times as necessary.

Power calibration

- 6 Measure power of stimulation laser at the 16X water-immersion Nikon objective output. This was done at least once at beginning of each cohort.