


Sep 24, 2019

Hornwort sporophyte induction - Bonn

DOI

dx.doi.org/10.17504/protocols.io.4x5gxq6



Eftychis Frangedakis¹

¹University of Cambridge



Eftychis Frangedakis

University of Cambridge, Plant Sciences

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.4x5gxq6

External link: <https://doi.org/10.1186/s12870-015-0481-x>

Protocol Citation: Eftychis Frangedakis 2019. Hornwort sporophyte induction - Bonn. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.4x5gxq6>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: June 29, 2019

Last Modified: September 24, 2019

Protocol Integer ID: 25309



Materials

G1910 - Sigma

Gelzan™ CM

Gelrite®

Synonym: Agar substitute gelling agent, Gellan Gum

BCD medium

Stock B (1L)

25g $\text{MgSO}_4 \times 7\text{H}_2\text{O}$

Stock C (1L)

25g KH_2PO_4

pH to 6.5 with KOH

Stock D (1L)

101g KNO_3

1.25g $\text{FeSO}_4 \times 7\text{H}_2\text{O}$

Stock CaCl_2 (1L)

14.7g $\text{CaCl}_2 \times 2\text{H}_2\text{O}$

Trace element solution (1L)

55mg $\text{CuSO}_4 \times 5\text{H}_2\text{O}$

614mg H_3BO_3

55mg $\text{CoCl} \times 6\text{H}_2\text{O}$

25mg $\text{NaMoO}_4 \times 2\text{H}_2\text{O}$

55mg $\text{ZnSO}_4 \times 7\text{H}_2\text{O}$

389mg $\text{MnCl}_2 \times 4\text{H}_2\text{O}$

28mg KI

BCD working solution (1L)

10ml Stock B

10ml Stock C

10ml Stock D

10ml Stock CaCl_2

1ml Trace element solution

[for plate: 8g agar]

- 1 Grow small thallus fragments for two weeks in petri dishes on Knop or BCD media at pH 5.7 and containing 0.7% (w/v) Gelzan (A in figure).
- 2 Transfer plants to Magenta pots on Knop medium or BCD medium (B and C in figure) at pH 5.7, containing 0.7% (w/v) Gelzan. Add 2mL of sterile water into the pot using a pipette.
- 3 Place pots in a Panasonic MLR-352 Versatile Environmental Test Chamber (or similar growth chamber / tissue culture room) at 21°C, 12 h of light and 12 h of dark, 1500 lux light intensity.
- 4 After approximately one month, antheridia start to appear. Add another 1-2 mL of sterile water into the pot using a pipette (D in figure).
- 5 After one more month sporophytes emerge (E in figure).
- 6

