



Mar 01, 2023

## How to make Tol2 mRNA

DOI

<https://dx.doi.org/10.17504/protocols.io.8epv5jo94l1b/v1>

FishFloorUCL<sup>1</sup>

<sup>1</sup>University College London

FishFloorUCL



Francois Kroll

Sorbonne Université

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.8epv5jo94l1b/v1>

Protocol Citation: FishFloorUCL 2023. How to make Tol2 mRNA. **protocols.io**  
<https://dx.doi.org/10.17504/protocols.io.8epv5jo94l1b/v1>

**License:** This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working



**Created:** March 01, 2023

**Last Modified:** March 22, 2023

**Protocol Integer ID:** 77877

**Keywords:** tol2 mrna, transcription kit, transcription,  $\mu$ l

## Abstract

These are instructions to make highly concentrated ( $> 1000 \text{ ng}/\mu\text{L}$ ) Tol2 mRNA. Note, the in-vitro transcription kit (mMESSAGE mMACHINE) is not cheap. 5–6 reactions like the protocol suggests come to  $\sim 100\text{--}120\text{€}$ , so please be thrifty with it.

## Linearise the Tol2 plasmid

- 1 Find the Tol2 plasmid from the Wilson lab freezer. It is #1151 and is labelled Tol2.

The concentration (measured on Qubit) is ~ 690 ng/μL (as of the Eppendorf in 2022).

Thaw it.

- 2 We linearise the plasmid using the restriction enzyme NotI HiFi.

On ice, prepare 6 reactions in 0.5 mL Eppendorfs.

Each reaction is (pipet in this order):

- Tol2 plasmid 2.9 μL
- NotI HiFi 2 μL
- CutSmart buffer 10× 10 μL
- nuclease-free H<sub>2</sub>O 91.1 μL

Total = 100 μL

(Add the enzyme last).

- 3 Place the Eppendorfs in the 37°C water bath.

Incubate for at least 1 hour. Digesting for more than 1 hour may be beneficial (cf. <https://international.neb.com/tools-and-resources/usage-guidelines/restriction-endonucleases-survival-in-a-reaction>).

- 4 Clean (remove the enzyme etc.) the reaction using the QIAquick PCR Purification Kit. Similar column-based kits are probably OK.

Follow the protocol in the box or online.

Use a single column.

For the first step (adding 5 volumes of PB buffer): transfer each 3 reactions in a 2 mL Eppendorf to obtain two 2 mL Eppendorfs, each containing 300 μL of reaction.

Add 1500 μL (5 volumes) of PB buffer. Vortex and spin down (with the benchtop centrifuge).

Add ~ 700 μL on the column and spin down following the kit's instructions (big centrifuge). Discard flow-through. Repeat until you got all 600 μL of reaction (the two 2 mL Eppendorfs) through the column.

- 5 At the last step of the kit's protocol, elute in 30  $\mu\text{L}$  nuclease-free  $\text{H}_2\text{O}$ .

Measure the concentration on the Qubit (dsDNA BR or HS kit). We expect the concentration to be 200–400  $\text{ng}/\mu\text{L}$  \*.

#### Note

\* Logic is: each reaction had ~ 2,000 ng of plasmid (2.9  $\mu\text{L}$  of 690  $\text{ng}/\mu\text{L}$ ) and we did 6 reactions, so 12,000 ng total. If you lost 0% on the column, you will get 12,000 ng in 30  $\mu\text{L}$ , i.e. 400  $\text{ng}/\mu\text{L}$ . If you lost 50% on the column, you will get 6,000 ng in 30  $\mu\text{L}$ , i.e. 200  $\text{ng}/\mu\text{L}$ .

## In-vitro transcription

- 6 Next step is to turn the linearised Tol2 plasmid (DNA) into mRNA. Find the SP6 mMESSAGE mMACHINE kit.

! Make sure you are using the SP6 version.

In my experience, each reaction generates 3,000–4,000 ng mRNA. Therefore, to get ~ 12  $\mu\text{L}$  of > 1,000  $\text{ng}/\mu\text{L}$ , you will need at least 3 reactions. This is assuming a good output and losing nothing on the column during clean-up, so I try to do 5–6 reactions for critical injections where I need high integration rate directly in F0 injected embryos.

Prepare each reaction on ice in a small 0.5 mL Eppendorf.

Each reaction is:

- 2 $\times$  NTP/CAP 10  $\mu\text{L}$
- 10 $\times$  reaction buffer 2  $\mu\text{L}$
- linearised Tol2 plasmid \*\*\*
- enzyme mix 2  $\mu\text{L}$
- nuclease-free  $\text{H}_2\text{O}$  \*to 20  $\mu\text{L}$ \*

\*\*\* Calculate from the concentration found in the previous step so that you use 1  $\mu\text{g}$ . e.g. I found the concentration to be 343  $\text{ng}/\mu\text{L}$ , so I use 2.91  $\mu\text{L}$  here.

- 7 Place the Eppendorfs in the 37°C water bath.

Incubate for 2 hours.

There is no benefit in incubating for more than 2 hours (cf. reaction time course: <https://www.thermofisher.com/order/catalog/product/AM1340#:~:text=mMESSAGE%20mMACHINE%E2%84%A2%20kits%20are,structure%20at%20the%205'%20end.>). Should be OK if reaction lasts a bit longer if for whatever reason you cannot stop it in time.

- 8 Add 1  $\mu\text{L}$  TURBO DNase in each reaction (included in the kit). This removes the DNA, i.e. the linearised plasmid.

Incubate in the 37°C water bath for 15 min.

- 9 Clean-up the mRNA using the ZYMO RNA Clean & Concentrator kit.

Follow the kit's instructions. Use a single column.

At the end, elute in 13  $\mu\text{L}$  nuclease-free H<sub>2</sub>O. Add the 13  $\mu\text{L}$  directly on the column's filter and wait for a few minutes at room temperature.

Heating the nuclease-free H<sub>2</sub>O to 70°C prior to adding it to the water may also help with retrieving as much as possible mRNA from the column, but I have never tried.

(You can always dilute later if you find that the Tol2 mRNA is highly concentrated, so best to elute in a small volume.)

- 10 Measure the concentration on Qubit (BR or HS RNA kit). Qubit can measure low concentrations, so no need to waste precious Tol2 mRNA. You can e.g. dilute 10 $\times$ :
  - 0.4  $\mu\text{L}$  Tol2 mRNA
  - 3.6  $\mu\text{L}$  nuclease-free H<sub>2</sub>Oand measure that concentration. Just remember to multiply the concentration you find by 10.

You are expecting (hoping) > 1,000 ng/ $\mu\text{L}$ . For reference, if each of 6 reactions produced 3000 ng and nothing is lost on the column, you would expect ~ 1384 ng/ $\mu\text{L}$  (6  $\times$  3000 ng, eluted in 13  $\mu\text{L}$ ).

- 11 Aliquot the ~ 13  $\mu\text{L}$  in 3–4 small 0.5 mL Eppendorfs and store in the –80°C freezer.