

Dec 19, 2025

hPBMCs thawing and plating (no specific down-stream assay)

DOI

<https://dx.doi.org/10.17504/protocols.io.n2bvj1ow5vk5/v1>

Rebecca Wallings¹, Malu Tansey¹

¹Indiana University



Janna Jernigan Posey

Indiana University

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.n2bvj1ow5vk5/v1>

Protocol Citation: Rebecca Wallings, Malu Tansey 2025. hPBMCs thawing and plating (no specific down-stream assay). protocols.io <https://dx.doi.org/10.17504/protocols.io.n2bvj1ow5vk5/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 19, 2025



Last Modified: December 22, 2025

Protocol Integer ID: 235487

Keywords: ASAPCRN, human pbmcs from cyrovial, plating human pbmc, hpbmc, stream assay, cyrovial

Funders Acknowledgements:

The Michael J. Fox Foundation for Parkinson's Research (MJFF) and the Aligning Science Across Parkinson's (ASAP)

Grant ID: ASAP-020527

Abstract

General protocol for thawing and plating human PBMCs from cyrovials.

Materials

- Water bath at 37 degrees (use on in main lab)
- Dry ice
- 1 × 50mL falcon tube per cryovial, with 25mL of thawing media at 37 degrees
- MACs buffer at 37 degrees
- Thawing media at 37 degrees
- Plating media at 37 degrees
- Tissue culture plasticware: 24-well plates, p1000 pipette tips, p20 pipette tips, stripettes, aspirating stripettes, trypan blue, countess slides, eppendorfs

Protocol

- 1 Retrieve cryovials for plating from Mr W and place on dry ice.
- 2 Make sure you have 1 falcon with thawing media per cryovial ready in the hood, pre-warmed to 37 degrees.
- 3 Place cells in floaty and place in water bath and allow to thaw for approx. 90 seconds.
- 4 Check to see if thawed, a small amount of ice remaining is fine.
- 5 Pour contents of cryovial into falcon containing thawing media.
- 6 Centrifuge falcons at room temp for 10 minutes at 300 x g.
- 7 Aspirate supernatant and resuspend pellet in 10mL MACs buffer – make sure resuspend pellet first in 1mL of warmed MACs with a p1000 pipette to ensure pellet is broken up, then add 9mL MACs with stripette.
- 8 Take 10uL of cell suspension and place in Eppendorf.
- 9 Centrifuge falcons 300 x g 10 mins at room temp.
- 10 Whilst falcons are centrifuging proceed with cell count: mix 10uL of trypan blue with 10uL of cell suspension, place 10uL in both chambers of a countess slide and count using “hPBMC cryo” programme on countess.
- 11 Note cell count per mL (countess accounts for 1:1 trypan blue dilution, so this number is accurate already) and live cell % (aiming for 80+ % viability).
- 12 Calculate volume of plating media needed for each pellet to obtain 1×10^6 cells per mL.



- 13 Example: countess count = 5×10^5 per mL. $5 \times 10^5 \times 10$ (as cells were in 10mL) = 5×10^6 cells total. Therefore resuspend cells in 5mL for 1×10^6 per mL.
- 14 For each patient sample, plate 500uL of cell suspension per well of a 24-well plate, and plate 1 well per treatment required for experiment.
NOTE: This can be adjusted based on downstream assay being used. E.g. flow cytometry = 1 million/mL and plate 500uL in 24 well plate (results in 500,000 per well), whilst RNA needs more cells so plate 2mL in 6 well plate (results in 2 million per well).
- 15 For each patient, label wells with patient ID and treatments.
- 16 Place in incubator for MINIMUM of 2 hours to rest and proceed with downstream assay.