



Aug 15, 2024

Version 1

## hPSC Passaging and Propagation on laminin521 V.1

DOI

[dx.doi.org/10.17504/protocols.io.81wgbz9e3gpk/v1](https://dx.doi.org/10.17504/protocols.io.81wgbz9e3gpk/v1)

Gist Gist Croft<sup>1</sup>, Regine Tipon<sup>1</sup>, Niraj Sawarkar<sup>2</sup>, Sigi Benjamin<sup>2</sup>

<sup>1</sup>New York Stem Cell Foundation; <sup>2</sup>NYSCF



Sigi Benjamin

nyscf

### Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



**DOI:** <https://dx.doi.org/10.17504/protocols.io.81wgbz9e3gpk/v1>

**Protocol Citation:** Gist Gist Croft, Regine Tipon, Niraj Sawarkar, Sigi Benjamin 2024. hPSC Passaging and Propagation on laminin521. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.81wgbz9e3gpk/v1>

**Manuscript citation:**

Ruzo A, Croft GF, Metzger JJ, Galgoczi S, Gerber LJ, Pellegrini C, Wang H Jr, Fenner M, Tse S, Marks A, Nchako C, Brivanlou AH. Chromosomal instability during neurogenesis in Huntington's disease. *Development*. 2018 Jan 29;145(2):dev156844. doi: 10.1242/dev.156844. PMID: 29378824.



**License:** This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

**Protocol status:** Working

**We use this protocol and it's working**

**Created:** August 14, 2024

**Last Modified:** August 15, 2024

**Protocol Integer ID:** 105309

**Keywords:** ASAPCRN, PSC, iPSC, hESC, substrate, laminin521, edta, CEPT, epithelial monolayer culture, using laminin521, laminin, confluent single cell layer, cell line, laminin521, gentle cell dissociation reagent, cell

**Funders Acknowledgements:**

**Alining Science Across Parkinson's**

**Grant ID:** ASAP-000472

## Abstract

hPSC Passaging and propagation using Laminin521 and EDTA. Laminin actively supports survival, prevents spontaneous differentiation, and increases plating efficiency, thus net expansion, and enables 100% confluent single cell layer epithelial monolayer culture. EDTA-based (Gentle Cell Dissociation Reagent) minimizes stress and enhances viability. Cell lines require little or no adaptation.

## Attachments



[EDTA \(Gentle Cell\) P...](#)

95KB



## Guidelines

**Survival** in as single cells is supported by ROCKi at 10uM, which blocks apoptosis downstream of loss of cell-to-cell E-cadherin signalling. We routinely use 20uM at plating to ensure maximal survival. Thiazovivin 1uM is equivalent to ROCKi 10uM. CEPT cocktail is superior to either. CloneR (StemCell Tech) is best for single cell cloning survival but ingredients and mechanism is not published. Standard spit is 1:10, i.e. **50K/cm<sup>2</sup>**, e.g. ~0.5M cells/6wp well). ROCKi or equivalent lasts for 2 days. On day 2 after passage if cells have not reached self supporting quorums (>5 cells/colony) re-feed with ROCKi supplemented medium

**Seeding density:** split cells at 1:10 or 1:20 bby area to a new well or count and seed at 0.05M cells/cm<sup>2</sup>. Unlike other substrates cells can be seeded at >0.05M/cm<sup>2</sup> without ROCK inhibitor Y-27632 (ROCKi) since they will migrate to form colonies before they die. Lower seeding density does requires ROCKi even on laminin.

**IMPORTANT:** to ensure **normalized cell seeding density**, it is important to mix cell suspension immediately (<30sec) before removing aliquot for counting, and again immediately before seeding. Single cells soft-pellet quickly by gravity. If seeding takes too long, periodically mix full volume of cell suspension. If seeding multiple plate types, adjust **master seeding stock** to a fixed concentration and calculate the volume to add per well type well.

### Confluent density:

Cells on laminin 521--unlike other substrates--can be brought to 100% confluency without differentiation. They will pack in at 0.5M cells/cm<sup>2</sup> in a completely flat epithelial sheet with apical pole up, continue dividing, but not differentiate, resembling a homogeneous honeycomb appearance. We typically passage at 100% confluence or just before, but can passage up to 4 days post-confluence as long as no differentiaton is observed (it should not be).

**Feeding volume and frequency:** Seed cells in at least 2.5 ml/6wp well. during the first half week, if using StemFlex, days may be skipped. As cells proliferate, increase feeding volume and feed daily. After 1 week, feed volume is 7-8ml/6wp well and 8-9ml if 100% confluent.

**Passage interval:** generally ~1/week. see confluent density note above.

### Other medias:

This protocol is compatible with any common hPSC medium, conditioned medium, MTsr1 (plus), conditioned medium, etc. Increase feed frequency and do not skip-days feeding if media is not additionally pH buffered and has stabilized FGF (like StemFlex)

### Other substrates

We have not tested other suppliers laminin521 substrates. Passaging protocol works well for substrates other than laminin521 (geltrex or vitronectin). ROCKi or equivalent is required for single cell passage on these substrates. Passage before colonies are too large ( crowding in the center, or vertical growth) or begin to touch (can spur differentiation).

## Materials

	MATERIAL	SUPPLIER	CATALOG NUMBER	USE
	Gentle Cell Dissociation Reagent	STEMCELL Technologies	100-0485	Cell dissociation
	Human Recombinant Laminin 521	Biolamina	LN521-05	Plate Coating
	Y-27632 dihydrochloride	R&D Systems	1254/50	ROCK inhibition
	StemFlex	Thermo	A33493	PSC media
	CEPT cocktail kit	Tocris	7991	PSC survival

## Troubleshooting



## Plate Coating

- 1 **PREPARE LAMININ WORKING SOLUTION:** Dilute laminin521 stock (100ug/ml) to 5 ug/mL in PBS with Mg<sup>2+</sup> and Ca<sup>2+</sup>, (1:20 dilution). **IMPORTANT:** Mg and Ca cations in PBS required for laminin function, do not use PBS without Mg<sup>2+</sup> and Ca<sup>2+</sup>



### 1.1 VOLUMES OF WORKING SOLUTIONS REQUIRED FOR DIFFERENT PLATES

REAGENT	4wp (0.5 mL/well)	6wp (1.5 mL/well)	24wp (0.25 mL/well)
Laminin521	0.11 mL = 110 uL	0.5 mL = 500 uL	0.4 mL = 400 uL
PBS(Final Volume for 1 plate with spare)	2.2 mL	10 mL	8 mL

- 2 Coat laminin521: minimum 2 hrs at 37 degrees. Plates may also be coated at 4 degrees sealed 1-2 days, but warm in the incubator 1 h before use.

\*\*\* Used laminin521 may be stored at 4 degrees for 1 month and reused 3 x.

**IMPORTANT:** Re-supplement with 1:40 (2.5 ug/mL) fresh laminin at each re-use.

## Cell Passaging

- 3 Aspirate medium
- 4 Wash with PBS -/- (with 25mM glucose, optional)
- 5 Add **500**-1000 uL (per 24wp well) / 1.5 - 2.0 mL (per 6wp well) of Gentle Cell Dissociation Reagent (GCDR) or Accutase
- 6 Incubate at 37 degrees till cells slough with gentle touch, 10-30 min for GCDR, usually 5-10 min for accutase. (check every 10 min by moving the plate gently)
- 7 Pipet with p1000 to remove cells (>90%), and transfer to a 15 mL tube containing DMEM





- 8 Gently triturate to single cells (and homogenate for cell counting/seeding)
- 9 Rinse well 1x with fresh 1mL GCDR and pool to 15ml (repeat if necessary)
- 10 Spin  $300 \times g$  for 5 minutes
- 11 Aspirate and resuspend in 1 ml PBS -/- (+glucose if available) or GCDR (**steps 8-10 are optional**)
- 12 Top up to appropriate counting volume (5 mL, ~1M cells/mL), homogenize and remove 20  $\mu$ L for count - prepare counting sample in an eppendorf tube before applying to the slide.
- 13 Spin  $300 \times g$  for 5 min
- 14 Count while spinning
- 15 Resuspend at appropriate [ ]\*\* in StemFlex + Suppl  
  
\*\* standard split ration for weekly passage is 1:10-1:20 =0.025-0.05M cells/cm<sup>2</sup>
- 16 If using laminin521, collect/save coating solution right before seeding and store at 4 degrees. Do not let substrate dry.\*\*\*



## Protocol references

**Chromosomal instability during neurogenesis in Huntington's disease.** Albert Ruzo<sup>1</sup>, Gist F Croft<sup>1</sup>, Jakob J Metzger<sup>1,2</sup>, Szilvia Galgoczi<sup>1</sup>, Lauren J Gerber<sup>1</sup>, Cecilia Pellegrini<sup>1</sup>, Hanbin Wang Jr<sup>1</sup>, Maria Fenner<sup>1</sup>, Stephanie Tse<sup>1</sup>, Adam Marks<sup>1</sup>, Corbyn Nchako<sup>1</sup>, Ali H Brivanlou  
Development. 2018 Jan 29;145(2):dev156844. doi: 10.1242/dev.156844.

### **CEPT cocktail**

Chen Y., Tristan C.A., Chen  
L., Jovanovic V.M., Malley C., Chu P.H., Ryu S., Deng T., Ormanoglu P.,  
Tao D., et al. A versatile polypharmacology platform promotes  
cytoprotection and viability of human pluripotent and differentiated  
cells. Nat. Methods. 2021;18:528–541. doi: 10.1038/s41592-021-01126-2.