

Dec 20, 2024

 HSQC_CT_13C_aro.nan

DOI

<https://dx.doi.org/10.17504/protocols.io.14egn689yl5d/v1>



NAN KB¹, Alex Eletsy², John Glushka²

¹Network for Advanced NMR (NAN); ²University of Georgia



NAN-KB UGA

Network for Advanced NMR

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.14egn689yl5d/v1>

Protocol Citation: NAN KB, Alex Eletsy, John Glushka 2024. HSQC_CT_13C_aro.nan. **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.14egn689yl5d/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: April 12, 2024

Last Modified: December 20, 2024

Protocol Integer ID: 98127

Keywords: protein nmr, assignment, backbone, hsqc, carbon-13, constant-time, C13, side-chain, aromatic 2d 13c trosy, 2d constant time 13c hsqc, time 13c hsqc, 13c resonance, signals for each aromatic carbon, aliphatic 13c, aromatic side chain spin system identification, aromatic 3d spectra, hsqc, resolution in the aromatic region, required isotope labeling, chain aromatic 1h, isotope labeling, aromatic carbon, spin system identification routine sample screening, 13c, aromatic region, aromatic 3d, aromatic 1h polarization, monitoring chemical shift perturbation study, pulse sequence with sensitivity enhancement, 3d spectra, proton pair, spectrometer, sensitive probe, spin system identification, chemical shift perturbation studies due to ligand

Funders Acknowledgements:

NSF

Grant ID: 1946970

Disclaimer

This protocol is part of the knowledge base content of **NAN: The Network for Advanced NMR** (<https://usnan.nmrhub.org/>)

Specific filenames, paths, and parameter values apply to spectrometers in the NMR facility of the Complex Carbohydrate Research Center (CCRC) at the University of Georgia.

Abstract

This protocol describes running a 2D constant time ^{13}C HSQC pulse sequence with sensitivity enhancement, gradient coherence selection, and water flipback. It is optimized for high-resolution in the aromatic region. This produces a 2D phase-sensitive ^{13}C - ^1H correlation spectrum that displays signals for each aromatic carbon-proton pair.

Required isotope labeling: U- ^{15}N , ^{13}C or U- ^{13}C . Not suitable for samples with uniform ^2H labeling due to starting with aromatic ^1H polarization. Application to samples at natural abundance is possible but would require high concentration, sensitive probes and lengthy signal averaging.

Optimal MW is less than 25 kDa. Larger systems may benefit from aromatic 2D ^{13}C TROSY.

Field strength preference: High B_0 fields are preferred for better signal dispersion and minimization of strong coupling artifacts.

This pulse sequence can be used for:

- resonance assignment of side-chain aromatic ^1H and ^{13}C resonances (in His, Phe, Trp and Tyr)
- aromatic side chain spin system identification
- routine sample screening and stability monitoring
- chemical shift perturbation studies due to ligand binding, paramagnetic relaxation or pseudo-contact shifts
- optimization of aliphatic ^{13}C offset and spectral width for use with aromatic 3D ^{13}C -edited NOESY/TOCSY, (H)CCH-COSY, and other experiments.
- anchoring of aromatic 3D spectra (^{13}C -edited NOESY/TOCSY, (H)CCH-COSY, etc.) during interactive visual analysis

It uses the pulseprogram 'hsqcctetgpsisp_aro.nan' (hsqc= heteronuclear single quantum correlation; ct=constant time; et=echo-antiecho; si=sensitivity enhanced; sp=shaped 180 pulses) modified from the original Bruker library sequence.

Note: specific parameter values illustrated below may differ depending on the facility and spectrometer.

Guidelines

The number of directly acquired points (**2 TD**) should be set so the acquisition time $t_{2,\max}$ (**2 AQ**) is between ~50 ms (for larger proteins ~25 kDa) and ~120 ms (for smaller proteins). Longer times may cause excessive probe and sample heating during ^{13}C decoupling, and resolve undesirable $^3\text{J}_{\text{H,H}}$ splittings.

Total constant-time delay is equal to $2 \times \text{D23}$, and is normally set to multiples of $(1/J_{\text{CC}})$. Here $^1J_{\text{CC}}$ is assumed to be 57 Hz, though the true splittings may vary slightly from site to site. The most common settings are $(1/J_{\text{CC}})$ and $(2/J_{\text{CC}})$, corresponding to **D23** of 0.0088 and 0.00166 s, respectively. In the former case peak signs will depend coupled spin multiplicity (even or odd number of covalently linked ^{13}C spins), while in the latter case all peaks have the same sign.

Since the F1(^{13}C) is a constant-time dimension, number of acquired points **1 TD** should be set to the largest possible value for optimal S/N and resolution. The maximum **1 TD** value is $2 \times \text{int}(d20/\text{in}20)$, where $d20$ is the decremented delay, and $\text{in}20$ is the decrement value. Both $d20$ and $\text{in}20$ are computed internally according to pulse program code and their values are shown in ased view. Note that $\text{in}0 = \text{in}20 = 1/(2 \times \text{SW}(^{13}\text{C}))$. In TopSpin 4.x the ased command would produce a warning if **1 TD** is set too large and reset the internal $td2$ parameter to the maximum possible value.

NUS sampling is usually not required for 2D experiments, since time savings are small, unless running multiple 2D experiments. If using NUS keep sampling amount ~50%.

Before start

A sample must be inserted in the magnet either locally by the user after training or by facility staff if running remotely.

This protocol requires a sample is locked, tuned/matched on ^1H , ^{13}C and ^{15}N channels, and shimmed. At a minimum, ^1H 90° pulse width and offset O1 should be calibrated and a 1D proton spectrum with water suppression has been collected according to the **protocol PRESAT_bio.nan**.

It is recommended to calibrate ^1H carrier offset, ^1H H_2O selective flip-back pulse, as well as ^1H , ^{13}C , and ^{15}N 90° pulse widths using the "Optimization" tab of **bioTop**. Alternatively, ^1H 90° pulse width and offset can be calibrated using other methods, such as **pulsecal** or **calibo1p1**. Additional parameters, like ^{15}N and ^{13}C offsets and spectral width can be either optimized or manually entered in the "Optimization" tab of **BioTop**.

Note that since BioTop optimizations are saved in the dataset folder, all experiments should be created under the same dataset name when using BioTop for acquisition setup.

Familiarize yourself with the general workflow for NMR study of a protein sample is outlined in protocol "Aquisition Setup Workflow, Solution NMR Structural Biology".

Create aromatic ¹³C CT-HSQC experiment

- 1 Start with existing Dataset containing 1D proton data in EXPNO 1.
(see step 1)
- 1.1 Click on Acquire → 'Create Dataset' button to open dataset entry box.
or type **edc** command.

Dataset: The Name defaults to the one used for the previous 1D dataset.

If a different name is desired, fill this out.

The EXPNO should increment by +1

Directory should be the same as preliminary 1D.

The Title text box will be the same as for EXPNO 1.


Edit to designate an HSQC pulse program and add other details as appropriate.

- 1.2 Select starting parameter set:
Check 'Read parameterset' box, and click Select.

For standard NAN parameter sets, change the Source directory at upper right corner of the window:

Source = /opt/NAN_SB/par

Select HSQC_CT_13C_ARO_xxx.par, where xxx=900,800 or 600*.




* parameter set name may differ depending on spectrometer

- 1.3 Click OK at bottom of window to create the new EXPNO directory.
It will be the active experiment in the acquisition window and should now be listed on your data browser.

1.4 **Tune Carbon (and Nitrogen) channels if not done.**

Return to the Acquire menu and click Tune (or type atma on command line).



This will start tuning of the nitrogen channel, then the carbon channel, followed by a retuning of the proton channel (which should not change).

Load pulse calibrations: use getprosol (step 2.1) or bioTop (step 2.2)

2

Note

Loading the HSQC_CT_13C_ARO_xxx.par parameter set enters the default parameters into the experiment directory. While a good starting point, they may not be fully optimal or accurate for your particular sample or spectrometer hardware. The probe- and solvent-specific parameters, specifically the ^1H 90° pulse length, and possibly the ^{13}C and ^{15}N 90° pulse lengths, along with other dependent pulse widths and powers may need to be updated.

There are two ways of automatically updating an entire range of experimental parameters. The first is using **getprosol** command (step 2.1), which only updates pulse widths and power levels without altering other parameters, such as spectral widths and offsets. This method is suitable for reproducing existing experiments or parameters sets with minimal variations.

The second method utilizes the **BioTop** module of TopSpin (step 2.2), and can load additional experimental parameters, such as spectral widths, offsets, and number of time-domain points. These additional parameters are set according to calibrations or definitions within the 'Optimization' tab of the BioTop GUI and the corresponding XML description files (*bt_hsqcctetfpgps_i_aro.nan.xml* in this case). This method has a lower dependency on the particular settings of the starting parameter set, and is suitable for setting up experiments from scratch. With this method nearly all important acquisition parameters can be optimized for a particular sample, and then applied consistently to multiple NMR experiments with a single command.

2.1 Loading pulse widths and power levels with **getprosol**:

Use the calibrated proton P1 value obtained from the proton experiment (protocol PRESAT_bio.nan) and note the standard power level attenuation in dB for P1 (PLW1); otherwise type calibo1p1 and wait till finished.

Then execute the getprosol command:

getprosol 1H [calibrated P1 value] [power level attenuation for P1 (PLdB1)]

e.g. getprosol 1H **9.9** -13.14.

Where for example, the calibrated **P1=9.9** at power level -13.14 dB attenuation

This also loads default ¹⁵N and ¹³C pulse widths and power levels from the PROSOL table, and are assumed to be sufficiently accurate.

⇒ [go to step #2.3](#) If not using BioTop

2.2 Loading experimental parameters from **BioTop**:

If you previously performed parameter calibrations using the "Optimizations" tab of the BioTop GUI, or entered parameters manually in the "Optimizations" tab, you can simply type **btprep** at the command line.

This is equivalent to calling **getprosol** with all ¹H, ¹³C and ¹⁵N optimized parameters followed by additional parameters loaded from "Optimizations" tab in BioTop GUI. In this particular case these parameters are based on the *bt_hsqcctetgpsisp_ali.nan.xml* description file: ¹H offset in Hz (**O1**), ¹H spectral width (**2 SW**), ¹⁵N offset in ppm (**O3P**), half constant-time (CT) delay (**D23**), ¹³C spectral width (**1 SW**), ¹³C max acquisition time (**1 AQ**), ¹³C offset in ppm (**O2P**).

Inspect and adjust parameters

2.3 The default parameters from HSQC_CT_13C_ARO_xxx.par will provide an aromatic 2D 13C CT-HSQC spectrum of a typical protein sample collected with the traditional sampling scheme (i.e. not using non-uniform sampling 'NUS'). Often the only parameters to change will be **NS = number of scans** in order to increase the signal to noise, and **1 TD** (TD in F1) that changes the number of increments (points) and hence the digital resolution in the 13C dimension.

Check the experiment time (type expt or click on 'clock' icon) after any change.

These and additional parameters can be accessed and changed on the parameter windows seen below.

- 2.4 Select the 'Acqparams' tab to display acquisition parameters. Two display modes can be selected, the full display mode (click on the '**A**' icon or type **eda**), or pulse program-specific mode (click on the 'pulse' icon, or type **ased**). The former gives you access to all parameters and provides an overview of all spectral dimensions at once, while the latter is useful because it only displays acquisition parameters used in the pulse sequence and can be parsed sequentially as a checklist.

First examine the specific dimension parameters in Acqparams 'eda' mode (click 'A' icon). This view gives a comprehensive list of parameters covering two dimensions, F2(1H) and F1(13C) in columns.

Experiment	F2	F1	Frequency axis
Experiment PULPROG: hsqcctetgpsi_aro.nan ... E Current pulse program AQ_mod: DQD Acquisition mode FnTYPE: traditional(planes) nD acquisition mode for 3D etc. FnMODE: Echo-Antiecho Acquisition mode for 2D, 3D etc. TD: 1024 Size of fid DS: 64 Number of dummy scans NS: 4 Number of scans TD0: 1 Loop count for 'td0' TDav: 0 Average loop counter for nD exper			
Width			
SW [ppm]	12.6292	40.0000	Spectral width
SWH [Hz]	11363.636	9051.226	Spectral width
IN_F [μsec]		110.4875	Increment for delay
AQ [sec]	0.0450560	0.0151361	Acquisition time
FIDRES [Hz]	22.194603	66.067345	Fid resolution
FW [Hz]	240000000.000		Filter width
Nucleus 1			
NUC1	1H	13C	Observe nucleus
O1 [Hz]	4229.65	28281.54	Transmitter frequency offset
O1P [ppm]	4.701	125.000	Transmitter frequency offset
SFO1 [MHz]	899.7942297	226.2806565	Transmitter frequency
BF1 [MHz]	899.7900000	226.2523750	Basic transmitter frequency
Nucleus 2			

Parameters to check:

- **FnTYPE** - 'traditional planes' or 'non-uniform sampling' - see step 3 below
- **SW[ppm]** - F2(1H) ~12-15 ppm (defined in bioTop); F1(13C) ~40 ppm
- **O1** - ¹H H₂O offset in Hz (calibrated with BioTop or calibo1p1, **O1P** will be around 4.7 ppm depending on temperature)
- **O2P** - ¹³C aromatic offset (~125 ppm)
- **O3P** - ¹⁵N offset (~130 ppm, for decoupling ¹⁵N in Trp and His rings)
- **2 TD** - Number of ¹H time domain real points (~1024-2048, preferably 2^N, keep **2 AQ** at ~50-120 ms)
- **1 TD** - Number of ¹³C time domain real points (int(2*d20/in20) for maximum signal-to-noise)
- **DIGMOD** - 'baseopt' (zero 1st order phase correction)

2.5 Then examine the parameters in the pulse program-specific '**ased**' mode (click on the 'pulse' icon). Most parameters are also accessible in the '**eda**' mode (step 2.4 above). However, the '**ased**' mode allows more convenient access to individual parameters within arrays, such as delays, pulse widths, constants, etc. It also displays parameter values computed internally within a pulse sequence, and provides context description from the relevant pulse program comment lines.

- **CNST2** - effective one-bond $^1J_{CH}$ coupling value (≥ 165 Hz for aromatic groups); used to calculate INEPT transfer delays
- **D1** - 1- 2 sec
- **D23**: constant-time delay 8.8 ms
- **P1** - 1H 90° high power pulse (calibrated with calibo1p1 or BioTop)
- **P3** - ^{13}C 90° high power pulse (calibrated with BioTop)
- **P21** - ^{15}N 90° high power pulse (calibrated with BioTop)

General		General			
Channel f1	PULPROG	hsqcctetgpsl_aro.nan	...	E	Pulse program for acquisition
Channel f2	TD	1024			Time domain size
Channel f3	SWH [Hz, ppm]	11363.64	12.6292		Sweep width in acquisition direction
Gradient channel	AQ [sec]	0.0450560			Acquisition time
F1 indirect dimension	RG	101			Receiver gain
	DW [μsec]	44.000			Dwell time
	DE [μsec]	18.00			Pre-scan-delay
	CNST2	167.0000000			= J(XH)
	d0 [sec]	0.0000030000			Incremented delay (2D)
	D1 [sec]	1.000000000			Relaxation delay: 1.5 * T1
	d4 [sec]	0.0014970060			1/(4)*XH
	d11 [sec]	0.0299999993			Delay for disk I/O
	D16 [sec]	0.000200000			Delay for homospoil/gradient recovery
	d20 [sec]	0.0075862501			= d23
	D23 [sec]	0.008800000			d23 = T : 8.8 or 17.6 msec
	DELTA1 [sec]	0.0004635916			DELTA1=p29+d16-p1*0.78+de+4u
	DELTA2 [sec]	0.0006832561			DELTA2=d4-larger(p2,p4)/2-p19-d16
	DELTA3 [sec]	0.0086529749			DELTA3=d23-d0-p4/2-larger(p2,p21*4.5)
	DELTA4 [sec]	0.0006970061			DELTA4=d4-p19-d16
	DS	64			32
	F1CNT	2			F1CNT = min(2 , td1)
	in0 [sec]	0.0000552437			1/(2 * SW(X)) = DW(X)
	in20 [sec]	0.0000552375			= in0
	INF1 [μsec]	110.4875			1/SW(X) = 2 * DW(X)
	NS	4			4 * n



Channel f1		
SFO1 [MHz]	899.7942297	
O1 [Hz, ppm]	4229.65	4.701
NUC1	1H Edit...	
P1 [μsec]	10.780	
p2 [μsec]	21.56	
PLW1 [W, dB]	20.627	-13.14
Frequency of ch. 1		
Frequency of ch. 1		
Nucleus for channel 1		
F1 channel - 90 degree high power pulse		
F1 channel - 180 degree high power pulse		
F1 channel - power level for pulse (default)		
Channel f2		
SFO2 [MHz]	226.2806565	
O2 [Hz, ppm]	28281.54	125.000
NUC2	13C Edit...	
CPDPRG 2	garp ... E	
P3 [μsec]	13.750	
p4 [μsec]	27.50	
PCPD2 [μsec]	50.00	
PLW2 [W, dB]	51.394	-17.11
PLW12 [W, dB]	3.8844	-5.89
Frequency of ch. 2		
Frequency of ch. 2		
Nucleus for channel 2		
File name for cpd2		
F2 channel - 90 degree high power pulse		
F2 channel - 180 degree high power pulse		
F2 channel - 90 degree pulse for decoupling sequence		
F2 channel - power level for pulse (default)		
F2 channel - power level for CPD/BB decoupling		
Channel f3		
SFO3 [MHz]	91.1867817	
O3 [Hz, ppm]	11852.74	130.000
NUC3	15N Edit...	
P21 [μsec]	28.950	
PLW3 [W, dB]	69.542	-18.42
Frequency of ch. 3		
Frequency of ch. 3		
Nucleus for channel 3		
Pulse 21		
F3 channel - power level for pulse (default)		
Gradient channel		
GSNAM 1	GMSG10 1.00 ... E	

Acquire and Process Data

- 3 Type '**expt**' to calculate the expected run time.

 [go to step #2.3](#) If necessary to re-adjust parameters

Type '**rga**' or click on 'Gain' in Topspin Acquire menu to execute automatic gain adjustment.

Type '**zg**' or click on 'Run' in Topspin Acquire menu to begin acquisition.

- 3.1 You can always check the first FID by typing 'efp' to execute an exponential multiplied Fourier transform. It will ask for a FID #, choose the default #1. You can evaluate the 1D spectrum for amide proton signal to noise and water suppression.

To take a look at the 2D, wait for ≥ 16 FIDS and then click on 'Proc.Spectrum' on the Topspin Process menu. This will execute an automated processing macro.



Although the resolution will be poor, you can evaluate the signal to noise (S/N) and whether the ^{15}N offset and spectral width are suitable.

This can be repeated at any time as additional FIDs are acquired.

Protocol references

J.Cavanaugh, W.Fairbrother, A.Palmer, N.Skelton: Protein NMR Spectroscopy: Principles and Practice. Academic Press 2006 ; Hardback ISBN: 97801216449189, eBook ISBN: 9780080471037

G.W. Vuister & A. Bax, J. Magn. Reson. 98, 428-435 (1992))

A.G. Palmer III, J. Cavanagh, P.E. Wright & M. Rance, J. Magn. Reson. 93, 151-170 (1991)