



Feb 24, 2020

HuBMAP: Embedding Fixed Frozen OCT Samples

DOI

dx.doi.org/10.17504/protocols.io.basniede

Marda Jorgensen¹, Jerelyn Nick¹

¹University of Florida

Human BioMolecular Atlas Program (HuBMAP) Method Development Community

Tech. support email: Jeff.spraggins@vanderbilt.edu



Marda Jorgensen

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.basniede

Protocol Citation: Marda Jorgensen, Jerelyn Nick 2020. HuBMAP: Embedding Fixed Frozen OCT Samples. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.basniede>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working

Created: December 20, 2019

Last Modified: April 02, 2020

Protocol Integer ID: 31278

Abstract

The purpose of this Standard Operating Procedure (SOP) is to outline procedures for the OCT embedding of HuBMAP frozen fixed specimens.



Guidelines

- Managers and supervisors - are responsible for making sure that technicians are properly trained and equipment and facility are maintained in good working order.
- Laboratory personnel - are responsible for reading and understanding this SOP and related documents and to perform these tasks in accordance with the SOPs.

Materials

MATERIALS

☒ KimWipes **Fischer Scientific**

☒ Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International (Avantor) Catalog #25608-930**

☒ 16% Paraformaldehyde **Fisher Scientific Catalog #15710**

☒ D-Sucrose (Molecular Biology) Fisher BioReagents **Fisher Scientific Catalog #BP220-1**

☒ PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher Scientific Catalog #BP399-1**

☒ Tissue Tek Cryomold (25mmx20mmx5mm) **Catalog #25608-916**

☒ 2-Methylbutane **Fisher Scientific Catalog #O3551-4**

☒ Ice / Dry Ice Bucket (EVA Foam) **Fisher Scientific Catalog #03-395-152**

STEP MATERIALS

☒ Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International (Avantor) Catalog #25608-930**

☒ PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher Scientific Catalog #BP399-1**

16% Paraformaldehyde

PBS 10x Solution

D-Sucrose

OCT

Tissue Tek Cryomold

Methylbutane

Dry Ice

Ice Bucket



Protocol materials

☒ 16% Paraformaldehyde **Fisher Scientific Catalog #15710**

☒ D-Sucrose (Molecular Biology) Fisher BioReagents **Fisher Scientific Catalog #BP220-1**

☒ PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher Scientific Catalog #BP399-1**

☒ PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher Scientific Catalog #BP399-1**

☒ KimWipes **Fischer Scientific**

☒ Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International (Avantor) Catalog #25608-930**

☒ Ice / Dry Ice Bucket (EVA Foam) **Fisher Scientific Catalog #03-395-152**

☒ Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International (Avantor) Catalog #25608-930**

☒ Tissue Tek Cryomold (25mmx20mmx5mm) **Catalog #25608-916**

☒ 2-Methylbutane **Fisher Scientific Catalog #O3551-4**

☒ PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher Scientific Catalog #BP399-1**

☒ Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International (Avantor) Catalog #25608-930**

Safety warnings

- ! ■ Use physical safety precautions when working with sharps (disposable blades).

Before start

- Ensure you have proper scalpel blades, forceps, and your personal preference of gauzes/wipes.
- Embedding can be a messy process, to protect your clothes it is best to wear a lab coat or apron.
- Gloves are highly suggested to protect your fingers from the spectrum of hot and cold one encounters during the process.

Reagent Preparation

1h

1

Prepare **1x PBS**

PBS is required to prepare reagents needed for this procedure. Dilute 10x PBS 1:10 in deionized water to make working 1x solution.



PBS Phosphate Buffered Saline 10X Solution Fisher BioReagents **Fisher**
Scientific Catalog #BP399-1

2

Prepare a **4% paraformaldehyde solution**.

Safety information

Paraformaldehyde is toxic. Use gloves, wear a lab coat and work in a fume hood.

2.1

Add 60ml of 100mM PBS to a sample cup.



60 mL PBS

Note

Prepare containers to accomodate number and size of tissues using atleast a 1:20 tissue to fixative volume ratio.

2.2

Open 2 -10mL ampules of 16% paraformaldehyde stock.
Add contents to the 1x PBS in the sample cup, invert to mix.



20 mL 16% Paraformaldehyde Stock

2.3

Label containers, solution expires one week after being made. Store at 4° C if not immediately used.

- 1 L 1x PBS

3d 8h

- 30m

PREVIEW



Document



NAME

SOP Appendix for Thymus

CREATED BY

Marda Jorgensen

PREVIEW

Document



NAME

SOP Appendix for Lymph Node

CREATED BY

Marda Jorgensen

PREVIEW

-

- 5 Place tissue into pre-labelled cassettes maintaining orientation.

15m

- 6 Transfer the tissue into a 120mL specimen container, pre-filled with 80mL freshly prepared fixative (4% paraformaldehyde in PBS).

1d

Fix specimen for 24 hours at room temperature under mild agitation using a rocker.

24:00:00 In Fixative on Rocker

Room temperature

- 7 Drain fixative (collect as hazardous waste) and wash the tissue **three (3) times** for 10 minutes each in 1X PBS at a volume that generously covers. Invert several times during this process, or return to rocking platform.

1h

🕒 00:30:00

- 8 Infiltrate and cryoprotect the tissue with 30% sucrose.

2d

🌡️ 4 °C in Sucrose/PBS Mixture

🕒 72:00:00

- 8.1 Combine equal volumes of 30% sucrose and 1x PBS to form 15% Sucrose. Replace PBS wash with 15% sucrose and place at 4° C for 8 hours

- 8.2 Drain 15% sucrose and replace with 30% sucrose. Equilibrate tissue 4° C for at least 48 hours.

Tissue is stable in sucrose for at least one week.

Freeze the Blocks

30m

- 9 Prepare a pre-labeled Cryomold

1m

Equipment

Cryomold®

NAME

Histology

TYPE

Tissue Tek®

BRAND

25608-916

SKU

<https://us.vwr.com/store/product/4639407/tissue-tek-cryomold-molds-adapters-sakura-finetek>

LINK


25mm x 20 mm x 5mm

SPECIFICATIONS

and fill it half way with OCT compound.



Tissue-Tek® O.C.T. Compound, Sakura® Finetek **VWR International**
(Avantor) Catalog #25608-930

- 10 Remove the tissue from the cassette using forceps and quickly touch the tissue to a Kimwipe to remove external sucrose droplet from tissue surface. 1m
- 11 Place the tissue into the OCT-containing Cryomold, maintaining original tissue orientation. 1m
- 12 Using forceps, push the tissue lightly to the bottom of the Cryomold to secure it. 1m
- 13 Prepare dry ice/methylbutane slurry for freezing OCT blocks.
- 13.1 Place 1-2 inches of dry ice pellets into the bottom of an ice bucket or styrofoam box.
- 13.2 Add enough 2-methylbutane solution to cover the dry ice by roughly 5mm .
- 13.3 Place lid on freezing chamber and allow methylbutane to chill.
Chamber is ready when fog dissipates and bottom of bucket becomes visible.
- 14 Freeze the tissue in the Cryomold by resting it on the surface of the methylbutane slurry. As the OCT inside the Cryomold begins to freeze, lightly push the tissue into the bottom of the mold one last time.
- 15 Add additional OCT to cover the tissue completely and fill the Cryomold. Allow the tissue to equilibrate in the mold for ten (10) minutes. 20m
 00:10:00 In Cryomold
- 16 When the OCT Cryomold is completely frozen and opaque, wrap the mold containing the tissue in a pre-labeled aluminum foil square, and store at -80° C in a freezer rack box.



Note

Enter block location in Freezer Log for future retrieval