

Jan 27, 2023

Human axillary lymph node fine-needle aspirate sample processing and cryopreservation

DOI

dx.doi.org/10.17504/protocols.io.j8nlkw36wl5r/v1

Jacqueline HY Siu¹, Katrina M Pollock¹, calliope.dendrou¹

¹University of Oxford



Jacqueline HY Siu

University of Oxford

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.j8nlkw36wl5r/v1>

Protocol Citation: Jacqueline HY Siu, Katrina M Pollock, calliope.dendrou 2023. Human axillary lymph node fine-needle aspirate sample processing and cryopreservation. **protocols.io** <https://dx.doi.org/10.17504/protocols.io.j8nlkw36wl5r/v1>

License: This is an open access protocol distributed under the terms of the **Creative Commons Attribution License**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: Working

We use this protocol and it's working



Created: January 13, 2023

Last Modified: January 27, 2023

Protocol Integer ID: 75197

Keywords: FNA, lymph node, LN, fine-needle aspirate, human axillary lymph node, needle aspirate sample, needle aspirate samples of the adult, needle aspirate sample processing, cell suspension, cyropreservation this protocol, cyropreservation, cell yield

Funders Acknowledgements:

Chan Zuckerberg Initiative

Grant ID: 2021-239944

Abstract

This protocol was used to generate a single-cell suspension from fine-needle aspirate samples of the adult human axillary lymph nodes and cryopreserved in CS10 for < 6 months. Cell yield was highly dependent on the participant; however, viability was >90% prior to freezing and between 80-90% after thawing (trypan blue dye exclusion).

Guidelines

- All steps should be performed at 4°C unless stated otherwise.

Troubleshooting

Safety warnings

- ❗ Sample preparation should be carried out in a Class II microbiological safety cabinet in a designated Containment Level 2 blood handling facility. Centrifuge steps should be performed in a designated blood-handling centrifuge with aerosol-tight inner lids.



Before start

Just prior to starting:

- Pre-cool centrifuge to 4°C
- Pre-cool cryopreservation vials in the -20°C
- Pre-cool cell freezing container if necessary

Reagents:




- **R10 Media** is made up with RPMI 1640 with 25 mM HEPES (Cat No: R5886, Sigma), 10% heat inactivated fetal bovine serum (HI-FBS) (Cat No: F4135, Sigma), 1% Pen/Strep (Penicillin-Streptomycin 10,000 Units/mL (Cat No: 15140-122, Gibco), and 2 mM of L-Glutamine (Cat No: 25030-024, Gibco). Sterile bottle 500 mL filter system 0.22 µm (Cat No: 430758/ 430769 Corning or Cat No: SEGPU0538/ SEGPU0545, Merck) or 0.22 µm syringe filters. Store at 4°C when not in use for a maximum of 1 month.
- **Cell Wash Buffer** is made up with PBS with 2% human AB serum and 2 mM EDTA. Store at 4°C when not in use.



Sample collection & transfer


- 1 Collect lymph node FNA samples in sterile filtered R10 media in 15 mL Falcon tubes. Transfer samples in an appropriate container and keep the tubes at 4°C using cooled gel packs.

Sample processing

- 2 Upon sample arrival, top up tubes with cold RPMI + 5% human AB serum (serum has been previously heat inactivated and sterile filtered) as required to ensure even volumes. 1m
- 3 Centrifuge at 400xg at 4°C for 10 min.  400 x g, 4°C, 00:10:00 10m
- 4 Remove the supernatant, resuspend in 5 mL of red blood cell (ACK) lysis buffer (Gibco, Cat: A10492-01). Incubate at room temperature for 5 min (check colour), up to 10 min max. Top up with cell wash buffer (PBS + 2% FBS, 2 mM EDTA).  Room temperature 10m
- 5 Centrifuge at 400xg at 4°C for 10 min. Wash again with cell wash buffer (PBS + 2% FBS, 2 mM EDTA).  400 x g, 4°C, 00:10:00 10m
- 6 Adjust final volume to 1 ml of cold RPMI + 5% AB serum. 1m
- 7 Take a 10µl aliquot of the cell suspension and dilute with 10µl Trypan Blue. Load this mixture onto a haemocytometer and perform a cell count. 5m

Sample cryopreservation

10m

- 8 Label the tubes and chill at -20°C for ~ 10 minutes. The minimum number should be ~500,000 cells (to be frozen in 100 µl; 100 µl is the minimum freezing volume). Aliquot cells such that there is a maximum of ~1 million cells per cryovial.
- 9 Add up to 10 ml of cold RPMI + 5% hAB serum to the tubes. Centrifuge at 400xg for 10 min at 4°C.  400 x g, 4°C, 00:10:00 10m



- 10 Resuspended samples in CS10 medium at a concentration of $1 \times 10^6/100 \mu\text{l}$. 1m
- 11 Once fully mixed, aliquot $100 \mu\text{l}$ of the sample into the chilled cryotubes. 1m
- 12 Transfer the cryotubes into appropriate cell freezing container. Ensure all slots of the MrFrosty are filled, using the filler vials if necessary. 2m
- 13 Transfer the samples into liquid nitrogen. Ideally this transfer should be performed within 24 hr, but may be extended to a maximum of 72 hr. 🕒 Overnight