

Sep 29, 2025

Human iPSC-derived MgBr assembloids

 In 1 collection

DOI

<https://dx.doi.org/10.17504/protocols.io.3byl46x7ogo5/v1>

Maria Jose Perez J.¹

¹Institut Imagine

Team Deleidi



Federico Bertoli

ASAP

Create & collaborate more with a free account

Edit and publish protocols, collaborate in communities, share insights through comments, and track progress with run records.

Create free account

OPEN  ACCESS



DOI: <https://dx.doi.org/10.17504/protocols.io.3byl46x7ogo5/v1>

Protocol Citation: Maria Jose Perez J. 2025. Human iPSC-derived MgBr assembloids. **protocols.io**
<https://dx.doi.org/10.17504/protocols.io.3byl46x7ogo5/v1>

License: This is an open access protocol distributed under the terms of the **[Creative Commons Attribution License](#)**, which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited



Protocol status: Working

We use this protocol and it's working

Created: September 25, 2025

Last Modified: September 29, 2025

Protocol Integer ID: 228165

Keywords: derived mgbr assembloids human ipsc, mgbr assembloids human ipsc, human ipsc, ipsc

Funders Acknowledgements:

ASAP

Abstract

Human iPSC-derived MgBr assembloids

- 1 Seed 9,000 iPSCs into each well of a Nunclon™ Sphera™ 96-well microplate (Cat. number 174925; Thermo Fisher Scientific) to generate embryoid bodies (EBs).
- 2 Culture EBs from day 0 to day 5 in dual SMAD inhibition medium consisting of mTeSR™ Plus supplemented daily with 1% NEAAs, 1% GlutaMAX, 100 μM β-mercaptoethanol, 200 nM LDN 193189 hydrochloride (LDN), 10 μM SB, and 50 μM Y-27632 (ROCK inhibitor).
- 3 On day 6, switch to cortical induction medium consisting of 1:1 DMEM-Ham's F-12 and Neurobasal medium supplemented with 1% NEAAs, 1% GlutaMAX, 1% N2, 1% B27 without vitamin A, 100 μM β-mercaptoethanol, 20 ng/mL FGF2, and 20 ng/mL EGF.
- 4 On day 7, embed EBs in Matrigel at a 3:2 ratio of Matrigel to medium to form the embedding mixture.
- 5 Wash EBs with fresh medium, embed in Matrigel in a 6-well ultralow-attachment plate, and incubate the Matrigel-EB mixture for 30 minutes at 37°C.
- 6 Maintain organoids in cortical induction medium until day 14.
- 7 From days 14 to 30, maintain organoids in differentiation medium consisting of 100% Neurobasal medium supplemented with 1% NEAAs, 1% GlutaMAX, 1% N2, 1% B27 without vitamin A, 1% PS, 100 μM β-mercaptoethanol, 2.5 μg/mL insulin, 200 μM ascorbic acid, 20 ng/mL BDNF, and 1 mM dcAMP.
- 8 On day 20, dissociate Matrigel from organoids and place plates on an orbital shaker to increase oxygenation and improve medium distribution.
- 9 On day 30, transfer individual brain organoids into wells of a Nunclon™ Sphera™ 96-well microplate to generate iPSC-derived MgBr assembloids.
- 10 On day 10 of in vitro differentiation, seed 2×10^5 predifferentiated iPSC-derived microglia on top of each organoid in a final volume of 200 μL assembloid maturation medium.
- 11 Prepare assembloid maturation medium with 100% Neurobasal medium supplemented with 1% NEAAs, 1% GlutaMAX, 1% N2, 1% B27, 1% PS, 100 μM β-mercaptoethanol, 20 ng/mL BDNF, 1 mM dcAMP, 10 ng/mL GDNF, 100 ng/mL IL-34, 100 ng/mL M-CSF, and 10 ng/mL GM-CSF.
- 12 Change maturation medium twice daily for 3 days without removing microglia in suspension.



- 13 Transfer assembloids to a 6-well ultralow-attachment plate after 3 days.

- 14 After 2 additional days, place MgBr assembloids on an orbital shaker and maintain in maturation medium until day 45 (15 DPI) or day 65 (35 DPI) of differentiation.