



Jul 09, 2018

🌐 Human Kidney / Tumour Tissue Disaggregation for Single Cell RNA Sequencing (10x Genomics platform)

DOI

dx.doi.org/10.17504/protocols.io.qtrdwm6

Kevin Loudon¹, John Ferdinand¹, Alexandra Riding¹, Menna Clatworthy¹

¹University of Cambridge



John Ferdinand

University of Cambridge

OPEN  ACCESS



DOI: dx.doi.org/10.17504/protocols.io.qtrdwm6

Protocol Citation: Kevin Loudon, John Ferdinand, Alexandra Riding, Menna Clatworthy 2018. Human Kidney / Tumour Tissue Disaggregation for Single Cell RNA Sequencing (10x Genomics platform). **protocols.io**

<https://dx.doi.org/10.17504/protocols.io.qtrdwm6>

License: This is an open access protocol distributed under the terms of the [Creative Commons Attribution License](#), which permits unrestricted use, distribution, and reproduction in any medium, provided the original author and source are credited

Protocol status: In development

We are still developing and optimizing this protocol

Created: June 08, 2018

Last Modified: July 09, 2018

Protocol Integer ID: 12881

Materials

STEP MATERIALS

⊗ Dead Cell Removal Kit **Miltenyi Biotec** Catalog #130-090-101

⊗ Dead Cell Removal Kit **Miltenyi Biotec** Catalog #130-090-101



Protocol materials

⊗ Dead Cell Removal Kit **Miltenyi Biotec Catalog #130-090-101**

⊗ Dead Cell Removal Kit **Miltenyi Biotec Catalog #130-090-101**

⊗ Dead Cell Removal Kit **Miltenyi Biotec Catalog #130-090-101**

Tissue Preparation

- 1 Take dissected tissue (renal cortex, medulla or tumour) and weigh tissue (typical biopsy size used 0.5 - 1 gram)
- 2 Pour approximately 2-3 mLs of "Digest Mix" onto sample in 10cm³ petridish

Note

PREPARATION OF DIGEST MIX

Ingredients

- (1) RPMI alone
- (2) Liberase TM (Sigma Aldrich)
- (3) DNase (Sigma Aldrich)

For 50mLs of RPMI add:

- > 625 microlitres of Liberase (Stock solution 2.5mg/mL)
- > 250 microlitres of DNase (Stock solution 0.05mg/mL)

- 3 Using a razor blade mince into small pieces approximately 2mm³.
- 4 Transfer tissue into a gentleMACS C tube and add further 3-4 mLs of Digest mix.

Note

GentleMACS C tube by Miltenyi Biotec (Cat.130-096-334)

- 5 Place in shaking incubator at 37°C for 30 minutes.
- 6 Homogenise sample in GentleMACS tube using program "Spleen 4" and "Lung 2" on GentleMACS dissociator.
- 7 Pass through a 100µm cell strainer with of a 2.5ml syringe plunger and wash through with cold running buffer.

Note

PREPARATION OF RUNNING BUFFER**Ingredients (for 1 litre)**

- (1) 1L PBS
- (2) 5ml BSA (from reagent diluent kit)
- (2) 4ml 0.5M EDTA

- 8 Centrifuge in a bench top centrifuge at 2000 RPM for 10 minutes and CAREFULLY remove the supernatant.
- 9 If sample is contaminated with red blood cells an additional red cell lysis step can be taken.
- 10 To ensure optimal yield for 10X Genomics single cell platform, a live cell enrichment step is required - this was performed using Miltenyi 'Dead Cell Removal Kit' (Please see manufacturers instructions for further details).

Note

LIVE CELL ENRICHMENT (Miltenyi - Dead Cell Removal Kit)**Ingredients**

- (1) Dead Cell removal Kit - Miltenyi (Order No. 130-090-101)
- (2) MACS Column (LS or MS)

In brief for MACS colum LS

- (1) Use LS column for 10^8 dead cells or 10^9 total cells.
- (2) Remove supernatant completely following previous steps
- (3) Resuspend pellet in 100 μ L of '*Dead Cell Removal MicroBeads*' per approximately 10^7 total cells.
- (4) Incubate 15 minutes at room temperature (20–25 °C).
- (5) Rinse column with 1x binding buffer as per manufacturers instructions.
- (6) Apply cell suspension in 1-10mLs of binding buffer and collect the effluent as the **NEGATIVE** cell population (i.e the live cells).
- (7) Wash cells with PBS for 5 minutes at 1500rpm.



Dead Cell Removal Kit Miltenyi Biotec Catalog #130-090-101

- 11 Count the cells and resuspend the live cell suspension in appropriate volume of PBS for the 10X application.